
· HEINE EL10 LED Examination Light

HEINE EL10 LED Examination Light

DEUTSCH	3
ENGLISH	9
FRANÇAIS	15
ESPAÑOL	21
ITALIANO	27
SVENSKA	33
NEDERLANDS	39
NORSK	45
SUOMI	51
PORTUGUÊS	57

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Please read and follow these instructions for use and keep them for future reference.

General Conditions of Warranty

Instead of the statutory warranty time period of 2 years, HEINE will grant a guarantee of 5 years from the date of the consignment of the goods ex works, concerning its equipment (excluding disposables, e.g. bulbs, single-use articles, and rechargeable batteries). Date of consignment means that HEINE hands the goods over to the transport carrier, freight forwarder or any other person designated by the Customer for the transport of the goods without loading the collecting vehicle.

The guarantee covers irrevocable workmanship, on condition of the proper use of the equipment and the observation of the operating instructions. During the warranty- and guarantee time period, errors and deficiencies arising on the equipment will be rectified free of charge, in so far as such are evidenced by defective materials, processing and/or constructional errors. Should buyer complain of a material deficiency during the warranty time period, then the onus of proof is always to be on the orderer, that the product was defective already upon receipt of the goods. The statutory warranty and the guarantee do not apply to loss or damage caused by wear and tear, negligent use, the non-employment of original HEINE components and/or spares (in particular bulbs, as these have been especially developed for HEINE instruments in accordance with the following criteria: colour temperature, useful service life, safety, optical quality and performance. The statutory warranty and the guarantee do not apply to interventions by persons not authorised by HEINE or when the operating instructions are not observed by the customer. Any modification of a HEINE product with parts or additional parts which do not conform to the original HEINE specification will invalidate the warranty for the correct function of the product and further invalidate any warranty claims which result from such a change or modification. Further claims, in particular claims for replacement of loss or damage, which are experienced otherwise than directly on the HEINE product itself, are hereby excluded.

Warnings and Safety Information

Caution! Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow, foreground color black.)

Note! Note indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

Intended use

HEINE EL10 LED is a mains operated medical examination light for illuminating the body's surfaces and cavities during medical examinations in closed examination rooms.

The device is intended only for medical examinations in which failure of the light will not pose a risk to the patient. A main source of lighting must always be switched on.


Product overview

Gooseneck:

- 1 Gooseneck (flexible component)
- 2 Gooseneck (rigid component)
- 3 Gooseneck (flexible component)

Light head:

- 4 Hand grip
- 5 Brightness adjuster and ON/OFF switch (illumination)
- 6 Focussing ring

Light foot:

- 7 ON/OFF switch (system)
- 8 Connection socket for the power cable
- 9 Installation cover
- 10 Connection for the mounting systems

Installation instructions

When mounting the accessories, please follow the enclosed mounting instructions for the mounting systems.

Caution! To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

To ensure all-pole disconnection from the mains at all times, the ME device must be installed so that the power cord is accessible and disconnectable.

When using the examination light with a stand, clamp or wall mount, ensure the examination light is standing or stably attached.

Check that the attached examination light is stable and secure before each use.


Use the accessories supplied with the examination light to ensure the examination light is held securely and properly in the respective mount.

When mounting the examination light on a stand, ensure that the requirements of IEC 60601-1 relating to the risk of instability as well as casters and wheels are observed.

When moving the light, ensure that the power cable is not taut. Unplug the power cable when moving the light over longer distances.





When mounting the light on the wall, ensure that the subsurface is able to bear the mechanical stress.

The wall plugs supplied with the product are universal wall plugs and suitable for most building materials (e.g. concrete, solid brick, plasterboard or gypsum cardboard). Before mounting the device, check whether you need to use a special wall plug for the subsurface of your wall. Before mounting the device on the wall, check that there are no cables or pipes running under the planned drill holes that could get damaged.

For mounting the mains adaptor to the wall, clamp or stand mount, push up the cover (9) and take it off. Use the screws supplied with the product to fasten the examination light to the respective mounting system (clamp or stand mount) or use the wall plugs supplied to mount the light on the wall.

The screws for the mounting systems are located in the corners of the mounting unit.

1 Clamp mount:

Please use the screws enclosed with the product for mounting the light on the light foot.

This can be fastened to an even surface (e.g. table plate) with a maximum thickness of 40 mm or to a circular tube with a maximum diameter of 30 mm. Please check that the mounted light is securely in place before using the light.

2 Stand mount:

This mount system is designed for tubes with a diameter of 25 mm.

Please use the screws enclosed with the product for mounting the light on the light foot.

3 Wall mounting:

The rubber mat can be used as a drilling template for marking the holes on the wall and is intended to level out uneven areas on the wall for the mains adaptor.

Please use the screw and universal wall plug combination supplied with the product if the subsurface is suitable to do so.

Use

⚠ The HEINE EL10 LED must not be used for eye examinations or as a surgical light.

The HEINE EL10 LED is intended for use in closed spaces. It is not designed for mobile use or for use outdoors. Improper use can result in damage to the device or pose a risk to the user.

Prior to each use, check that the HEINE EL10 LED is functioning properly as intended. Do not use the device if visible damage is apparent.

The device does not have an application part and therefore care must be taken to ensure that the examination light does not come into contact with the patient!

HEINE EL10 LED must not be used in strong magnetic fields such as within the vicinity of MRI equipment

The HEINE EL10 LED must not be used in oxygen-enriched environments.

Fire hazard! The HEINE EL10 LED must not be used in an environment of flammable anaesthetic agents (class AP) or with oxidising agents (class APG).

Risk of burns! Always keep a minimum distance of 20 cm from the light.

Risk of burns! Do not cover the examination light.


Ensure that the device is handled properly and carefully.

Do not let the power cable become taut as this could damage the device or pose a tripping hazard.

Only after the device has been securely mounted, should the power cable supplied with the device be connected by plugging it into the mains connection socket (8). To use the device, switch it on by pressing the power switch (7) so that the switch is in the ON position I and the switch lights up green.

Positioning

⚠ To avoid contamination, position the examination light prior to the examination as well as during the examination using the assistance.


The swivel arm of the HEINE EL10 LED examination light has two flexible (1, 3) and one rigid section (2).

To position the light precisely, hold the light head by the hand grip (4) and the gooseneck by the flexible part (1) and bend the gooseneck into the desired position. Fine adjustments can be made via the hand grip (4).

Always use the indicated handling points to move or position the examination light.


To prevent any glare when the device is switched on, reduce the brightness to a minimum before switching off the examination light or switch it off using the brightness adjuster.

Do not bend the device more than 90 degrees.

Operation

⚠ The optimal working distance is 30 cm. The minimum distance of 20 cm must not be exceeded.


Brightness control:

To individually adjust the illuminance, turn the brightness adjuster (5) to the right (see figure).

To switch on the light source, turn the brightness adjuster to the ☺ symbol.

To switch it off, turn the adjuster until you see the symbol. ☹

The actual illuminance is indicated on the indicator bar

Focus adjustment:

To individually adjust the size of the light field (see figure) turn the focussing ring (6).

The smallest light field is indicated by the symbol * and the largest light field by the symbol ●.

Hygienic Cleaning / Processing

⚠ To prevent cross-contamination, the touched contact surfaces and any other contaminated surfaces must be cleaned and disinfected after each use.

The examination light must only be cleaned and disinfected once it has cooled down and is disconnected from the power supply. The device must be disconnected from the mains supply for this. It must be switched off at the power supply and unplugged from the mains.

The examination light may only be cleaned and disinfected by a member of the hygiene personnel or a person instructed by the hygiene personnel.

Occupational safety must be assured in all process steps, e.g. by wearing appropriate protective clothing, goggles and gloves, and ensuring room air quality.

Only agents that have been tested and approved by HEINE Optotechnik with regard to material compatibility and effectiveness of disinfection may be used.

If other procedures / disinfectants are used the fundamental effectiveness of the procedure must be proven as part of a successful validation process.

Spray disinfection is not permitted due to the lack of effectiveness and associated health hazards, and as this may also damage the product.


The application times, environmental temperatures and concentrations of the disinfectants vary and must be applied in accordance with the manufacturer's instructions. The manufacturer's instructions must be observed.

No sharp, pointed or abrasive objects should be used for cleaning, as these always risk damaging the surfaces.

Reprocessing instructions (according to DIN EN ISO 17664:2004):

Classification according to KRINKO:	Non-critical
Spaulding Classification USA:	Non-critical
Place of use:	No special requirements.
Storage and transport:	No special requirements. It is recommended to reprocess an instrument as soon as possible after its use.
Cleaning preparations:	No special requirements. The device does not need to be disassembled for cleaning.

Cleaning:	<p>Coarse dirt or contamination must be removed after use with a disposable cloth / paper towel to prevent residue drying on the device.</p> <p>Dirt that has dried on the device makes the cleaning more difficult and may also cause material changes.</p>
Disinfection:	<p>According to the recommendation of the Commission for Hospital Hygiene and Infectious Disease Prevention (KRINKO) at the Robert Koch Institute (RKI) and of the Federal Institute for Drugs and Medical Devices (BfArM):</p> <p>Only disinfectants with proven antibacterial, antiviral and fungicidal effectiveness are to be used.</p> <p>German Federal Health Gazette 2012 - 55:1244-1310; DOI 10.1007/s00103-012-1548-6</p> <p>According to the FDA recommendation in the USA, the disinfection must be high-level disinfection.</p> <p>USA: CidexOPA Solution EU: Mikrobac Tissues, Incides Tissues</p> <p>In principle, alcohol-based (2-propanol) instrument disinfectants with a CE mark are suitable. We have tested and validated the material compatibility with Mikrobac Tissues (BODE Chemie GmbH), Incides N from (ECOLAB GmbH) and CidexOPA Solution.</p> <p>Please follow the instructions for use of the disinfectant you use.</p> <p>Procedure for disinfecting the light head:</p> <p>Carefully wipe the surfaces, applying the disinfectant to all parts:</p> <ol style="list-style-type: none"> 1. Wipe the surface of the light housing and the control elements in a cross-wise direction with a disinfectant wipe. 2. Wipe the indentations of the control elements lengthways with a disinfectant wipe. 3. Wipe the hand grip area in a cross-wise direction with a disinfectant wipe. 4. Wipe the protective glass in circular movements with a disinfectant wipe. <p>Procedure for disinfecting the light foot:</p> <p>Carefully wipe the surfaces, applying the disinfectant to all parts.</p> <p>Removing chemical residues:</p> <p>At the end of the application time, wipe over with a disposable cloth.</p>
Drying:	Allow the housing surfaces to dry fully at room temperature.
Maintenance:	Not applicable
Checks and functional testing:	<p>The examination light must be checked for impurities (e.g. encrustations, coatings) and damage and if necessary cleaned and disinfected again, or, in case of damage, taken out of service.</p> <p>The function check involves switching on the light and checking that the control elements are working.</p>
Packaging:	Not applicable
Sterilisation:	Not permitted!
Storage:	Not applicable
Additional information:	Not applicable

Maintenance


- The HEINE EL10 LED examination light does not require any maintenance. The LED lamp cannot be replaced. In case of a fault or defect, please contact HEINE.

Replacing fuses

- Switch off the examination light and remove the power cable from the device. The fuses are located in the fuse holder of the panel mount connector. Use only the fuses specified in the technical data as replacement fuses.

Changing the protective glass


- Once the protective glass has been replaced, check it regularly to ensure it is correctly and securely fitted.

First, remove the cover using a screwdriver on the lateral projections. Unscrew the holder using the tool provided. The protective glass can now be replaced. Fit the new protective glass in the reverse sequence of steps. Tighten the threaded ring with moderate force.

Technical safety inspection

- The power supply unit and power cable should be inspected for damage and that they are in good working order at regular intervals.
- The technical safety inspection must be carried out every two years in accordance with the generally accepted rules of technology.

- For the electrical inspection, the installation cover must first be removed.
- The testing point is located on the left hand side and is labelled accordingly.

Service

- The device may not be modified in any way.
- For repair, send the device to HEINE or an authorized service partner only, to avoid the risk of electrical shocks, burns or injuries due to improper repair work or handling.

This product does not require regular service.

Disposal

- The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

Technical data

Mains power:	100–240 V AC, 50/60 Hz
Power consumption:	33 VA
Fuse sets:	2 x T 2A H 250 V
Life time of the lamp:	approx. 30,000 hours
Minimum distance from light source to patient:	20 cm
Working distance from light source to patient:	30 cm
Central illuminance E_c :	$\geq 45,000 \text{ lx}^*$
Light field diameter d_{10} (at a 30 cm working distance):	approx. 8.5 – 14 cm
Overall irradiance:	Typically 115 W/m^2
Colour temperature:	$\geq 3000 \text{ K}$
Colour rendering index R_a :	≥ 85
Special colour index R_g :	Typically 58
<i>Environmental conditions for operation:</i> Temperature: Relative humidity: Air pressure:	$+10 \text{ }^\circ\text{C}$ to $+35 \text{ }^\circ\text{C}$ 10 % to 75 % 700 hPa to 1060 hPa
<i>Environmental conditions for storage and transport:</i> Temperature: Relative humidity: Air pressure:	$-20 \text{ }^\circ\text{C}$ to $+50 \text{ }^\circ\text{C}$ 10 % to 95 % 500 hPa to 1060 hPa
Protection class:	I
Operation altitude:	3000 m
Protection provided by enclosures (IP code):	IP20
Risk group (according to DIN EN 62471)	1 (low risk)
Weight incl. wall mount:	approx. 2.6 kg
Weight of clamp mount:	approx. 0.2 kg
Length of flexible arm:	97 cm

* smallest light field and at a 30 cm working distance, depending on t

General instructions

 Condensation can form when moving from a cold to a warm environment. Wait until the unit has warmed up and is absolutely dry before using it. Proper functioning is only guaranteed if original spare parts and accessories from HEINE are used.

The expected life cycle amounts with designated use and the observation of warning and safety information as well as the maintenance instructions up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

Electromagnetic Compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

	This device is intended for use by medical professionals in the electromagnetic environment specified below. The user of this device should assure that it is used in such an environment
	The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission reduced electrical immunity of the medical equipment.
	The ME device may not be stacked directly near or used directly beside other devices. If the ME device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

Guidance and manufacturer's declaration – electromagnetic emissions

The EUT is intended for use in the electromagnetic environment specified below. The customer or the user of the EUT should assure that it is used in such environment.

Emission test	Compliance	Electromagnetic environment – Guidelines
RF emissions CISPR11	Group 1	EUT uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The EUT is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. Warning: This device is intended only for use by medical professionals. This is a device of class A CISPR 11 in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the MEG or restrict the connection to the site.
Harmonic Emissions IEC 61000-3-2	Class A	Symmetrical three-phase EUTs and other EUTs.
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Passed	

Guidance and manufacturer declaration - Electromagnetic immunity

The EUT is intended for use in the electromagnetic environment specified below.
The customer or the user of the EUT should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guideline
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	± 2 kV for mains cables ± 1 kV for input and output lines	The supply voltage quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U _T , (>95% dip in U _T) for 1/2 period 40% U _T , (60% dip in U _T) for 5 periods 70% U _T , (30% dip in U _T) for 25 periods <5% U _T , (>95% dip in U _T) for 5 seconds	< 5% U _T , (>95% dip in U _T) for 1/2 period 40% U _T , (60% dip in U _T) for 5 periods 70% U _T , (30% dip in U _T) for 25 periods <5% U _T , (>95% dip in U _T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ME device requires continued operation during power mains interruptions, it is recommended that the EUT be powered by a UPS (uninterruptible power supply) or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Comment: U_T is the a.c. supply voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The EUT is intended for use in the electromagnetic environment specified below. The customer or the user of the EUT should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	3 V eff	Portable and mobile RF communication equipment should be used no closer to any part of the EUT, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	3 V/m	<p>Recommended separation distance: $d = 3,5/3 \cdot \text{SQRT}(P/W)$ $d = 3,5/3 \cdot \text{SQRT}(P/W)$ 80 MHz to 800 MHz $d = 7/3 \cdot \text{SQRT}(P/W)$ 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a** Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.
 To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EUT is used exceeds the applicable RF compliance level above, the EUT should be observed to verify normal operation.
 If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the EUT.

- b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

Recommended separation distances for portable and mobile RF communication equipment and the EUT

The EUT is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the EUT can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the EUT as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$	80 MHz to 800 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$	800 MHz to 2,5 GHz $d = 7/3 \cdot \text{SQRT}(P)$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3





For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

	The CE mark indicates that the product complies with the European medical device directive 93/42/EEC.
	Catalogue- or reference number
	Serial number
	Manufacturer
	Date of manufacture
	Product bearing this symbol may not be disposed of together with general household waste, but instead requires separate disposal according to local provisions. (European Waste Electrical and Electronic Equipment Directive, WEEE)
	Temperature limits in °C for storage and transport
	Temperature limits in °F for storage and transport
	Humidity limitation for storage and transport
	Pressure limitation for storage and transport
	Fragile, handle with care!
	Keep dry!
	Follow instructions for use! (Background: blue, foreground: white.)
	Unique Device Identification

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