HEINE QUALITY
MADE IN GERMANY

∴ HEINE Hand-held Indirect Ophthalmoscope

HEINE Optotechnik GmbH & Co. KG
Dornierstr. 6 · 82205 Gilching · Germany
E-Mail: info@heine.com · www.heine.com

MED 3814 2020-07-15

V-200.00.052

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HEINE Hand-held Indirect Ophthalmoscope

Please read and follow these instructions for use of and keep them for future reference.

Intended Use

The HEINE Hand-held Indirect Ophthalmoscope is intended to examine media (cornea, aqueous humor, lens, vitreous body) and the retina of the eye. The devices contain viewing optics and a battery-powered illumination unit. It should only be operated by qualified medical personnel within a professional healthcare facility.

For U.S. only:


⚠ Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

Warnings and Safety Information

⚠ CAUTION! Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

🔧 NOTE! This symbol indicates valuable advice in terms of set up, operation or maintenance, as applicable. Notes are important, but not related to hazardous situations.

Product overview

Designations:

- | | |
|-------------------|----------------------|
| 1 Loupe | 8 Extension piece |
| 2 Prism | 9 Holder |
| 3 Orbital support | 10 Friction clutch |
| 4 Marking | 11 Scale |
| 5 Sleeve | 12 Adjustment screw |
| 6 Connector | 13 Viewing apertures |
| 7 Adjustment pin | 14 HEINE handle |

Setting up

To set up the instrument, plug the instrument head on the HEINE handle (14).

Operation

Switch the illumination on, push the sleeve (5) down as far as it will go and illuminate a white surface from a distance of approx. 50 cm. The size of the illuminated area can be continuously adjusted by turning the sleeve (5). To switch on the red-free filter, push the sleeve (5) up to the green marking (4).

Monocular examination

Attaching the loupe

The loupe makes examination easier for presbyopes and gives a magnified fundus image: Slide the holder (1) onto the cylindrical upper part of the instrument. The loupe should be fitted behind and parallel to the vertical face of the prism. The beam should not traverse the loupe.

Using the orbital support

The orbital support (3) helps the user to steady the instrument during the examination. Spectacle wearers can use the extension piece which is delivered with the instrument. The extension piece (8) should be inserted sideways into the orbital support so that the guides (9) engage in the slit provided.

Binocular examination

Preparation

Release the friction clutch (10) on the binocular attachment. Remove the loupe from the instrument and attach the binocular attachment as illustrated. Tighten the friction clutch by rotating the lever provided until the attachment sits firmly on the instrument.

Adjusting the optics

Turn the knob (12) to adjust the P.D. to the appropriate value shown on the scale (11). With the instrument switched on, select the desired spot size. Look through the binocular attachment without accommodating at a light-coloured surface at a distance of 30–40 cm. A small object, such as a pencil, must be in focus. (+2 diopter lenses are built into the viewing apertures (13). Spectacle wearers should view through the "long-range" part of bi-focal spectacles). The fine adjustment of the P.D. can then be made by centring the light-spot in each viewing aperture by closing one eye at a time. Two replacement plain lenses are included with the instrument. If required the eye-pieces can be unscrewed and 2D lenses exchanged for the above.

The setup and operation of the HEINE handles are described in a separate instruction for use.

Hygienic reprocessing

The instruction is available:

- at www.heine.com
- in a paper version which you can request from the address listed

Maintenance

The instruments do not require maintenance or service.

Changing the light source

 Allow the device to cool down before changing the bulb.

- Remove the instrument from the handle (14) and pull the bulb out of the socket (6).
- Clean the glass of the new bulb with a soft cloth. Replace the bulb in the socket so that the guide (7) engages in the cutout of the socket.

Service

The device has no components serviceable by the end-user.

General Notes

 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

If you don't use the device for a longer period of time, please remove the batteries in advance.

The expected life cycle, when the device is normal used and the warning and safety information as well as the maintenance instructions are observed, is up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

General Warnings

 Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not use the device outside.

Light exposure hazard

 Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures.

The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

While no acute optical radiation hazards have been identified for direct or indirect ophthalmoscopes, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum intensity should not exceed 12 minutes (with the use of HEINE A.R. 16D/Ø54mm A.R. lens).

Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	Inside hospitals except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Performance features of the ME system that have been determined to be essential to the performance	None
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Hand-held indirect ophthalmoscope, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	See attachment
Test	Compliance
RF emissions CISPR11	Group 1 Class B
Harmonic Emissions	Passed*
Voltage Fluctuations/Flicker	Passed*

*n/a: "Not applicable" in the internally powered mode

Technical specification

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Nominal voltage	3.0 V – 3.7 V
Nominal current	Typ. 440 – 760 mA
Protection class	internal power supply
Device classification according ISO 15004-2	Group: 2
XHL Xenon spare bulb	X-002.88.050 (3,5 V)

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