HEINE QUALITY
MADE IN GERMANY

∴ HEINE SIGMA® 250 (M2)

CE

 HEINE Optotechnik GmbH & Co. KG
Dornierstr. 6 · 82205 Gilching · Germany
E-Mail: info@heine.com · www.heine.com
MED 113379 2020-08-18

V-200.00.066

HEINE SIGMA 250 (M2)

DEUTSCH	3
ENGLISH	7
FRANÇAIS	11
ESPAÑOL	15
ITALIANO	19
SVENSKA	23
NEDERLANDS	27
DANSK	31
SUOMI	35
PORTUGUÊS	39

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Please read and follow these instructions for use of and keep them for future reference.

Intended Use

The indirect Ophthalmoscope HEINE SIGMA 250 (M2) is a battery powered indirect ophthalmoscope for medical professionals, containing illumination and viewing optics intended to examine the media (cornea, aqueous, lens, vitreous) and the retina of the eye. It should only be operated by qualified medical personnel within a professional healthcare facility.

For U.S. only:

Federal law restricts this device to sale by or on the order of a Physician or Practitioner!

Warnings and Safety Information

CAUTION! Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

NOTE! This symbol indicates valuable advice in terms of set up, operation or maintenance, as applicable. Notes are important, but not related to hazardous situations.

Product overview


Name of positions

HEINE SIGMA 250 (M2):

- 1 Pressure clasp
- 2 Release
- 3 Aperture lever
- 4 Filter lever
- 5 Filter blue/yellow, diffuser
- 6 Stereopsis control
- 7 Parallax control
- 8 Beam control
- 9 Teaching mirror accessory
- 10 Dust cover
- 11 PD
- 12 Ocular
- 13 2D-Lens

S-FRAME:

- 14 Swivel mount
- 15 Temple / cord retainer
- 16 Nose rest

Setting up

To set up the HEINE SIGMA 250 (M2), couple the connection cable with the HEINE mPack mini. The setup and operation of the HEINE mPack mini are described in a separate instruction document. The spectacle frame has a vertical swivel mount (14), onto which the HEINE SIGMA 250 (M2) can be fitted from above by means of the pressure clasp (1). To remove the HEINE SIGMA 250 (M2), press the release (2) and pull them up and out of the mounting slot. Attach the bulb cord to the frame temple by means of the cord retainers enclosed (15). The retaining cord can be simply pushed on to the ends of the temples. We recommend that spectacle-wearers have purchase the clip-in correction frame and have it fitted with their prescription.

Operation

Light spot apertures

With lever (3) you can select two different-sized spots.

Filters

With the filter lever (4) you can add a red-free filter to any light spot. The optional blue or yellow filter or diffuser accessory (5) can be push-fitted to the front of the instrument housing.

Stereopsis and Parallax

The stereopsis control (6) below the instrument body adjusts the observation beam to suit the size of the patient's pupil (forward position close to patient: large pupil). The parallax control (7) adjusts the illumination beam to suit the size of the patient's pupil (upper position: large pupil). With both controls on the small setting it is possible to examine in stereopsis with 2 mm pupils. We recommend large stereopsis and small parallax settings for the examination of the periphery.

Lighting direction

The rotary beam control (8) can be used to swivel the beam vertically to improve reflex reduction and retinal illumination.

Teaching mirror

The optional teaching mirror accessory (9) can be push-fitted from below to the front of the housing.

Basic setting

Remove the dust cover (10).

Connect the power cord to your power source and switch on the illumination.

Carry out the following procedure:

Aperture lever (3) – large spot

Filter lever (4) – white light

Stereopsis control (6) – large pupil

Parallax control (7) – large pupil

Beam control (8) – middle of the range of adjustment

PD (11) – select your PD or set to an average position

Adjusting the optics

Put on the spectacle and check to ensure that the eyepieces are at the correct height. Remove the instrument and bend the nose rest (16) carefully to adjust the height. Bend the temple (15) and if necessary attach the retaining cord to the frame and adjust so that the instrument fits comfortably.

Setting the optics

Move the HEINE SIGMA 250 (M2) as close as possible to your eyes and look at the light spot at a distance of 30 – 40 cm. A small object like a pencil must be clearly-focused. Adjust the HEINE SIGMA 250 (M2) so that the light spot appears centred vertically in your view.

Adjust the PD setting by viewing the light spot alternately with the left and right eye and sliding the ocular (12) so that the spot is central. If you look through the ocular with both eyes, there should be no diplopic images.

Correct adjustment of the optics is particularly important when examining through small pupils. Once set, the adjustments need only be changed if another examiner uses the instrument.

Each instrument is supplied with separate plano lenses which can be used to replace the + 2D lenses if preferred. Change the ocular lenses only in a clean environment in order to avoid dust within the device.

Setting the brightness

Use the control dial on the HEINE mPack mini to adjust the brightness of the HEINE SIGMA 250 (M2):

Position 1: maximum brightness

Position 0: OFF

Hygienic reprocessing

The instruction is available:

- at www.heine.com

- in a paper version which you can request from the address listed

Maintenance

The device is maintenance free.

Changing the light source

With the HEINE SIGMA 250 (M2) the LED cannot be changed.

Service

The device has no components serviceable by the end-user.

General Notes

 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

The expected life cycle, when the device is normal used and the warning and safety information as well as the maintenance instructions are observed, is up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

General Notes and Warnings

 Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anesthetic environments).

Do not modify the device.

Use only original HEINE spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not use the device outside.

Light exposure hazard

 Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures. This device should be used with filters that eliminate UV radiation (<400 nm) and, whenever possible, filters that eliminate short-wavelength blue light (<420 nm).

The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

While no acute optical radiation hazards have been identified for direct or indirect ophthalmoscopes, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.


Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum intensity should not exceed 150 minutes (with the use of HEINE A.R. 16D/Ø54mm A.R. lens).

Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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Technical specification

Nominal voltage	3 V DC
Nominal current	max. 190 mA
Environmental conditions for operation	+10 °C to +35 °C 30 % to 90 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	-10 °C to +55 °C 10 % to 95 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-40 °C to +70 °C 10 % to 95 % rel. humidity 500 hPa to 1060 hPa
Device classification according to ISO 15004-2	Group 2
Protection Class	Internally powered

List of HEINE SIGMA 250 (M2) accessories

Teaching Mirror	C-000.33.036
Blue filter	C-000.33.313
Yellow filter	C-000.33.314
Diffusor	C-000.33.315

Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	Inside hospitals except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Performance features of the ME system that have been determined to be essential to the performance	None
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the SIGMA 250 (M2), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	See attachment
Test	Compliance
RF emissions CISPR11	Group 1 Class B
Harmonic Emissions	Passed*
Voltage Fluctuations/Flicker	Passed*

*n/a: „Not applicable“ in the internally powered mode

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