HEINE SIGMA 150® HEINE SIGMA 150® M2



Indirektes binokulares Brillenophthalmoskop

CE 93/42/EWG/CEE/EEC

HEINE

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. + 49 (0) 81 52 / 38 - 0 · Fax + 49 (0) 81 52 / 38 - 202
E-mail: info@heine.com · www.heine.com
med 6296 2 / 1.09

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DEUTSCH

HEINE SIGMA 150® / HEINE SIGMA 150® M2 Indirektes binokulares Brillenophthalmoskop

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Before using the ophthalmoscope read these instructions carefully and keep them in a safe place for future reference.

In accordance with local regulations this product should be disposed of as an electronic device separately.

Safety information

Range of applications: The HEINE SIGMA 150® is designed solely for examination of the eye, particularly the retina.

Safety in use: Use only those devices as power supplies which are designed for medical applications. Use as little light as possible and always switch off the power supply after use.

Do not leave switched on for longer than 10 minutes. After such use leave the unit switched off for at least 10 minutes.

Spectacle frame

The spectacle frame has a vertical swivel mount (1), onto which the optics can be fitted from **above** by means of the special fixing tongue (2).

To remove the optics, press the release (3) and **pull them up** and out of the mount. Attach the bulb cord to the frame temple by means of the cord retainer enclosed.

The retaining cord can be simply pushed on to the ends of the side-pieces. We recommend that spectacle-wearers have their own correction fitted to the frame.

Controls for examination and illumination

Light spot apertures: With lever (6) you can select one of two different-sized spots.

Filters: With the filter lever (7) you can add a red-free filter to any light spot. The optional blue or yellow filter accessory (13) can be push-fitted to the front of the instrument housing.

Stereopsis and Parallax: The **stereopsis control** below the instrument body adjusts the **observation beam** to suit the size of the patient's pupil (forward position close to patient: large pupil). The **parallax control** (8) adjusts the **illumination beam** to suit the size of the patient's pupil (upper position: large pupil).

Please note: With both controls on the small setting it is possible to examine in stereopsis with 2 mm pupils. We recommend **large stereopsis** and **small parallax** settings for the examination of the **periphery**.

The rotary beam control (5) can be used to swivel the beam vertically to improve reflex reduction and retinal illumination.

Teaching mirror: The optional teaching mirror accessory (14) can be push-fitted from below to the front of the housing.

Setting-up

- Remove the dust cover (15)
- Connect the power cord to your power source and switch on the illumination.
- Carry out the following procedure:
Aperture lever (6) – large spot
Filter lever (7) – white light
Stereopsis control – large pupil
Parallax control (8) – large pupil
Beam control (5) – middle of the range of adjustment
PD (4) – select your PD or set to an average position

Adjusting the optics

- Put on the spectacle ophthalmoscope and check to ensure that the eyepieces are at the correct height. Remove the instrument and bend the nose rest (16) carefully to adjust the height.
- Attach the retaining cord to the frame so that the instrument fits comfortably.
- Observe the light spot on a level surface at a distance of 30–40 cm. Adjust the PD setting for the left and right eye by closing first one then the other and centralising the light spot image in your field of view.

Important: The eyepieces contain 2D lenses. This enables normal-sighted, emmetropic patients to view the retinal image without effort. If required, these can be replaced by plano lenses supplied with the instrument. To do this, unscrew the eyepieces (4).

Changing the bulb

Important Note: This instrument will only function properly if an original HEINE replacement bulb is used.

- Switch off and allow the instrument to cool down.
- Push the bulb compartment cover (9) up and out.
- Shake the used bulb out of its housing (11).
- Wipe the new bulb clean with a soft cloth and insert it into the compartment (11). The pin (12) must coincide with the groove (10).
- Replace the cover (9).

Care of the instrument

Clean glass surfaces with clean, soft cloth using alcohol or spectacle glass cleaner, if necessary. Read the instructions for cleaning the HEINE S-Frame®.

Technical specification

Weight of the optics: 90 g
Halogen bulb: 6 V/5 W
Color temperature: ca. 3000 K
Working life: > 30 hours

Accessories and spare parts

S-Frame	C-000.33.036
Teaching mirror	C-000.33.302
Blue filter	C-000.33.313
Yellow filter	C-000.33.314
HEINE XHL® Xenon halogen spare bulb	X-004.88.093

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Warranty

Instead of the usual 2-year guarantee, we guarantee this device for 5 years from the date of delivery from our factory (this excludes consumables such as bulbs, tips and batteries). We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product during the guarantee period, the buyer must prove that any such defect was present at the time the product was delivered. The usual legal guarantee and warranty do not apply to incorrect use, the use of non-original HEINE parts or accessories. (Especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance.) Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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