# HEINE EN 100® Ausbaueinheit

CE 93/42/EWG/CEE/EEC

## HEINE

HEINE Optotechnik GmbH & Co. KG  
Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49(0)81 52/38-0 · Fax +49(0)81 52/38-202  
E-mail: info@heine.com · www.heine.com  
med 5001 1/2.10

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### HEINE EN 100® Ausbaueinheit

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
## Instructions ENGLISH

### HEINE EN 100® Extension Unit

 Before operating the HEINE EN 100® Extension Unit read these instructions carefully and keep them in a safe place for future reference.

 Applied Part BF  In accordance with local regulations this product should be disposed of as an electronic device separately.

 AC powered

 Fuse  WITH RESPECT TO ELECTRICAL SHOCK, FIRE, AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA C22.2 No. 601.01 28KY

### Contents

HEINE EN 100® Extension Unit  
Mounting plate  
3 Fixing plugs U6  
3 Screws 4x35  
Instrument locks with instructions S-277.04.06

### Safety Information

**Range of applications:** The HEINE EN 100® Extension Unit is designed for use exclusively as a power source for medical diagnostic instruments with 3,5 V bulbs and should only be operated by medically-trained staff.

**General Safety:** This unit is maintenance-free. Do not open the housing. Repairs should only be carried out by authorised personnel.

**Safety in use:** Do not operate this unit in the presence of an explosion hazard or in a strong magnetic field.

We can only guarantee the correct performance of this product when it is used with the HEINE EN 100® and HEINE EN 100-1®.  
It is possible to use a maximum of two handles for continuous use.

### Assembly and Disassembly

The HEINE EN 100® Extension unit should only be fixed to the wall by means of the mounting plate supplied. Please refer to the enclosed assembly instructions.

- Disconnect the mains cord of the EN 100/EN 100-1 and detach the unit from the mounting plate.
- Hold the mounting plate for the extension unit at the same height and right next to the EN 100/EN 100-1 mounting plate and mark the drill-holes.
- Use a 6 mm drill to drill holes at least 35 mm deep and insert the fixing plugs.
- Screw the mounting plate in place with a suitable screwdriver.
- Remove the blanking plate on the right-hand side of EN 100/EN 100-1 housing and re-mount the EN 100/EN 100-1 on its mounting plate, making sure that all three pins engage when pressed in and down.
- Check all three pins again.
- Connect the wires from the EN 100 Extension unit to the EN 100/EN 100-1 and push them together, at the same time inserting the wires into the free space in the housing. Locate the unit with all three pins and press in and down to attach it to the plate.

To disassemble, disconnect the mains cord and remove the unit from the mounting plate by pressing upwards.

**Please note:** The new version of the HEINE EN 100® extension unit allows you to mount several units sequentially. In this case, handles can also be used both simultaneously or independently of each other.

### Cleaning

Disconnect the mains cord!  
Clean and disinfect using 70% Isopropyl alcohol.

### Specification

Output voltage:	3,5 V (can be selected by the user)
Output current:	max. 800 mA
Operation:	Continuous with intermittent loading
Classification:	Type BF
Dimensions:	110 x 205 x 85 mm
Weight:	440 g
Ambient temperature:	+10 °C to +40 °C
Rel. Humidity:	10 to 75 %
Air pressure:	700 to 1060 hPa
Protection class:	IP X0

**Applicable standards:** EN 60601-1  
EN 60601-1-2  
UL 60601-1  
CAN/CSA C22.2 No. 601.01

### Warranty

Instead of the usual 2-year guarantee, we guarantee this device for 5 years from the date of delivery from our factory (this excludes consumables such as bulbs, tips and batteries). We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product during the guarantee period, the buyer must prove that any such defect was present at the time the product was delivered. The usual legal guarantee and warranty do not apply to incorrect use, the use of non-original HEINE parts or accessories. (Especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance.) Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.  
For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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 CA

 Fusible  WITH RESPECT TO ELECTRICAL SHOCK, FIRE, AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA C22.2 No. 601.01 28KY

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