## Mode d'emploi

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#### Warranty

Instead of the usual 2-year guarantee, we guarantee this device for 5 years from the date of delivery from our factory (this excludes consumables such as bulbs, tips and batteries). We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product during the guarantee period, the buyer must prove that any such defect was present at the time the product was delivered. The usual legal guarantee and warranty do not apply to incorrect use, the use of non-original HEINE parts or accessories. (Especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance.) Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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### HEINE Retinoscope Trainer

#### General Information

The HEINE Retinoscope Trainer is not intended for use on patients and is not a medical product according to directive 93/42 EWG/CEE/EEC.

The instrument is intended as a training device for the technique of retinoscopy. It can also be used as a reference device to check the function of instruments used for refraction measurement, provided the relevant optical principles are applied. It is not designed for ophthalmoscopy, instead of this, we recommend the Ophthalmoscopy Trainer C-000.33.010.

We are unable to give an introduction to Retinoscopy and its application in refraction in these instructions. We recommend the published literature on this subject.

#### Technical Description

The optics of the trainer include an achromatic lens with a pupil aperture and an orange-coloured curved shell which can be moved along the optical axis.

##### Technical data:

Weight	600g
Achromat	f = 32mm
Range of refraction	-7 bis + 6 dpt
Accuracy	+/-0,1 dpt at 0dpt*)
Pupil size:	2 to 8mm

\*)See below: Retinoscopy.

#### Using the instrument

**Setting up the trainer:** Adjust the body of the trainer so that the optical axis corresponds to your viewing axis, (angle can be adjusted by loosening screw1).

**Selecting refraction:** Rotate the housing to select a value as shown on the scale(2). Select the desired value carefully, particularly around the 0 dpt position. View the scale from directly above to reduce the chance of error.

**Selecting Pupil Size:** Six different pupil sizes can be selected by rotating the selector wheel (3). Ensure that the desired value clicks into place. The selected pupil size is indicated (4).

**Using trial lenses:** A slot is provided in front of the lens to accommodate two trial lenses. The scale (6) gives an indication of the cyl. Angle

**Retinoscopy:** Please ensure that you view along the optical axis of the trainer to avoid getting a false result. You can check this by making sure that the two reflexes generated by the retinoscope on the front and back side of the trainer lens co-incide. In general, the higher the refraction, the less precise is the result.

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