HEINE VIDEO OMEGA 2 C®

– mit Kamera –

CE 93/42/EWG/CEE/EEC

HEINE

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49(0)8152/38-0 · Fax +49(0)8152/38-202
E-mail: info@heine.com · www.heine.com
med 7596 1/2.08

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Warranty

We guarantee this device (except bulbs) for five years from date of purchase. We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product, the buyer must prove that any such defect was present at the time the product was delivered. This guarantee excludes any defects which occur because of wear and tear, incorrect use, the use of non-original HEINE parts, (especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance). Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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Instructions for use

HEINE VIDEO OMEGA 2 C® – with Camera –

Read these instructions carefully before using the HEINE VIDEO OMEGA 2 C® and keep them in a safe place for future reference.

Safety information

Range of applications: The HEINE VIDEO OMEGA 2 C® is used exclusively for eye examination, particularly retinal examination. A video image of the eye is shown on a monitor.

Safety in use: Use only power sources which are designed for use in medical applications. This also applies to any accessories and the monitor used with the system.

Video installation

A complete video installation consists of the following components:

1. HEINE VIDEO OMEGA 2 C® Ophthalmoscope, 6 V with C-mount
2. CCD camera with C-mount and 3 m cord
3. Power Source HEINE EN 50® or mPack
4. Multi-coated anti-reflex ophthalmoscope lens HEINE A.R.20
5. Two scleral depressors (large/small)
6. Spare 6 V bulb for ophthalmoscope
7. Desk-top video processor incl. monitor cord and BNC/cinch adaptor for colour monitor with BNC or cinch connector
8. Carrying case for items 1. to 6.
9. Colour monitor and video recorder/processor (optional).

Important note: Colour monitor and video recorder/processor are not included in the specification as delivered.

Setting-up

1. Instructions:

First, please read and follow the special instruction leaflets for **HEINE OMEGA, HEINE EN 50® and Video camera** which are supplied with this unit. Compared with the HEINE OMEGA 100, the **HEINE VIDEO OMEGA 2C** differs in the following ways: There is no control to select different-sized light spots or to select filters or to vary convergence and parallax.

2. Connect the various components. Please take care not to kink the connecting cords as this can cause damage.

3. Switch on the HEINE EN 50®, camera processor and monitor.

4. Adjust and set up the ophthalmoscope, lens and HEINE EN 50® just as you would with the standard OMEGA (see special instruction leaflets).

5. Focus the monitor image by adjusting the lever (2) on the left of, and below the camera.

No further adjustment is necessary.

Important note: The monitor image is reversed left to right compared with the image seen by the examiner.

Adjusting the camera

The camera is adjusted before delivery for an average examination distance (a) of 350 mm. The focus lever (2) adjust the focal point by ca. ± 100 mm. If your usual examination distance varies significantly and the monitor image is out of focus, or if the camera needs to be adjusted after being fitted to another instrument, proceed as follows:

1. **Connect up the individual components** and switch on the HEINE EN 50®, camera processor and monitor.

2. Check the adjustment

• Rest the HEINE VIDEO OMEGA 2 C® on a table and place a vertical flat object e.g. a book at the desired examination distance (a).

Please note: When determining the desired examination distance, please bear in mind that the camera must be adjusted to the virtual image produced by the ophthalmoscope lens.

As a general rule, the examination distance equals the distance from the examiner's eyes to the ophthalmoscope lens minus 100 mm.

• Adjust the lever (2) to the mid-position.

• Check the focus without regard to whether the monitor image is vertical or not. If the image is focused, go on to point 4. If not, move the object until focused and go on to point 3.

3. Adjust the examination distance

• Switch off the camera processor and unscrew the cord from the camera head.

• Loosen the camera in the C-mount (3) by unscrewing it.

• Unscrew the focus lock-ring (1) so that the camera can be adjusted.

• If the object is too close to the camera, screw the camera into the C-mount (3) by about 3 turns. (If not close enough, the reverse applies).

• Screw the cord back in and switch on the camera processor.

• Repeat point 2. and 3. and vary amount of turns applied.

4. Adjust the position of the image

• Rotate the camera by not more than half a turn until an upright image reversed left to right is obtained.

• Lock the camera in this position with the focus lock ring (1).

• The **vertical position** of the image can be set by the screw (4) located on the upper part of the optics unit (VERT.ADJUST).

• In rare cases it may be necessary to adjust the **horizontal position** of the image. This can be set by the screw (5) located at the front of the optics unit (HORIZ.ADJUST).

Using the camera with instruments other than OMEGA

The C-mount (3) makes it possible to fit the camera to other similarly-equipped instruments. Switch off the camera processor first, remove the camera cord and unscrew the C-mount adaptor including the camera.

With other instruments it may be necessary to re-adjust the camera focus. To do this, please follow the instrument manufacturer's instructions.

Technical Specification

Ophthalmoscope

Modified HEINE OMEGA 100 with 6 V halogen bulb

Power Supply

HEINE EN 50® with continuous brightness control, 0 to 6 V/10 W maximum.

Lens

Ophthalmoscopy lens HEINE A.R.20, aspheric with special multi-coating.

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