# Projektor HEINE HK 7000®

## HEINE

**HEINE Optotechnik GmbH & Co. KG**  
Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49(0)8152/38-0 · Fax +49(0)8152/38-202  
E-mail: info@heine.com · **www.heine.com**  
**med 4001**    1/10.08

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### Warranty

Instead of the usual 2-year guarantee, we guarantee this device for 5 years from the date of delivery from our factory (this excludes consumables such as bulbs, tips and batteries). We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product during the guarantee period, the buyer must prove that any such defect was present at the time the product was delivered. The usual legal guarantee and warranty do not apply to incorrect use, the use of non-original HEINE parts or accessories. (Especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance.) Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded. For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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**EN 60 601-1, EN 60 601-1-2**

**IEC 60601-1 second edition, UL 60601-1 first edition**

## Instructions

	Before using the HK 7000 projector read these instructions carefully and keep them in a safe place for future reference.
--	--

	Hot surface		Fuse
	Applied Part BF		In accordance with local regulations this product should be disposed of as an electronic device separately.
	AC powered		

#### Contents:

Projector
Mains cord
Spare bulb
2 spare fuses
Instructions

#### Safety Information

**Range of applications:** The fiber optics projectors HK 7000 is designed exclusively for use as a light source with a fiber optics cable for medical diagnostic instruments and medical examination lights. The unit is designed for continuous use.

Never use the projector for transillumination because of risk of skin burning.

**General safety:** This unit needs no maintenance. Do not remove the housing. With the exception of a bulb change, repairs may only be carried out by authorized persons.

**Safe installation:** Set up the projector out of reach of non-authorized personnel horizontally so that free circulation of air on all sides is possible. Do not place on soft furniture, a cloth or carpet. Do not put inflammable objects in the light exit (1) or the bulb compartment (6). Do not look into the light exit (1) and do not open the lid (5) when the projector is switched on.

**Safety in use:** The mains voltage on the identification plate must correspond to the local mains voltage. Do not use this unit in the presence of an explosion hazard or in a strong magnetic field.

#### Assembly instructions

Instructions reference S 209.50.02 are provided with the appropriate assembly kit for wall – or wheeled-stand mounting.

#### Instructions for installation and use

- Connect the unit to the mains with the mains cord.
- Plug in the fiber optics cable to the light exit (1) until the cable plug engages with a click.
- Switch on at (3), the green indicator will light up.
- Set the desired brightness with the control (2).
- The switch (4) can be used to select the normal maximum brightness (for long bulb life) or increased maximum brightness (reduced bulb life).
- If possible, switch off after use to prolong bulb life.

HEINE fiber optics cables are equipped with an automatic connector: Hold the handle in one hand and use the other to insert the 3 male prongs of the connector into the guides of the handle-head. Press the instrument into the handle until it clicks into place.

To remove the instrument, hold the handle in one hand and rotate the knurled upper ring anti-clockwise. With the other hand remove the instrument by pulling gently.

#### Changing the bulb

**Please note: The performance of this unit can only be guaranteed if genuine HEINE bulbs are used.**

- Remove the mains cord plug.
- Hot bulb allow the unit to cool down.
- Open the lid (5) and swing the bulbholder upwards and out of the bulb compartment (6).
- Pull the bulb out of the bulbholder.
- Insert the new bulb into the bulbholder, taking care not to touch the glass with bare hands. Push the bulb in as far as it will go and avoid bending the contact pins.
- Close the lid (5). The bulbholder will return automatically to its correct position.
- Re-connect the mains cord.
- Switch off the projector when not in use to prolong bulbs life.

#### Changing a fuse

- Remove the mains cord.
- The fuses are located in the fuse compartment of the mains cord socket. Use only replacement fuses as specified in the technical specification.

#### Cleaning and Sterilization

**Projectors:**

Always remove the mains cord! Always clean the unit after it has cooled down. To clean, use only a soft cloth and soapy water or similar non-aggressive solutions.

**HEINE Fiber optics cables:**

Gas sterilization or wiping clean with alcohol or soapy water are recommended. Do not autoclave or soak in any solutions.

The optical surfaces at both ends of the cable can be cleaned by a cotton wool swab soaked in alcohol.

<b>Technical data</b>	<b>Fuse</b>
Mains voltage (see nameplate)	230/240V/1.6A/50-60 Hz 2 x 1.6 A slow 100V/3.5A/50-60 Hz 2 x3.15 A slow 120V/2.5A/50-60 Hz 2 x2.5 A slow/ 250V AC 127V/2.5A/50-60 Hz 2 x2.5 A slow
Power consumption:	280VA
Bulb:	24V / 150 W
Classification:	Type BF
Dimensions:	217 x 177 x 125 mm
Weight:	3400 g
Ambient temperature:	+ 10° to +40 °C
Relative humidity:	10 to 75 <span> </span> %
Air pressure:	700 to 1060 hPa
Protection class:	IP 00

#### Accessories

XHL Spare halogen bulb 150 W	<b>Y-096.15.102</b>
Spare fuses	<b>Y-096.15.110</b>
Wall support	<b>Y-096.15.105</b>
Wheeled floor stand	<b>Y-096.50.002</b>
Standard fiber optics cable	<b>Y-003.99.518</b>

#### Standards applicable:

**EN 60 601-1, EN 60 601-1-2**

**IEC 60601-1 second edition, UL 60601-1 first edition**

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