HEINE BETA® Griffe

HEINE BETA® Handles

HEINE Optotechnik GmbH & Co. KG
Kientalsr. 7 · 82211 Herrsching · Germany
Tel. +49(0)81 52/38-0
Fax +49(0)81 52/38-202
E-Mail: info@heine.com · www.heine.com
med 1196 1/3.10

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DEUTSCH

HEINE BETA® Griffe

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ENGLISH

HEINE BETA® handles

Before using the HEINE BETA® handle read these instructions carefully and keep them in a safe place for future reference.

Safety Information Symbols

In accordance with local regulations this product should be disposed of as an electronic device separately.

Range of applications: BETA handles are designed exclusively for use with medical diagnostic instruments with bulb illumination.

General safety:

– Do not open the bottom insert or the transformer.
 – All repairs other than a change of batteries should be carried out by qualified persons.

Safety in use:

- Handles should not be used in rooms where an explosion hazard is present.
- The handles fitted with rechargeable batteries should not be used when being charged in the charger and should be switched off before recharging.
- The performance of this instrument can only be guaranteed if genuine HEINE rechargeable batteries are used.
- The voltage shown on the identification plate of the mains transformer or bottom insert of the BETA R 3,5V handle must be identical to the local mains voltage.

Application

Automatic Connector (AV):

– Attaching the instrument:

– Attach an instrument head by inserting it from above so that the connector lugs mate with the grooves in the handle. Press the head into the handle until it engages with an audible click.

– Removing the instrument:

– Remove the head by twisting the knurled ring with thumb and forefinger of the hand holding the handle. Pull off the head with the other hand.

Twist-Lock Connector (TL):

– Attaching the instrument:

– To connect an instrument head, hold the handle in one hand, match up its mounting lugs with the recesses in the head, press in and turn the handle clockwise until securely attached.

– Removing the instrument:

– To remove the head rotate it anti-clockwise and pull it off.

ON/OFF switch and brightness control: Press the button and rotate the black ring clockwise.

Using HEINE BETA® handles

BETA handles consist of three basic components:

- Handle module
 - Rechargeable 3,5 V battery or batteries
 - Bottom insert
- The handle module (sleeve with head) is the same for all handles. The battery and bottom insert determine the handle type. By changing the bottom insert and battery the handle can be converted from one model to another.

Changing the battery:

- Unscrew the bottom insert from the handle module.
- Remove the batteries or rechargeable battery.
- Insert the new battery.
- Screw the bottom insert into the handle module.

1. HEINE BETA® battery handle:

Use alkaline "C"-size batteries. With the BETA extension piece a third battery can be added to give 3,5 V operation (see accessories).

2. HEINE BETA® NT rechargeable handle:

You can fit either 3,5 V NiMH or Li-Ion rechargeable batteries.

3. HEINE BETA® cord handle:

The cord handle contains no battery. It may only be connected to 6 V max. supply.

4. HEINE BETA® R 3,5 V rechargeable handle:

The bottom insert incorporates a cap to protect the plug pins. This should be screwed on after charging. You can fit either 3,5V NiMH or Li-Ion rechargeable batteries.

Technical specification

Mains supply: 230 V/50–60 Hz
 See bottom insert: 120 V/60Hz
 Power consumption: max. 2 VA
 Safety class: Double insulated
 Recharging current: max. 75 mA

5. HEINE BETA® TR 3,5 V rechargeable handle:

Use only the charging transformer supplied with the handle. During charging the green LED lights up to indicate that a charging current is flowing. You can fit either 3,5 V NiMH or Li-Ion rechargeable batteries.

Technical specification

Mains supply: 230 V/50–60 Hz
 See bottom insert: 120 V/60Hz
 Power consumption: max. 5 VA
 Safety class: Double insulated
 Thermal fuse: in transformer
 Recharging current: max. 75 mA

3,5 V rechargeable battery instructions

Spontaneous discharge can occur when the battery is not recharged over a period of weeks. In this case a 14-hour charge is enough.

Replace the battery when a 14-hour charge is no longer sufficient to guarantee enough working capacity.

Maintenance

(Cleaning / Disinfection / Sterilising)

- The components can be wiped with a cloth moistened with water, cleaning solution or disinfectant.
- Ensure that the handle is dry after cleaning.
- Moisture must not be allowed to enter the rechargeable battery or the bottom insert.
- BETA R 3,5 V and BETA TR 3,5 V rechargeable handles should be disconnected from the mains supply before cleaning. BETA NT 3,5 V handles must be first removed from the charger.

Accessories and spare parts

Handle module	X-000.99.010
Bottom insert for battery handle	X-000.99.119
Extension piece for third battery	X-000.99.020
Bottom insert for cord handle	X-000.99.213
Bottom insert for R 3,5 V handle	X-002.99.377
Bottom insert for TR 3,5 V handle	X-002.99.385
Bottom insert for NT 3,5 V handle	X-002.99.412
3,5 V rechargeable battery (NiMH)	X-002.99.382
3,5 V rechargeable battery Li-Ion (Lithium Ion)	X-007.99.381
Main transformer pour poignée TR 3,5 V handle	X-000.99.308
Handle rest for TR 3,5 V handle	X-000.99.307

Applicable Standards:

IEC 60 601-1 (Electrical Medical Devices – Part 1: General conditions for safety including the relevant performance data.)

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