HEINE LAMBDA 100®

CE 93/42/EWG/CEE/EEC

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49 (0) 81 52 / 38 - 0
Fax +49 (0) 81 52 / 38 - 2 02
E-Mail: info@heine.com · www.heine.com
med 6896 1/10.11

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HEINE LAMBDA 100® Retinometer

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Instructions for use

HEINE LAMBDA 100® Retinometer

Read these instructions carefully before using the HEINE LAMBDA 100® retinometer and keep them in a safe place for future reference.

In accordance with local regulations this product should be disposed of as an electronic device separately.

Safety information

Range of applications: The HEINE LAMBDA 100® retinometer is only used to test the visual acuity of the eye.

Safety in use: Use only power sources which are designed for use in medical applications.

Please note: The performance of this instrument can only be guaranteed if genuine HEINE rechargeable batteries or alkali manganese batteries are used.

Terms used:

- 1 Brow rest
- 2 Light exit
- 3 Grid orientation selector
- 4 Visual acuity selector
- 5 Bulb

How LAMBDA 100 works

The HEINE LAMBDA 100® works on the principle of the Maxwellian view: A micro-aperture is illuminated by a halogen bulb through a red filter and imaged by an optical system into the patient pupil. The optical system consists basically of two lenses between which optical grids with variable spacing can be positioned in the parallel beam which passes through them. The resulting diffraction forms a circular test pattern with equally-spaced red and black lines on the retina. The distance between the lines corresponds to that of the Snellen E (Visus 1 = 33 lines / degree of visual angle). The orientation of the lines can be selected by means of a prism in 45° steps. Because the beam in the pupillary plane is very narrow (a few tenths of a mm), a tiny "window" in the opacity of the lens is enough to allow the light to pass through for a successful examination.

General Instructions for acuity-testing with the Retinometer

Testing potential acuity with retinometers has been the method of choice for some years especially where opacities are present. It is usually possible to get very reliable results provided a few basic preconditions are observed:

With very dense opacities below an optotype Visus of 0,1 the reliability of the estimate is reduced, the post-operative visus is often better than that obtained by the retinometer.

With amblyopic eyes, the tendency is more towards false positives because of the size of the test pattern.

With ametropia of more than 6 Dpt. the results can be improved by examining the patient while he is wearing his glasses. The same comments apply to high astigmatism.

Further factors which can have a negative influence on the results are nystagmus, tremor, communication problems, senility or debility of the patient.

Generally speaking, degenerative conditions can cause considerable problems, which can only be overcome if the examiner has the necessary experience of this technique.

How to use the instrument

The instrument is very easy to use: With the lever (3) select the grid orientation in 45° steps and then choose one of the six acuity values with the selector wheel (4). Use the brow rest (1) to support the retinometer on the patient's forehead so that the distance to the patient's cornea is correct and to hold the instrument steady. The light from the light exit can only be seen if the room is darkened. (Support the brow rest on a sheet of white paper. With the instrument switched on you will see three small spots of light, the distance between which changes when the acuity setting is changed). Please read the instructions for use supplied with your battery handle or power supply to see how to switch on and change the brightness of the light etc.

Examination

Always dim the room light before examining. Eyes with opaque media should always be dilated. This is not necessary where no opacity is present. Explain to the patient as clearly as possible what he will see with the help of the retinometer instruction card supplied with the set. Use the better eye first to get the patient used to the instrument.

- Make sure that the patient is in a comfortable position. He should look into the far corner of the darkened room.
- Switch on the retinometer. If preferred you can set the desired brightness by shining the light into your own eye. (It may be necessary to change the brightness during the examination).
- Use the lever (3) to choose the grid orientation and select a low visual acuity value with the selector wheel (4).
- Rest the retinometer upright against the patient's forehead and direct the red light beams onto the pupil. At the normal reading distance look alongside or over the top of the retinometer (there is a "window" built into the brow rest) and observe the points of light on the patient's cornea.
- Move the beam over the pupil by slightly turning or tilting the retinometer until the patient recognizes the grid with the red and black pattern. Agree with the patient how he will indicate the grid orientation (with a hand-movement or verbally). If he mentions circular spots or distorted lines, this could be an indication of a malfunction of the macula.
- Change the grid orientation by rotating the lever (3) with your free hand without moving the retinometer.
- Increase the visual acuity value until the patient is no longer able to recognize the direction of the grid.

The last value which he was able to recognize represents the potential acuity.

Changing the bulb

Please note: The performance of this unit can only be guaranteed if genuine HEINE bulb are used.

Remove the retinometer from the handle and pull the bulb (5) out and down. Push the new bulb in as far as the stop.

Cleaning and maintenance

Use a soft cloth with a little alcohol or cleaning fluid to clean the instrument body. Clean the light exit (2) with a cotton wool bud and alcohol.

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