**Instruction for use**  
**HEINE SIGMA® 250 (M2)**

# **Instruction for use**

## **HEINE SIGMA 250 (M2)**

<b>DEUTSCH</b>	<b>3</b>
<b>ENGLISH</b>	<b>7</b>
<b>FRANÇAIS</b>	<b>11</b>
<b>ESPAÑOL</b>	<b>15</b>
<b>ITALIANO</b>	<b>19</b>
<b>SVENSKA</b>	<b>23</b>
<b>NEDERLANDS</b>	<b>27</b>
<b>DANSK</b>	<b>31</b>
<b>SUOMI</b>	<b>35</b>
<b>PORTUGUÊS</b>	<b>39</b>

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Please read and follow these instructions for use of and keep them for future reference.

### Intended Use

The indirect Ophthalmoscope HEINE SIGMA 250 (M2) is a battery powered indirect ophthalmoscope for medical professionals, containing illumination and viewing optics intended to examine the media (cornea, aqueous, lens, vitreous) and the retina of the eye. The product may only be used by qualified medical professionals and in professional healthcare facilities.

#### For U.S. only:

**Federal law restricts this device to sale by or on the order of a Physician or Practitioner!**

### Warnings and Safety Information

**CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

**NOTE!** This symbol indicates valuable advice. This notes are important, but no related to hazardous situations.

### Product overview


#### Name of positions

#### HEINE SIGMA 250 (M2):

- 1 Pressure clasp
- 2 Release
- 3 Aperture lever
- 4 Filter lever
- 5 Filter blue/yellow, diffuser
- 6 Stereopsis control
- 7 Parallax control
- 8 Beam control
- 9 Teaching mirror accessory
- 10 Dust cover
- 11 PD
- 12 Ocular
- 13 2D-Lens

#### S-FRAME:

- 14 Swivel mount
- 15 Temple / cord retainer
- 16 Nose rest

## **Setting up**

To set up the HEINE SIGMA 250 (M2), couple the connection cable with the HEINE mPack mini. The setup and operation of the HEINE mPack mini are described in a separate instruction document. The spectacle frame has a vertical swivel mount (14), onto which the HEINE SIGMA 250 (M2) can be fitted from above by means of the pressure clasp (1). To remove the HEINE SIGMA 250 (M2), press the release (2) and pull them up and out of the mounting slot. Attach the bulb cord to the frame temple by means of the cord retainers enclosed (15). The retaining cord can be simply pushed on to the ends of the temples. We recommend that spectacle-wearers have purchase the clip-in correction frame and have it fitted with their prescription.

## **Operation**

### **Light spot apertures**

With lever (3) you can select two different-sized spots.

### **Filters**

With the filter lever (4) you can add a red-free filter to any light spot. The optional blue or yellow filter or diffuser accessory (5) can be push-fitted to the front of the instrument housing.

### **Stereopsis and Parallax**

The stereopsis control (6) below the instrument body adjusts the observation beam to suit the size of the patient's pupil (forward position close to patient: large pupil). The parallax control (7) adjusts the illumination beam to suit the size of the patient's pupil (upper position: large pupil). With both controls on the small setting it is possible to examine in stereopsis with 2 mm pupils. We recommend large stereopsis and small parallax settings for the examination of the periphery.

### **Lighting direction**

The rotary beam control (8) can be used to swivel the beam vertically to improve reflex reduction and retinal illumination.

### **Teaching mirror**

The optional teaching mirror (9) can be push-fitted from below to the front of the housing.

### **Basic setting**

Remove the dust cover (10).

Connect the power cord to your power source and switch on the illumination.

Carry out the following procedure:

Aperture lever (3) – large spot

Filter lever (4) – white light

Stereopsis control (6) – large pupil

Parallax control (7) – large pupil

Beam control (8) – middle of the range of adjustment

PD (11) – select your PD or set to an average position

### **Adjusting the optics**

Put on the spectacle and check to ensure that the eyepieces are at the correct height. Remove the instrument and bend the nose rest (16) carefully to adjust the height. Bend the temple (15) and if necessary attach the retaining cord to the frame and adjust so that the instrument fits comfortably.

## Setting the optics

Move the HEINE SIGMA 250 (M2) as close as possible to your eyes and look at the light spot at a distance of 30 – 40 cm. A small object like a pencil must be clearly-focused. Adjust the HEINE SIGMA 250 (M2) so that the light spot appears centred vertically in your view.

Adjust the PD setting by viewing the light spot alternately with the left and right eye and sliding the ocular (12) so that the spot is central. If you look through the ocular with both eyes, there should be no diplopic images.

Correct adjustment of the optics is particularly important when examining through small pupils. Once set, the adjustments need only be changed if another examiner uses the instrument.

Each instrument is supplied with separate plano lenses which can be used to replace the + 2D lenses if preferred. Change the ocular lenses only in a clean environment in order to avoid dust within the device.

## Setting the brightness

Use the control dial on the HEINE mPack mini to adjust the brightness of the HEINE SIGMA 250 (M2):

Position 1: maximum brightness

Position 0: OFF

## Hygienic reprocessing

The instruction is available:

- at [www.heine.com](http://www.heine.com)

- in a paper version which you can request from the address listed

## Maintenance

The device is maintenance free.

## Changing the light source

With the HEINE SIGMA 250 (M2) the LED cannot be changed.

## Service

The device has no components serviceable by the end-user.

## General Notes

The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

The expected life cycle, when the device is normal used and the warning and safety information as well as the maintenance instructions are observed, is up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

Note to the user and/or patient:

All serious incidents that occur in connection with the product must be reported to HEINE Optotechnik GmbH & Co. KG and the member state's competent authority.

## General Notes and Warnings

Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anesthetic environments).

Do not modify the device.

Use only original HEINE spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not use the device outside.

## Light exposure hazard

 Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures. This device should be used with filters that eliminate UV radiation (<400 nm) and, whenever possible, filters that eliminate short-wavelength blue light (<420 nm).

The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

While no acute optical radiation hazards have been identified for direct or indirect ophthalmoscopes, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

ISO 15004-2: Group 2 (LED)

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum intensity should not exceed 150 minutes (with the use of HEINE A.R. 16D/ Ø54mm A.R. lens).

ANSI Z80.36-2016: Group 2 (LED)

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater is the risk of ocular damage. Exposure to light from this instrument when operated at maximum intensity will exceed the recommended maximum exposure (RME) of 2.2 J/cm<sup>2</sup>, unless additional action is taken by the user to minimize exposure, after \_\_\_ min. (see table exposure guideline). The risk of retinal injury at an exposure of 2.2 J/cm<sup>2</sup> is not high, but because some patients may be more susceptible than others, caution is advised if this radiant exposure value is exceeded. However, because of a significant risk of injury at exposures exceeding 10 J/cm<sup>2</sup>, the user should avoid exposures longer than \_\_\_ min. (see table exposure guideline).

## Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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**Technical specification**

Nominal voltage	3 V DC
Nominal current	max. 190 mA
Environmental conditions for operation	+10 °C to +35 °C 30 % to 90 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	-10 °C to +55 °C 10 % to 95 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-40 °C to +70 °C 10 % to 95 % rel. humidity 500 hPa to 1060 hPa
Device classification according to ISO 15004-2	Group 2
Protection Class	Internally powered

**List of HEINE SIGMA 250 (M2) accessories**

Blue filter	C-000.33.313
Yellow filter	C-000.33.314
Diffusor	C-000.33.315

**Options**

Teaching Mirror	C-000.33.036
-----------------	--------------

## Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	Inside hospitals except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Performance features of the ME system that have been determined to be essential to the performance	None
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the SIGMA 250 (M2), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	See attachment
<b>Test</b>	<b>Compliance</b>
RF emissions CISPR11	Group 1 Class B
Harmonic Emissions	Passed*
Voltage Fluctuations/Flicker	Passed*

\*n/a: „Not applicable“ in the internally powered mode

**Erläuterung der verwendeten Symbole**  
**Explanation of utilized symbols**  
**Explication des symboles utilisés**  
**Explicación de los símbolos utilizados**  
**Spiegazione dei simboli utilizzati**  
**Förklaring av symboler som används**  
**Verklaring van de gebruikte symbolen**  
**Forklaring af de anvendte symboler**  
**Käytettyjen symbolien selitys**  
**Explicação dos símbolos utilizados**

	<p>CE-Kennzeichnung kennzeichnet die Übereinstimmung mit der Europäischen Medizinprodukteverordnung (EU) 2017/745.</p> <p>The CE mark indicates that the product complies with the medical device regulation (EU) 2017/745.</p> <p>Le marquage CE indique que le produit est conforme au règlement (UE) 2017/745 relatif aux dispositifs médicaux.</p> <p>La marca CE indica que el producto cumple el Reglamento europeo sobre productos sanitarios (UE) 2017/745.</p> <p>Il marchio CE indica la conformità con il regolamento sui dispositivi medici (UE) 2017/745.</p> <p>CE-märkning markerar en överensstämmelse med förordningen om medicinska produkter (EU) 2017/745.</p> <p>CE-markering duidt de overeenstemming aan met de verordening betreffende medische hulpmiddelen (EU) 2017/745.</p> <p>CE-mærkningen angiver overensstemmelse med forordningen om medicinsk udstyr (EU) 2017/745.</p> <p>CE-merkintä tarkoittaa, että laite lääkinnällisiä laitteita koskevan asetuksen (EU) 2017/745 kanssa.</p> <p>O símbolo CE identifica a concordância com o regulamento de Dispositivos Médicos (EU) 2017/74.</p>
---	---

	<p>Katalog- oder Bestellnummer</p> <p>Catalogue- or order number</p> <p>Número de catalogue ou de commande</p> <p>Número de catálogo o de pedido</p> <p>Codice catalogo e di dell'ordine numero</p> <p>Katalog- eller Beställningsnummer</p> <p>Catalogus- of Bestelnummer</p> <p>Katalog- eller Ordrenummer</p> <p>Luetelo- tai viitenumero</p> <p>Número de catálogo ou pedido</p>
---	--

	<p>Hersteller</p> <p>Manufacturer</p> <p>Fabricant</p> <p>Fabricante</p> <p>Produttore</p> <p>Tillverkare</p> <p>Fabrikant</p> <p>Produceur</p> <p>Valmistaja</p> <p>Fabricante</p>
---	---

Herstelldatum  
Date of manufacture  
Date de fabrication  
Fecha de fabricación  
Data di produzione  
Tillverkningsdatum  
Productiedatum  
Produktionsdato  
Valmistuspäivä  
Data de fabricação

Getrennte Sammlung von Elektro- und Elektronikgeräten  
(Europäische WEEE Richtlinie)  
Product bearing this symbol may not be disposed of together with general household waste, but instead requires separate disposal according to local provisions.  
(European Waste Electrical and Electronic Equipment Directive, WEEE)  
Tri sélectif des appareils électriques et électroniques  
(Directive européenne DEEE)  
Desechado separado de aparatos eléctricos y electrónicos  
(Directiva Europea RAEE)  
Raccolta differenziata di apparecchi elettrici ed elettronici  
(direttiva europea RAEE)  
Separat insamling av elektriska och elektroniska apparater  
(det europeiska WEEE-direktivet)  
Gescheiden inzameling van elektrische en elektronische apparaten  
(Europese AEEA richtlijn).  
Separat indsamling af elektrisk og elektronisk udstyr  
(det europæiske WEEE-direktiv)  
Sähkö- ja elektroniikkalaitteille tarkoitettu erillinen keräyspiste  
(eurooppalainen WEEE-standardi)  
Colecção separada de aparelhos eléctricos e eletrónicos  
(Diretrizes Europeias WEEE)

Zulässiger Temperaturbereich in °C für Lagerung und Transport  
Temperature limits in °C for storage and transport  
Plage de température admise en °C pour le stockage et le transport  
Rango de temperatura permitida en °C para almacenar y transportar el producto  
Temperatura ammessa in °C per conservazione e trasporto  
Tillåtet temperaturintervall i °C för lagring och transport  
Toegestane temperaturen in °C voor opslag en transport  
Tilladt temperaturområde i °C ved opbevaring og transport  
Näyttää pakkauksen sallitun säilytys- ja kuljetuslämpötilan (°C)  
Limite de Temperatura em °C para armazenamento e transporte

Zulässiger Temperaturbereich in °F für Lagerung und Transport  
Temperature limits in °F for storage and transport  
Plage de température admise en °F pour le stockage et le transport  
Rango de temperatura permitida en °F para almacenar y transportar el producto  
Temperatura ammessa in °F per conservazione e trasporto  
Tillåtet temperaturintervall i °F för lagring och transport  
Toegestane temperaturen in °F voor opslag en transport  
Tilladt temperaturområde i °F ved opbevaring og transport  
Näyttää pakkauksen sallitun säilytys- ja kuljetuslämpötilan (°F)  
Limite de Temperatura permitida em °F para armazenamento e transporte

Zulässige Luftfeuchtigkeit für Lagerung und Transport  
Humidity limitation for storage and transport  
Humidité admise pour le stockage et le transport  
Humedad del aire permitida para almacenar y transportar el producto  
Umidità atmosferica ammessa durante il trasporto e la conservazione  
Tillåten luftfuktighet för transport och lagring  
Toegestane luchtvochtigheid voor opslag en transport  
Tilladt luftfugtighed ved opbevaring og transport  
Sallittu ilmankosteus kuljetuksen ja varastoinnin aikana  
Umidade do ar admissível para o armazenamento e transporte

Zulässiger Luftdruck für Lagerung und Transport  
Pressure limitation for storage and transport  
Pression atmosphérique admise pendant le transport et le stockage  
Presión de aire permitida para almacenar y transportar el producto  
Pressione atmosferica ammessa durante il trasporto e la conservazione  
Tillåten lufttryck för lagring och transport  
Toegestane luchtdruk voor opslag en transport  
Tilladt lufttryk ved opbevaring og transport  
Sallittu ilmanpaine kuljetuksen ja varastoinnin aikana  
Pressão do ar admissível para o armazenamento e transporte

Vorsicht Bruchgefahr!  
Fragile, handle with care!  
Fragile ! Manipuler avec soin  
Atención. Frágil.  
Attenzione: pericolo di rottura!  
Försiktigt! Risk för brott  
Voorzichtig, kans op breuk!  
Forsigtig, risiko for brud!  
Varo särkymisvaaraa!  
Frágil, manuseie com cuidado!

Trocken lagern!  
Keep dry!  
Conserver au sec !  
Conservar en un lugar seco!  
Evitare ambienti umidi!  
Förvaras torrt!  
Droog bewaren!  
Opbevarer tørt!  
Säilytetään kuivassa paikassa!  
Armazenar em ambiente seco!

Gebrauchsanweisung verbindlich befolgen.  
(Hintergrundfarbe: blau, Vordergrundfarbe: weiß)  
Follow instructions for use!  
(Background color: blue, foreground color: white.)  
Suivre le mode d'emploi.  
(Couleur de fond : bleu ; couleur du premier plan : blanc)  
Seguir obligatoriamente las instrucciones de uso.  
(Color de fondo: azul, color de primer plano: blanco)  
Attenersi obbligatoriamente alle istruzioni per l'uso.  
(Colore dello sfondo: blu, colore in primo piano: bianco)  
Bruksanvisningen ska alltid följas.  
(Bakgrundsfärg: blå, förgrundsfärg: vit)  
De gebruiksaanwijzing is bindend en dient gevolgd te worden.  
(achtergrondkleur: blauw, voorgrondkleur: wit)  
Følg altid brugsanvisningen.  
(Baggrundsfarve: Blå; forgrundsfarve: Hvid)  
Käyttöohjetta on noudatettava tarkasti.  
(Taustaväri: sininen, etualan väri: valkoinen)  
Siga as instruções de uso!  
(Cor de fundo: azul, cor de primeiro plano: branco)

**UDI**

Unique Device Identification

**MD**

Medical Device

**HEINE Optotechnik GmbH & Co. KG**

Dornierstr. 6, 82205 Gilching, Germany

[www.heine.com](http://www.heine.com)