**HEINE** QUALITY  
MADE IN GERMANY

## HEINE® NC2 Dermatoscope

 **HEINE Optotechnik GmbH & Co. KG**  
Dornierstr. 6 · 82205 Gilching · Germany  
E-Mail: [info@heine.com](mailto:info@heine.com) · [www.heine.com](http://www.heine.com)  
MED 113857      2020-06-02

V-200.00.068

# HEINE® NC2 Dermatoscope

DEUTSCH	3
ENGLISH	6
FRANÇAIS	9
ESPAÑOL	12
ITALIANO	15
SVENSKA	18
NEDERLANDS	21
DANSK	24
NORSK	27
SUOMI	30
PORTUGUÊS	33

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# HEINE® NC2 Dermatoscope

 Please read and follow these instructions for use of and keep them for future reference.

## Intended Use

### NC2 Dermatoscope

The HEINE® NC2 Dermatoscope is an internally powered medical examination light for non-invasive visual inspection of intact skin. The product can be used contactless or in contact with the skin of the patient. It should only be operated by qualified medical personnel within a professional healthcare facility.

### NC2 Dermatoscope with "Smartphone Adapter" option

The "Smartphone Adapter" option can be used to mechanically connect a smartphone to the NC2 dermatoscope. The NC2 App is designed for documentation through the recording, storage and display of images.

 **For U.S. only:**  
Federal law restricts this device to sale by or on the order of a Physician or Practitioner!


## Warnings and Safety Information

 **CAUTION!** This symbol indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** This symbol indicates valuable advice in terms of set up, operation or maintenance, as applicable. Notes are important, but not related to hazardous situations.

## Product overview

- 1 Instrument
- 2 Focusing ring
- 3 Contact plate
- 4 Adapter
- 5 Mounting case smartphone
- 6 Micro-USB connector
- 7 Charge indicator

## Setting up

Charge the instrument fully before first use.

## Charging

Use a USB cable to connect the USB connector (6) of the NC2 to a USB power supply (e.g. HEINE® E4-USB, PC or laptop). The battery indicator (7) flashes in green during the charging process. When the indicator is green permanently, the unit is fully charged. If the indicator is orange, you should recharge your device.

## Operation

To switch on the illumination, press the button on the instrument. (Press the button again to switch the illumination off.)

The NC2 is intended for a transient application of less than < 20 minutes with a 10 minutes break until the next application.

### *Contact mode (with contact plate)*

Place the contact plate (3) on the instrument such that it snaps into the recesses. The contact plate is held in place by magnetic force.

Place the instrument with the contact plate onto the skin area to be investigated. Because of the polarised illumination, you do NOT need any contact liquid such as immersion oil or disinfection spray. The image has 10-fold magnification.

### *Non-contact mode (without contact plate)*

Take the contact plate off and hold the unit at a distance of approximately 2 cm over the skin area to be investigated. The image has 6-fold magnification.

### *Adjusting the focus*

Your eye needs to be at a distance of approximately 2 cm from the focusing ring (2). Rotate the focusing ring until you get a sharp image.

### *Smartphone adapter*

To connect the NC2 to your smartphone, you need the adapter dock (5) and adapter (4) fitting your smartphone.

Place the smartphone in the adapter dock (5) and press it into the dock until it snaps into the brackets. Push the adapter (4) onto the adapter dock until it clicks into place. Rotate the focusing ring (2) outward to the limit stop such that the eyepiece focuser is moved out as far as possible. Carefully push the fixing bracket of the adapter over the eyepiece focuser until the bracket snaps into place and the adapter is firmly connected to the NC2.

The adapter can be rotated about the eyepiece focuser allowing the smartphone to be oriented in any direction (for right- and left-handed users).

Install the NC2 app on your smartphone. The app is available in the App Store\*.

To avoid any possible malfunctions, please deactivate all services (flight mode).

 Protecting patient data is users' sole responsibility (e. g. Protection of the smartphone with a pin, data transmission via a secure connection). Ensure that your smartphone is adequately protected against access by third parties.

An electrical conducting connection between smartphone, PC and a mains power source is not permissible during the recording of the picture.

 The pictures are stored in a file in the picture database of the smartphone. This file should only be viewed and changed with the NC2 app.

Carry out periodically a back-up of your data to avoid a loss. A back-up of the data is for example possible with the iCloud.

## Hygienic reprocessing

The instruction is available:

- at [www.heine.com](http://www.heine.com)
- in a paper version which you can request from the address listed

## Maintenance

The device is maintenance free.

## Service

The device has no components serviceable by the end-user.

\*App Store is a service mark of Apple Inc.

## General notes

 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

The expected life cycle amounts with designated use and the observation of warning and safety information as well as the maintenance instructions up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

## General warnings

 Use only CE marked USB power supplies (5 V) from reputable manufacturers that conform with the safety requirements from the IEC 60950-1 Information technology or IEC 60601-1 Medical electrical equipment.

Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

This product is not allowed to be entered or used in strong magnetic fields like MRI scanners.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Do not look directly into the light source to avoid dazzle from the intense light.

The NC2 is not suitable for eye examination.

## Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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## Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	Inside hospitals except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Performance features of the ME system that have been determined to be essential to the performance	None
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the NC2, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	See attachment
<b>Test</b>	<b>Compliance</b>
RF emissions CISPR11	Group 1 Class B
Harmonic Emissions	Passed*
Voltage Fluctuations/Flicker	Passed*

\*n/a: "Not applicable" in the internally powered mode

## Technical specification

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Power supply	5 V (± 10 %), 500 mA (Micro-B USB)
Power consumption	max. 1,2 W
Class	Charging: II Operation: internally powered
Charging time	typ. 2 h
Operating time	< 20 minutes with a 10 minutes break
Device classification according to IEC 62471	Group 2
Applied part	Type BF (contact mode)
Protection class	IP20

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