HEINE QUALITY
MADE IN GERMANY

⚡ HEINE SIGMA® 250 (M2)

CE

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49(0)81 52/38-0
Fax +49(0)81 52/38-202
E-Mail: info@heine.com · www.heine.com
MED 113379 **2017-04-13**


V-200.00.066

HEINE SIGMA® 250 (M2)

DEUTSCH	3
ENGLISH	7
FRANÇAIS	11
ESPAÑOL	15
ITALIANO	19
SVENSKA	23
NEDERLANDS	27
DANSK	31
SUOMI	35
PORTUGUÊS	39

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 Please read and follow these instructions for use of and keep them for future reference.

Intended Use

The indirect Ophthalmoscope HEINE SIGMA® 250 (M2) is a battery powered indirect ophthalmoscope for medical professionals, containing illumination and viewing optics intended to examine the media (cornea, aqueous, lens, vitreous) and the retina of the eye. The instrument may be used by professional medical staff and in professional healthcare institutions only.

 **For U.S. only:**
Federal law restricts this device to sale by or on the order of a Physician or Practitioner!

Warnings and Safety Information

 **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** Indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

Product overview


Name of positions

HEINE SIGMA® 250 (M2):

- 1 Pressure clasp
- 2 Release
- 3 Aperture lever
- 4 Filter lever
- 5 Filter blue /yellow, diffuser
- 6 Stereopsis control
- 7 Parallax control
- 8 Beam control
- 9 Teaching mirror accessory
- 10 Dust cover
- 11 PD
- 12 Ocular
- 13 2D-Lens

S-FRAME:

- 14 Swivel mount
- 15 Temple / cord retainer
- 16 Nose rest

Setting up

To set up the HEINE SIGMA® 250 (M2), couple the connection cable with the HEINE® mPack mini. The setup and operation of the HEINE® mPack mini are described in a separate instruction document.

The spectacle frame has a vertical swivel mount (14), onto which the HEINE SIGMA® 250 (M2) can be fitted from above by means of the pressure clasp (1). To remove the HEINE SIGMA® 250 (M2), press the release (2) and pull them up and out of the mounting slot. Attach the bulb cord to the frame temple by means of the cord retainers enclosed (15). The retaining cord can be simply pushed on to the ends of the temples. We recommend that spectacle-wearers have purchase the clip-in correction frame and have it fitted with their prescription.

Operation

Light spot apertures:

With lever (3) you can select two different-sized spots.

Filters:

With the filter lever (4) you can add a red-free filter to any light spot. The optional blue or yellow filter or diffuser accessory (5) can be push-fitted to the front of the instrument housing.

Stereopsis and Parallax:

The stereopsis control (6) below the instrument body adjusts the observation beam to suit the size of the patient's pupil (forward position close to patient: large pupil). The parallax control (7) adjusts the illumination beam to suit the size of the patient's pupil (upper position: large pupil). With both controls on the small setting it is possible to examine in stereopsis with 2 mm pupils. We recommend large stereopsis and small parallax settings for the examination of the periphery.

Lighting direction:

The rotary beam control (8) can be used to swivel the beam vertically to improve reflex reduction and retinal illumination.

Teaching mirror:

The optional teaching mirror accessory (9) can be push-fitted from below to the front of the housing.

Basic setting:

- Remove the dust cover (10).
- Connect the power cord to your power source and switch on the illumination.s
- Carry out the following procedure:
 - Aperture lever (3) – large spot
 - Filter lever (4) – white light
 - Stereopsis control (6) – large pupil
 - Parallax control (7) – large pupil
 - Beam control (8) – middle of the range of adjustment
 - PD (11) – select your PD or set to an average position

Adjusting the optics

Put on the spectacle and check to ensure that the eyepieces are at the correct height. Remove the instrument and bend the nose rest (16) carefully to adjust the height. Bend the temple (15) and if necessary attach the retaining cord to the frame and adjust so that the instrument fits comfortably.

Setting the optics

- Move the HEINE SIGMA® 250 (M2) as close as possible to your eyes and look at the light spot at a distance of 30 – 40 cm. A small object like a pencil must be clearly-focused. Adjust the HEINE SIGMA® 250 (M2) so that the light spot appears centred vertically in your view.
- Adjust the PD setting by viewing the light spot alternately with the left and right eye and sliding the ocular (12) so that the spot is central. If you look through the ocular with both eyes, there should be no diplopic images.
- Correct adjustment of the optics is particularly important when examining through small pupils. Once set, the adjustments need only be changed if another examiner uses the instrument.
- Each instrument is supplied with separate plano lenses which can be used to replace the + 2D lenses if preferred. Change the ocular lenses only in a clean environment in order to avoid dust within the device.

Setting the brightness

Use the control dial on the HEINE® mPack mini to adjust the brightness of the HEINE SIGMA® 250 (M2):

Position I: maximum brightness

Position 0: OFF

Hygienic Reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

Classification according to KRINKO: noncritical

Spaulding Classification USA: noncritical

 Allow the device to cool down before reprocessing.

Before cleaning, disconnect the device from the power source.

In the event of suspected contamination, carry out hygienic preparation of the instrument.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the agents and procedures mentioned below.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the instructions of the manufacturer of the reprocessing media.

Do not use spray or immersion disinfection, dripping wet or heavily foaming tissues, automated or ultrasonic reprocessing.

Procedure

Clean and disinfect the HEINE SIGMA® 250 (M2) and the HEINE® A.R. lens manually (wipe clean and wipe disinfect).

Recommended agents:

Cleaning agent: enzymatic (e.g. neodisher® MediClean)

Disinfectant agent: alcoholic (e.g. Incides® N) or agent hydrogen peroxide (e.g. PREempt® Wipes)

Changing the light source

With the HEINE SIGMA® 250 (M2) the LED cannot be changed.

Maintenance and Service

The instruments do not require maintenance or service.

We as manufacturer do not stipulate a safety inspection in accordance with § 6 Medical Device Operator Ordinance (MPBetreibV), ref. attachment 1, for the HEINE product.

The expected operational life is 7 years.

General Notes and Warnings

 Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anesthetic environments).

Do not modify the device.

Use only original HEINE spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not use the device outside.

General Notes

 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

Light exposure hazard

 Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures. This device should be used with filters that eliminate UV radiation (< 400 nm) and, whenever possible, filters that eliminate short-wavelength blue light (<420 nm). The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

While no acute optical radiation hazards have been identified for direct or indirect ophthalmoscopes, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum intensity should not exceed 150 minutes (with the use of HEINE® A.R. 16D/Ø54mm A.R. lens).

Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

Electromagnetic Compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

 This device is intended for use by medical professionals in the electromagnetic environment specified below. The user of this device should assure that it is used in such an environment.

The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission reduced electrical immunity of the medical equipment.

The device may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

The appendix contains the tables „Guidance and Manufacturer's declaration – Electromagnetic immunity“, technical data and an explanation of the used symbols.

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Appendix

Technical specification

Nominal voltage	3 V DC
Nominal current	max. 190 mA
Environmental conditions for operation	+10 °C to +35 °C 30 % to 90 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	-10 °C to +55 °C 10 % to 95 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-40 °C to +70 °C 10 % to 95 % rel. humidity 500 hPa to 1060 hPa
Device classification according to ISO 15004-2	Group 2
Protection Class	Internally powered

List of HEINE SIGMA® 250 (M2) accessories

Teaching Mirror	C-000.33.036
Blue filter	C-000.33.313
Yellow filter	C-000.33.314
Diffusor	C-000.33.315

Guidance and manufacturer's declaration – electromagnetic emissions

The SIGMA 250 (M2) is intended for use in the electromagnetic environment specified below.
The customer or the user of the SIGMA 250 (M2) should assure that it is used in such environment.

Emission test	Compliance	Electromagnetic environment – Guidelines
RF emissions CISPR 11	Group 1	The SIGMA 250 (M2) uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The SIGMA 250 (M2) is suitable for use in all establishments, excluding domestic establishments and those directly connected to the public mains network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	not applicable	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	not applicable	

Guidance and manufacturer declaration - Electromagnetic immunity

The SIGMA 250 (M2) is intended for use in the electromagnetic environment specified below.
The customer or the user of the SIGMA 250 (M2) should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guideline
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	not applicable	
Surge IEC 61000-4-5	± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth	not applicable	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T , (>95% dip in U_T) for 1/2 period 40% U_T , (60% dip in U_T) for 5 periods 70% U_T , (30% dip in U_T) for 25 periods <5% U_T , (>95% dip in U_T) for 5 seconds	not applicable	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Comment: U_T is the a.c. supply voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The SIGMA 250 (M2) is intended for use in the electromagnetic environment specified below. The customer or the user of the SIGMA 250 (M2) should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	not applicable	Portable and mobile RF communication equipment should be used no closer to any part of the SIGMA 250 (M2), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 3,5/ n/a * \text{SQRT} (P/W)$ $d = 3,5/3 * \text{SQRT} (P/W)$ 80 MHz to 800 MHz $d = 7/3 * \text{SQRT} (P/W)$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ME Equipment is used exceeds the applicable RF compliance level above, the ME Equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SIGMA 250 (M2).
- b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than n/a.

Recommended separation distances for portable and mobile RF communication equipment and the SIGMA 250 (M2)

The SIGMA 250 (M2) is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the SIGMA 250 (M2) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the SIGMA 250 (M2) as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3,5/ n/a * \text{SQRT} (P)$	80 MHz to 800 MHz $d = 3,5/3 * \text{SQRT} (P)$	800 MHz to 2.5 GHz $d = 7/3 * \text{SQRT} (P)$
0.01	n/a	0.12	0.23
0.1	n/a	0.38	0.73
1	n/a	1.2	2.3
10	n/a	3.8	7.3
100	n/a	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

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