**HEINE** QUALITY  
MADE IN GERMANY

## HEINE® NT4

CE

**HEINE Optotechnik GmbH & Co. KG**  
Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49(0)81 52/38-0  
Fax +49(0)81 52/38-202  
E-Mail: [info@heine.com](mailto:info@heine.com) · [www.heine.com](http://www.heine.com)  
**MED 112353      2016-06-30**

# HEINE® NT4

DEUTSCH	3
ENGLISH	7
FRANÇAIS	11
ESPAÑOL	15
ITALIANO	19
SVENSKA	23
NEDERLANDS	27
DANSK	31
NORSK	35
SUOMI	39
PORTUGUÊS	43

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# HEINE® NT4

 Please read and follow these instructions for use of and keep them for future reference.

## Intended Use

The HEINE charger NT4 is to be used exclusively for charging HEINE 4 series rechargeable handles. These must include a HEINE rechargeable battery.

 **For U.S. only:**  
Federal law restricts this device to sale by or on the order of a Physician or Practitioner!

## Warnings and Safety Information


 **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** Indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

## Product overview

## Setting up

### Primary plug:

The mains socket adapter of the power supply is interchangeable, this enables worldwide operation. The pictograms describe the procedure of changing the adapters.

Place the power supply in such a position where it is accessible and easily unplugged.

Do not place the device on soft surfaces or fabrics.

Apart from the appropriate charging handles, no foreign objects (e.g. paperclips) should be placed into the charging slots.

By plugging the mains transformer into the mains socket, the charger is ready for use.

The charger ensures that batteries are kept in good condition irrespective of handle size.

Instruments placed next to each other must not touch one another.

When using different sized handles, reduction rings can be used or exchanged to ensure an optimal fit for all handles at any time.

Reduction rings for:

BETA 4 NT

BETA 4 SLIM / F.O. 4 SLIM LED NT / F.O. 4 SLIM NT

Standard Laryngoscope handles do not require a reduction ring.

After use, the handle can be replaced into the charging slot and remain there. Handles cannot be overcharged.

The device must be disconnected from the mains supply when the adaptors are being fitted.

If the second charging compartment is not in use it can be sealed off by the dust cap which is supplied with the unit.

The pulsating charge indicator (light ring) around the charging slot, indicates that a charging current is flowing and the battery is charging. The charging process takes place in accordance with the JEITA\* temperature standard for Li-ion batteries. This can lead to short-term charging interruptions at ambient temperatures > 30 °C which optimizes the battery life.

(\* Japan Electronics and Information Technology Industries Association)

A constant light around the charge indicator indicates that a HEINE handle is present but no charging process is taking place.

The charge indicator does not light up when the charging slot is empty or a malfunction is present.

If a fully charged handle remains in the charger, there will be brief recharging cycles (< 5 min) this is visible from the flashing of the charge indicator. The handle is fully charged and ready for use.

The LED indicates the operational readiness of the power supply.

**⚠** The charger may only be operated with the mains voltage specified on the nameplate of the power supply.

Unplug the unit from the mains if not in use over long periods.

In the event of fluid ingress in the housing, the charger should be taken out of service and checked by an authorised person.

Follow the instructions for the charger and rechargeable handles.

When replacing the handles in the charging compartments, ensure that the internal contacts in the compartment are clean and that the instruments are switched off.

The rechargeable battery can only be charged when fitted to the handle.

The HEINE Charger is intended for installation in medical areas outside the patient environment (at least 1.5 metres from the patient or patient support pursuant to IEC 60601-1, see figure).


Fig1: The above figure defines the patient environment (minimum radius around the patient).

To prevent mechanical and electrical breakdown, never install or store the charger in the following areas:

- areas exposed to high moisture or condensation,
- areas exposed to extreme environmental conditions,
- areas subject to constant vibrations,
- areas severe temperature fluctuations,
- outdoors.

## Hygienic Reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

Classification according to KRINKO: non-critical

Spaulding Classification USA: noncritical

 Before cleaning, disconnect the device from the power source.

In the event of suspected contamination, carry out hygienic preparation of the instrument.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the cleansing agents and procedures named in these instructions for use.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the instructions of the manufacturer of the reprocessing media.

Do not use spray or immersion disinfection, dripping wet or heavily foaming tissues, automated or ultrasonic reprocessing.

 Place only reprocessed handles into the charger.

Hold the charger with the slots facing down when cleaning and disinfecting.

## Procedure

Clean and disinfect the device manually (wipe clean and wipe disinfect).

### Recommended agents

Cleaning agent: Neodisher® MediClean

Disinfectant agent:

EU: quaternary ammonia compounds (e.g. Microbac® tissues)

USA: Active ingredient orthophthalaldehyde (for example Cidex® OPA)

## Maintenance

The device is maintenance free. However, portable appliance testing must be carried out at least every 24 months by a qualified electrician.

 The power supply must not be opened (risk of electrical shock).

## Service

The device does not require service.

## General Warnings

 Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Turn the handle off during the charging process.

The device must not be placed near strong Magnetic fields, e.g. MRI units.

## General Notes

 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

Do not pull on the cord/cable to disconnect the power supply from the mains.

## Disposal

The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

## Electromagnetic Compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

**⚠** This is a device in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site.

The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission reduced electrical immunity of the medical equipment.

The device may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

The appendix contains following tables

- Guidance and Manufacturer's declaration – Electromagnetic immunity
- Technical specification
- Explanation of the used symbols

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**Guidance and manufacturer's declaration – electromagnetic emissions**

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – Guidelines</b>
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Passed	

**Guidance and manufacturer declaration - Electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – Guideline</b>
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	± 2 kV for mains cables ± 1 kV for input and output lines	The supply voltage quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% $U_T$ , (>95% dip in $U_T$ ) for 1/2 period 40% $U_T$ , (60% dip in $U_T$ ) for 5 periods 70% $U_T$ , (30% dip in $U_T$ ) for 25 periods <5% $U_T$ , (>95% dip in $U_T$ ) for 5 seconds	< 5% $U_T$ , (>95% dip in $U_T$ ) for 1/2 period 40% $U_T$ , (60% dip in $U_T$ ) for 5 periods 70% $U_T$ , (30% dip in $U_T$ ) for 25 periods <5% $U_T$ , (>95% dip in $U_T$ ) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Comment:  $U_T$  is the a.c. supply voltage prior to application of the test level.

**Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	3 V eff	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	3 V/m	<p><b>Recommended separation distance:</b>  <math>d = 3,5/3 * \text{SQRT}(P/W)</math>  <math>d = 3,5/3 * \text{SQRT}(P/W)</math> 80 MHz to 800 MHz  <math>d = 7/3 * \text{SQRT}(P/W)</math> 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

**Recommended separation distances for portable and mobile RF communication equipment and the device**

The device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3,5/3 * \text{SQRT}(P)$	80 MHz to 800 MHz $d = 3,5/3 * \text{SQRT}(P)$	800 MHz to 2,5 GHz $d = 7/3 * \text{SQRT}(P)$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**Technical specification**

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Mains supply	100 – 240 V AC / 50 – 60 Hz
Current consumption	max. 205 mA
Output	max. 6 V / 0,9 A
Fuse	Integral overload protection
Charging time	max. 6 h
Safety class	II
Protection class	IP 20

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### Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

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