HEINE QUALITY
MADE IN GERMANY

HEINE® iC2 Funduscope

 HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49(0)81 52/38-0
Fax +49(0)81 52/38-202
E-Mail: info@heine.com · www.heine.com
MED 113539 2018-03-07

HEINE® iC2 Funduscope

DEUTSCH	3
ENGLISH	7
FRANÇAIS	11
ESPAÑOL	15
ITALIANO	19
SVENSKA	23
NEDERLANDS	27
DANSK	31
NORSK	35
SUOMI	39
PORTUGUÊS	43

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HEINE® iC2 Funduscope

 Please read and follow these instructions for use of and keep them for future reference.

Intended Use

The HEINE® iC2 Funduscope is a hand-held, non-invasive, internally powered digital ophthalmoscope, which is intended for recording, saving and displaying pictures of the internal part of the eye, in particular the ocular media and the ocular background.

The device is intended for transient application. It should only be operated by qualified medical personnel within a professional health facility.

 **For U.S. only:**
Federal law restricts this device to sale by or on the order of a Physician or Practitioner!

Warnings and Safety Information

 **CAUTION!** This symbol indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** This symbol indicates valuable advice in terms of set up, operation, maintenance or repair, as applicable. Notes are important, but not related to hazardous situations.

Product overview

- 1 Controller for adjusting the dioptre compensation
- 2 Double-stage release button for switch on the device and for picture recording
- 3 Micro USB-port
- 4 Charger indicator
- 5 Data matrix code and serial number of the device
- 6 Mounting case smartphone
- 7 Protective cover

Setting up

Installation of the app

Install the iC2 app on your smartphone. The app is available in the AppStore.

To avoid any possible malfunctions, please use the smartphone without SIM-card.

Choose a mounting case, which is suitable for the used smartphone. Insert the smartphone into the appropriate mounting case so it is firmly seated. Make sure that the camera is in the correct position.

Setting up of the connection between the iC2 and the smartphone

The connection between the iC2 and the smartphone is possible, using the data matrix code or the choice of the available devices in the list.

- Start the iC2 App on your smartphone
- Switch on the iC2, using the release button (second pressure point)
- The description of the further steps for the connecting the smartphone with the iC2 is included in the app.
- The connection is via *Bluetooth*® wireless technology

Connecting the mounting case smartphone with the iC2

Insert the guide pins, attached to the mounting case, into the guide rails of the iC2 so it is firmly seated.

 Make sure that the mounting case is securely attached to the iC2.

Operation

- Use the iC2 only with the corresponding app.
- Adjust the brightness as low as necessary during the application and ensure that the light is switched off after each application.

The iC2 is controlled by operating elements located on the iC2 and in the iC2 app.

The device provides the following functions:

- Create a patient file / Recording of patient data
- Picture recording
- Filter selection
- Storage of pictures
- Three-level brightness setting via the setting menu of the app.

 Protecting patient data is users' sole responsibility (e. g. Protection of the smartphone with a pin, data transmission via a secure connection). Ensure that your smartphone is adequately protected against access by third parties.

An electrical conducting connection between smartphone, PC and a mains power source is not permissible during the recording of the picture.

-  The pictures are stored in a file in the picture database of the smartphone. This file should only be viewed and changed with the iC2 app.
- Carry out periodically a back-up of your data to avoid a loss. A back-up of the data is for example possible with the iCloud.

Switch off the device

Switch off the device after each application. This is possible through the following procedure:

- Closing the iC2 app causes the device to switch off automatically
- If there is no *Bluetooth*[®] connection between the iC2 and the app for five minutes and the double-stage release button isn't pushed, the device will switch off automatically.

Reset the device

If the device does not react, it can be reset by pushing the release button for 30 seconds. This procedure has no influence on the data of the iC2 app.

Application duration

The device can be used continuously in an ambient temperature of between 10 °C and 30 °C. In an ambient temperature of between 30 °C and 35 °C, the iC2 can be used for a brief period of < 15 minutes and with a break of 10 minutes before being used again.

Charging

Charging via USB

Please connect the USB cord, which is attached to the power supply, to the micro USB charging port (3).

Wireless charging

When using the iC2 charging station, insert the device into the charging station until it lies full-facedly.

-  Pay attention to the fact that no foreign objects are between the device and the charging station.
- When replacing the iC2 in the charger, ensure that the charger and the iC2 aren't contaminated.

The setting up and operation of the iC2 charging station is described in a separate instruction of use. Due to local regulations, the iC2 charging station may not be available in your country.

Charge Indicator

When the iC2 is charged via USB, the LED will flash green.

When the battery is fully charged, the LED will be constant green.


When using the iC2 charging station, the charging procedure is indicated by the charging station.

When the battery reaches a critical charge state, the LED will light up orange. The device must now be recharged.

Apart from the charge indicator LED, the charging status can also be seen in the settings via the submenu "iC2 Connection". In addition, the charge status is indicated by a bar graph display in the picture capture mode.

⚠ Please ensure that the device is only charged to 80 % in an ambient temperature of 35 °C to preserve the rechargeable battery.

Charge the device outside of the patient environment (1.5 m).

Hygienic reprocessing

The instruction is available:

- online at www.heine.com
- in a paper version which you can request from the address listed

Maintenance

The device is maintenance free.

Service

The device has no components serviceable by the end-user.

General notes

⚠ The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

The expected life cycle, when the device is normal used and the warning and safety information as well as the maintenance instructions are observed, is up to 5 years.

Beyond this period, the product may continue to be used if it is in a safe and good condition.

General warnings

⚠ Use only CE marked USB power supplies (5 V) from reputable manufacturers that conform with the safety requirements from the IEC 60950-1 Information technology or IEC 60601-1 Medical electrical equipment.

Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

This product is not allowed to enter or be used in areas with strong magnetic fields e.g. MRI scanners.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Store and use the device in dry and dust-free environments only!

Light exposure hazard

Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures.

The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

While no acute optical radiation hazards have been identified for the funduscope, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

Disposal

The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	Inside hospitals except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Performance features of the ME system that have been determined to be essential to the performance	None
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the iC2 Funduscope, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	As applicable, see attachment
Test	Compliance
RF emissions CISPR11	Group 1 Class B
Harmonic Emissions	Passed*
Voltage Fluctuations/Flicker	Passed*

*n/a: "Not applicable" in the internally powered mode

Technical specification

Environmental conditions for operation	+10 °C to +35 °C 30 % to 90 % rel. humidity 800 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for transport	-10 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Work distance	30 ± 5 mm
Nominal charging voltage	5 V (± 10 %), 500 mA (Micro-B USB) Wireless charging
Class	Charging mode: Class II Operation: internally powered
Device classification according to ISO 15004-2:	Group 1
Operational time	typ. 4 h
Rechargeable battery charging time	USB: typ. 5.5 h Wireless: typ. 4 h

THIS INSTRUMENT DOES NOT WITHSTAND THE TEMPERATURE CONDITIONS FOR STORAGE ACCORDING TO ISO 15004-1. DO NOT STORE THIS INSTRUMENT AT CONDITIONS UNDER WHICH THE TEMPERATURE MAY BE ABOVE 50 °C OR BELOW -10 °C.

Options

USB cord with E4-USB plug-in power supply	X-000.99.303
Mounting case smartphone	C-000.33.15x
iC2 Charger	X-002.99.209
Case for iC2	C-150.00.000

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