HEINE UniSpec® Instrument Head

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Please read and follow these instructions for use and keep them for future reference.

Intended Use The HEINE UniSpec® Instrument Head and the UniSpec® disposable tubes with fiberglass illumination are intended for examination and performance of procedure within the rectum and the anus. It should only be operated by qualified medical personnel within a professional health facility.

For U.S. only: Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

Warnings and Safety Information CAUTION! Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to serious extent. (Background color yellow; foreground color black). NOTE! This symbol indicates valuable advice in terms of set up, operation or maintenance, as applicable. (Notes are important, but not related to hazardous situations).

- Product overview (see figure) 1 Viewing window 2 Viewing screw 3 Hole for the optical magnification 4 Extension-piece 5 Connector for the Hygiene Filter 6 Instrument Head for UniSpec® Tubes 7 Hygiene Filter 8 Double bulb insufflator 9 Obturator 10 Tube with extension-piece 11 Swivel lens 12 Illumination adapter

Setting up Loosen the fastening screw (2) and open the viewing window (1). If the viewing window isn't necessary, it can be removed by extracting it upwards out of the instrument head.

Pull the obturator out of the tube. Insert the tube with the extension-piece (4) as far as possible into the instrument head and secure this by twisting to the right. Insert the obturator (9) through the instrument head into the tube (10).

The instrument head can be used with the rubber insufflation bulb (8). To do so, connect the rubber insufflation bulb with the hygiene filter (7) to the viewing window (5). Control the connection. If it isn't securely attached, fix it more tightly.

Risk of allergic shock: The orange rubber insufflation bulb contains natural rubber latex. Do not use it in patients with a latex allergy.

The black insufflation bulb is latex- and phthalate free.

Operation When the tube is inserted, remove the obturator through the instrument head. The viewing window is equipped with a face plate and can be conveniently opened and closed at any time.

Application duration The sigmoidoscopes, proctoscopes and anosopes are intended for transient application. The application time should not exceed 15 minutes.

Hygienic reprocessing The instruction is available: - at www.heine.com - in a paper version which you can request from the address listed

Maintenance The device is maintenance free.

Service The device has no components serviceable by the end-user.

General Notes The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE.

General Warnings Check the correct operation of the device before use! Do not use the device if there are visible signs of damage. Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

Disposal The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

The appendix contains following tables - Electromagnetic disturbances - Requirements and tests - Technical specification - Explanation of the used symbols

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
Electromagnetic disturbances – Requirements and tests	
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.	
Statement for the operational environments	Inside hospitals except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Performance features of the ME system that have been determined to be essential to the performance	None
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the UniSpec instrument head, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	As applicable, see attachment
Test	Compliance
RF emissions CISPR11	Group 1 Class B
Harmonic Emissions	Passed*
Voltage Fluctuations/Flicker	Passed*

*n/a: „Not applicable“ in the internally powered mode

Erläuterung der verwendeten Symbole
Explanation of utilized symbols
Explication des symboles utilisés
Explicación de los símbolos utilizados
Spiegazione dei simboli utilizzati

	CE-Kennzeichnung kennzeichnet die Übereinstimmung mit der Europäischen Medizinprodukterichtlinie 93/42 EWG The CE mark indicates that the product complies with the European medical device directive 93/42/EEC. Le marquage CE indique la conformité à la directive européenne 93/42/CEE relative aux dispositifs médicaux. El marcado CE indica la conformidad con la directiva europea 93/42 /CEE relativa a los productos sanitarios. Il marchio CE indica la conformità con la direttiva europea sui dispositivi medici 93/42 CEE
	Katalog- oder Bestellnummer Catalogue- or order number Numéro de catalogue ou de commande Número de catálogo o de pedido Codice catalogo e di dell'ordine numero
	Seriennummer Serial number Numéro de série Número de serie Numero di serie
	Hersteller Manufacturer Fabricant Fabricante Produttore
	Hersteldatum Date of manufacture Date de fabrication Fecha de fabricación Data di produzione
	Getrennte Sammlung von Elektro- und Elektronikgeräten (Europäische WEEE Richtlinie) Product bearing this symbol may not be disposed of together with general household waste, but instead requires separate disposal according to local provisions. (European Waste Electrical and Electronic Equipment Directive, WEEE) Tri sélectif des appareils électriques et électroniques (Directive européenne DEEE) Desechado separado de aparatos eléctricos y electrónicos (Directiva Europea RAEE) Raccolta differenziata di apparecchi elettrici ed elettronici (direttiva europea RAEE)
	Zulässiger Temperaturbereich in °C für Lagerung und Transport Temperature limits in °C for storage and transport Plage de température admise en °C pour le stockage et le transport Rango de temperatura permitida en °C para almacenar y transportar el producto Temperatura ammessa in °C per conservazione e trasporto
	Zulässiger Temperaturbereich in °F für Lagerung und Transport Temperature limits in °F for storage and transport Plage de température admise en °F pour le stockage et le transport Rango de temperatura permitida en °F para almacenar y transportar el producto Temperatura ammessa in °F per conservazione e trasporto
	Zulässige Luftfeuchtigkeit für Lagerung und Transport Humidity limitation for storage and transport Humidité admise pour le stockage et le transport Humedad del aire permitida para almacenar y transportar el producto Umidità atmosferica ammessa durante il trasporto e la conservazione

Technical specifications	
Environmental conditions for operation	+10 °C to +35 °C 10 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Protection class	IP 20

Suitable power/light source		
HEINE BETA Handle 3.5 V and the Illumination Adapter	X-002.98.002	
BETA 4 NT rechargeable handle	X-007.99.396	
BETA 4 USB Li-ion rechargeable handle	X-007.99.387	
BETA 4 SLIM NT rechargeable handle	X-007.99.474	
BETA L Li-ion rechargeable handle	X-007.99.395	
BETA NT rechargeable handle	X-002.99.411	
BETA SLIM rechargeable handle	X-002.99.471	
Lamp Handle with switch and XHL Xenon Halogen bulb	X-004.99.607	
Sigmoidoscope / Proctoscope KIT RE7000	E-096.15.501	

Accessories	
HEINE UniSpec® Disposable Anoscopes 25 pcs. / 85x20mm 100 pcs. / 85x20mm	E-003.19.925 E-003.19.911
HEINE UniSpec® Disposable Proctoscopes 25 pcs. / 130x20mm 100 pcs. / 130x20mm	E-003.19.825 E-003.19.811
HEINE UniSpec® Disposable Sigmoidoscopes 25 pcs. / 250x20m 100 pcs. / 250x20mm	E-003.18.825 E-003.18.811
Rubber insufflation bulb (orange)	E-003.18.105
Insufflation bulb (black)	E-003.18.105
Hygiene Filter (25 pcs.)	E-000.18.116
HEINE Illumination adaptor	X-002.98.002

Options	
Swivel lens (only suitable for sigmoidoscopes)	E-003.18.096
Swivel lens (only suitable for proctoscopes and anoscopes)	E-003.19.099
Viewing window	E-003.18.099

	Zulässiger Luftdruck für Lagerung und Transport Pressure limitation for storage and transport Pression atmosphérique admise pendant le transport et le stockage Presión de aire permitida para almacenar y transportar el producto Pressione atmosferica ammessa durante il trasporto e la conservazione
	Vorsicht Bruchgefahr! Fragile, handle with care! Fragile ! Manipuler avec soin Atención. Frágil. Attenzione: pericolo di rottura!
	Trocken lagern! Keep dry! Conservér au sec ! Conservar en un lugar seco! Evitare ambienti umidi!
	Gebrauchsanweisung verbindlich befolgen. (Hintergrundfarbe: blau, Vordergrundfarbe: weiß) Follow instructions for use! (Background color: blue, foreground color: white.) Suivre le mode d'emploi. (Couleur de fond : bleu ; couleur du premier plan : blanc) Seguir obligatoriamente las instrucciones de uso. (Color de fondo: azul, color de primer plano: blanco) Attenersi obbligatoriamente alle istruzioni per l'uso. (Colore dello sfondo: blu, colore in primo piano: bianco)
	Anwendungsteil des Typs BF Type BF applied part Partie appliquée de type BF Pieza de aplicación del tipo BF Applicazione di tipo BF
	Einmalgebrauch Single use only Usage unique Un solo uso Monouso
	Enthält oder Anwesenheit von Naturkautschuklatex Contains or presence of natural rubber latex Contient ou présente des traces de latex de caoutchouc naturel Contiene o presenta látex de caucho natural Presenza o assenza di lattice di gomma naturale
	Latexfrei Latex free Sans latex Sin látex Senza lattice
	Bei beschädigter Verpackung nicht verwenden Do not use if packaging is damaged Ne pas utiliser si l'emballage est endommagé No utilizar si el envase está dañado Non usare se la confezione è danneggiata
	Verwendbar bis (JJJJ-MM-TT oder JJJJ-MM) Use by (YYYY-MM-DD or YYYY-MM) Utilisable jusqu'au (AAA-MM-JJ ou AAA-MM) Se puede utilizar hasta (AAAA-MM-DD o AAAA-MM) Scadenza (AAAA-MM-GG oppure AAAA-MM)
	Unique Device Identification