**HEINE** QUALITY  
MADE IN GERMANY

## HEINE® LAMBDA 100

**CE**

**HEINE Optotechnik GmbH & Co. KG**  
Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49(0)81 52/38-0  
Fax +49(0)81 52/38-202  
E-Mail: [info@heine.com](mailto:info@heine.com) · [www.heine.com](http://www.heine.com)  
med 2214      2015-02-10

V-200.00.055

# HEINE® LAMBDA 100

DEUTSCH	3
ENGLISH	7
FRANÇAIS	11
ESPAÑOL	15
ITALIANO	19
SVENSKA	23
NEDERLANDS	27
DANSK	31
SUOMI	35
PORTUGUÊS	39

# HEINE® LAMBDA 100 Retinometer

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# HEINE® LAMBDA 100 Retinometer

Please read and follow these instructions for use of and keep them for future reference.

## Intended Use

The HEINE® LAMBDA 100 Retinometer is intended to examine the potential visual acuity of patients with cataracts or other lens opacities. The device contains viewing optics and a battery-powered illumination unit. It must be used only by qualified medical personnel. The devices are intended for short-term examinations.

## For U.S. only:

 Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

## Warnings and Safety Information

 **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** Note indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

## Product overview

- 1 Brow rest
- 2 Light exit
- 3 Grid orientation selector
- 4 Visual acuity selector
- 5 Bulb
- 6 HEINE Handle

## Setting up

To set up the HEINE® LAMBDA 100 Retinometer, screw the instrument head into the HEINE battery handle or plug it on the HEINE rechargeable handle.

Verify that the lamp voltage complies with the supply voltage of the handle.

The coloured marking on the bottom of the lamp shows you the lamp voltage:

White ring = HEINE XHL® 2.5V bulb  
only for use with the HEINE Battery handle

Red ring = HEINE XHL® 3.5V bulb  
only for use with the HEINE Rechargeable Handle

## Operation

Always dim the room light before examining. Eyes with opaque media should always be dilated. This is not necessary where no opacity is present. Explain to the patient as clearly as possible what he will see with the help of the retinometer instruction card supplied with the set. Use the better eye first to get the patient used to the instrument. Make sure that the patient is in a comfortable position. He should look into the far corner of the darkened room.

Switch on the retinometer. If preferred you can set the desired brightness by shining the light into your own eye. (It may be necessary to change the brightness during the examination).

Use the lever (3) to choose the grid orientation and select a low visual acuity value with the selector wheel (4).

Rest the retinometer upright against the patient's forehead and direct the red light beams onto the pupil. At the normal reading distance look alongside or over the top of the retinometer (there is a "window" (2) built into the brow rest (1) and observe the points of light on the patient's cornea.

Move the beam over the pupil by slightly turning or tilting the retinometer until the patient recognizes the grid with the red and black pattern. Agree with the patient how he will indicate the grid orientation (with a hand-movement or verbally). If he mentions circular spots or distorted lines, this could be an indication of a malfunction of the macula.

Change the grid orientation by rotating the lever (3) with your free hand without moving the retinometer.

Increase the visual acuity value until the patient is no longer able to recognize the direction of the grid. The last value which he was able to recognize represents the potential acuity.

General Instructions for acuity-testing with the Retinometer:

Testing potential acuity with retinometers has been the method of choice for some years especially where opacities are present. It is usually possible to get very reliable results provided a few basic preconditions are observed: With very dense opacities below an optotype visus of 0,1 the reliability of the estimate is reduced, the post-operative visus is often better than that obtained by the retinometer. With amblyopic eyes, the tendency is more towards false positives because of the size of the test pattern. With ametropia of more than 6 Dpt. the results can be improved by examining the patient while he is wearing his glasses. The same comments apply to high astigmatism. Further factors which can have a negative influence on the results are nystagmus, tremor, communication problems, senility or debility of the patient. Generally speaking, degenerative conditions can cause considerable problems, which can only be overcome if the examiner has the necessary experience of this technique. (A detailed presentation on the topic, including references, can be found in: Lachenmayer, Potentielle Sehschärfe bei Störungen der brechenden Medien [Potential visual acuity in case of breaking media impairments'], Quintessenz-Verlag Munich 1993).

The HEINE retinometer are intended for short-term examinations < 2 min. with a 15-minute break before the next application.

## Hygienic Reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

Classification according to KRINKO: non-critical

Spaulding Classification USA: non-critical

 Allow the device to cool down before reprocessing.

In the event of suspected contamination, carry out hygienic preparation of the instrument.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the agents and procedures mentioned below. Other procedures/reprocessing agents have to be validated.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the instructions of the disinfectant manufacturer.

Do not use spray or immersion disinfection, dripping wet or heavily foaming tissues. Do not use automated or ultrasonic reprocessing.

### Procedure

Clean and disinfect the retinometer manually (wipe cleaning and wipe disinfection).

Recommended agents:

Cleaning agent: Neodisher® MediClean

Disinfectant agent: alcoholic (e.g. Incides® Tissues)

### Changing the light source

 Verify that the lamp voltage complies with the supply voltage of the handle.

Remove the LAMBDA 100 from the HEINE handle. Unscrew the connecting piece and carefully remove the bulb (5).

Wipe the dome of the new bulb with a clean cloth and place the bulb in the retinometer.

Screw the connecting piece back on.

### Maintenance and Service

The instruments do not require maintenance or service.

### General Warnings

 Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

### General Notes

 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

If you don't use the device for a longer period of time, please remove the batteries in advance.

### Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

## **Electromagnetic Compatibility**

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

 This is a device in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site.

The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission reduced electrical immunity of the medical equipment.

The device may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

The appendix contains following tables:

- Guidance and Manufacturer's declaration – Electromagnetic immunity
- Technical specification
- Explanation of the used symbols

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## Appendix:

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – Guidelines</b>
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	Symmetrical three-phase devices and other devices.
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Passed	

<b>Guidance and manufacturer declaration - Electromagnetic immunity</b>			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – Guideline</b>
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	± 2 kV for mains cables ± 1 kV for input and output lines	The supply voltage quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT, (>95% dip in UT) for 1/2 period 40% UT, (60% dip in UT) for 5 periods 70% UT, (30% dip in UT) for 25 periods <5% UT, (>95% dip in UT) for 5 seconds	< 5% UT, (>95% dip in UT) for 1/2 period 40% UT, (60% dip in UT) for 5 periods 70% UT, (30% dip in UT) for 25 periods <5% UT, (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Comment: UT is the a.c. supply voltage prior to application of the test level.			

**Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	3 V eff	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	3 V/m	<p><b>Recommended separation distance:</b>  <math>d = 3,5/3 * \text{SQRT}(P/W)</math>  <math>d = 3,5/3 * \text{SQRT}(P/W)</math> 80 MHz to 800 MHz  <math>d = 7/3 * \text{SQRT}(P/W)</math> 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80Hz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

**Recommended separation distances for portable and mobile RF communication equipment and the device**

The device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3,5/3 * \text{SQRT}(P)$	80 MHz to 800 MHz $d = 3,5/3 * \text{SQRT}(P)$	800 MHz to 2,5 GHz $d = 7/3 * \text{SQRT}(P)$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Technical specification

Environmental conditions for operation	+10°C to +35°C 30% to 75% rel. humidity 700hPa to 1060hPa
Environmental conditions for storage	+5°C to +45°C 45% to 80% rel. humidity 500hPa to 1060hPa
Environmental conditions for transport	-20°C to +50°C 45% to 80% rel. humidity 500hPa to 1060hPa
Nominal voltage	3.0V – 3.7V
Nominal current	Typ. 440 – 760 mA
Protection class	internal power supply
Device classification according ISO 15004-2	Group:1
Applied part	Type BF
XHL Xenon spare bulb	#077 (2,5V) #078 (3,5V)

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### Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

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