Instruction for use **HEINE DELTAone A** **Dermatoscope**

Please read and follow these instructions for use of and keep them for future reference

Intended Use

The dermatoscope is an internally powered medical examination light for non-invasive transient contact and non-contact examination of the skin. In contact mode, skin showing a pathological change or a change due to a disease can be examined. The small contact plate is not suitable for use on mucosa. The examination of wounds is only permitted in non-contact mode. The contact plate must be removed for non-contact mode.

The product may only be used by qualified medical professionals and in professional healthcare facilities. The application is not limited to any patient demography.

DELTAone A with "Smartphone Adapter" option

The "Smartphone Adapter" option can be used to mechanically connect a smartphone with a respective app. The app is designed for documentation by recording, storing and displaying images.

Warnings and safety information


CAUTION! This symbol indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

NOTE! This symbol indicates valuable advice. This notes are important, but no related to hazardous situations.

Product overview


Setting up

Charge the instrument fully before first use.

Charging

Charge the device outside of the patient environment (1.5 m).

Charging via USB

Use a USB cable to connect the USB connector (4) of the dermatoscope to a USB power supply (e.g. HEINE E4-USB).

Charge indicator

A charge indicator (5) is located on the side or underside of the dermatoscope.

Instrument	Charger	Battery	Charge indicator
Switched off	connected	charging	Green, flashing
		full	green
Switched on	Not connected	full ($\geq 25\%$)	green
		empty ($< 25\%$)	orange
		residual capacity ($< 10\%$)	orange, flashing

If the device is switched off, the charge indicator remains on for a further 3 seconds and then goes off.

Do not look directly into the light source when checking the charge level.

Operation


To switch on the illumination, press the button (1) on the instrument. The device starts with the last selected brightness level and in polarisation mode. To switch off, press and hold the button (1) on the instrument for at least 2 seconds. Press the button (6) to adjust the brightness in three steps.

The device switches off automatically after 10 minutes if the polarisation buttons or the brightness control are not used.

Amber light configuration

By pressing the button (1) there are 3 light configurations possible.

There are 3 configurations to cycle through. If the last setting is reached, the cycle starts again with neutral white.

Contact mode (with contact plate)

Place the contact plate (3, 3a, 3b) on the instrument, making sure that it snaps into place.

Place the instrument with the contact plate onto the area of skin being examined. With DELTAone A, extend the contact plate (3, 3a, 3b) over the focusing ring (2).

An unused GoCover contact plate must be used to examine the mucous membranes. This contact plate must be disposed after use. Under no circumstances should the small contact plate be used for this use. Do not use the contact plate, which is damaged.

When carrying out examinations with polarised illumination, you do NOT need any contact liquid such as immersion oil or disinfection spray. In polarised mode, the polarisation indicator (8) lights up blue. The polarised mode is available in the configuration neutral white and amber.

Use the polarisation button (7) to switch to non-polarised mode. In this case, use a contact fluid such as immersion oil.

Do not use a contact fluid or immersion oil, which contain alcohol.
Do not use a disinfection agents as contact fluid, which contain alcohol.

The non-polarised mode is available in the configuration neutral white and mixed amber & white.

When using the small contact plate (3a), examination is only possible in non-polarised mode.

Adjusting the focus in contact mode

Your eye needs to be at a distance of approximately 2 cm from the focusing ring (2). Rotate the focusing ring until you see a sharp image. The image has 10-fold magnification.

Non-contact mode (without contact plate)

Take off the contact plate (3, 3a, 3b) and hold the device at a distance of approximately 1 cm above the area of skin being examined. Use the polarisation button (7) to switch to non-polarised mode.

Adjusting the focus in non-contact mode

Your eye needs to be at a distance of approximately 2 cm from the focusing ring (2).

With DELTAone A, the focus is adjusted via the distance between the skin and the optics. There is no need to rotate the focusing ring (2). The image has 10-fold magnification.

Smartphone Adapter

To connect the dermatoscope to a smartphone, you need a suitable adapter for your smartphone.

Place the adapter onto the focusing ring (2). The adapter is held in place by magnetic force.

Install the corresponding app on your smartphone. The app is available in the App Store*.

To avoid any possible malfunctions, please deactivate all other services (switch on flight mode).

*App Store is a service mark of Apple Inc.

Protecting patient data is users' sole responsibility (e. g. protection of the smartphone with a pin, data transmission via a secure connection). Ensure that your smartphone is adequately protected against access by third parties.

An electrical conducting connection between smartphone, PC and a mains power source is not permissible during the recording of the picture.

All images are saved exclusively on the smartphone in the app's sandbox. The user is responsible for the regular creation of a data backup to prevent data loss, which can happen for example through deleting the app or losing the smartphone. A data backup can be made via cable to another computer or via the network to the iCloud.

Hygienic reprocessing

The instruction is available:

- at www.heine.com
- in a paper version which you can request from the address listed

Maintenance

The device is maintenance free.

Service

The device has no components serviceable by the end-user.

General Notes

The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

The expected life cycle amounts with designated use and the observation of warning and safety information as well as the maintenance instructions up to 7 years for the DELTAone A. Beyond this period, the product may continue to be used if it is in a safe and good condition.

The GoCover can be used up to 3 years from the date of manufacture.

Note to the user and/or patient:

All serious incidents that occur in connection with the product must be reported to HEINE Optotechnik GmbH & Co. KG and the member state's competent authority.

Charge your device during longer storage periods to protect the battery from deep discharge.

Store and use the device only in a dry and dust-free environment.

General warnings

Use only CE marked USB power supplies (5 V) from reputable manufacturers that conform with the safety requirements from the IEC 60950-1 Information technology or IEC 60601-1 Medical electrical equipment.

Do not use power supplies where damage has been detected.

Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

This product is not allowed to be entered or used in strong magnetic fields like MRI scanners

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Do not look directly into the light source to avoid dazzle from the intense light.

The dermatoscopes are not suitable for eye examination.

Disposal

Dermatoscope

The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

GoCover

Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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Electromagnetic disturbances – Requirements and tests	
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.	
Statement for the operational environments	Inside professional healthcare facilities except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
	The supply voltage quality should be that of a typical hospital environment.
	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Performance features of the ME system that have been determined to be essential to the performance	None
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle	
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the considered device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
	Illumination failures may occur due to electrostatic discharge.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	EMC compatibility is only ensured if original HEINE spare parts, accessories and power sources are used as described in the chapter "Options". The EMC compatibility when using power sources from other manufacturers must be evaluated by the user.
Test	Compliance
RF emissions CISPR11	Group 1 Class B
Conducted emissions (EN 55011/CISPR 11)*	Passed
Radiated emissions (EN 55011/CISPR 11)	
Harmonic current emissions (IEC 61000-3-2)*	
Voltage changes, voltage fluctuations and flicker (IEC 61000-3-3)*	
Immunity	See attached immunity test levels

Immunity test levels		
Test	Test level	
	IEC 60601-1-2 test levels	Compliance test levels
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: ±8 kV Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV	
Radiated RF EM fields (IEC 61000-4-3)	3 V/m 80-2700 MHz 80 % AM at 1kHz	
Electrical fast transients / bursts (IEC 61000-4-4)*	± 2 kV 100 kHz repetition frequency	
Surges (IEC 61000-4-5)*	± 0,5 kV, ± 1 kV	
Conducted disturbances induced by RF fields (IEC 61000-4-6)*	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz, FM: ±5 Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m	
Power frequency magnetic fields (IEC 61000-4-8)	30 A/m; 50 Hz or 60 Hz	Test was not necessary, because there are no magnetic sensitive components.
Voltage dips (IEC 61000-4-11)*	0 % U _i ; 0,5 cycle; at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	Test was not performed, because it was not considered necessary. No risks were identified.
Short interruptions (IEC 61000-4-11)*	0 % U _i ; 250/300 cycle	

*n/a: "Not applicable" in the internally powered mode

Technical specification


Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Power supply	USB 2.0 Micro B: 5 V DC, 650 mA
Power consumption	max. 3.25 W
Class	Charging: II Operation: internally powered
Charging time	max. 75 min
Operating time	typ. 80 min at highest brightness in polarisation mode typ. 600 min at lowest brightness in non-polarisation mode
Device classification according to IEC 62471	Risk group 1
Applied part	Contact plate type BF

Options

USB cord with E4-USB plug-in power supply	X-000.99.303
Contact plate DELTAone	K-000.34.222
Small contact plate DELTAone	K-000.24.207
Mounting case smartphone	K-000.34.26x
Universal smartphone connector	K-000.34.270
GoCover	K-000.34.010

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