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HEINE otoscope G100 LED

ENGLISH

Please read and follow these instructions for use of and keep them for future reference.

Intended use

The HEINE Otoscope G100 LED is intended for the examination of the ear, the external auditory canal and the tympanum. The instrument has a swivel lens for optical magnification and an internal powered illumination unit. The device is intended for transient examination and can also be used for instrumentation with non-active medical devices. The product may only be used by qualified veterinarian and in professional healthcare facilities.

⚠️ For U.S. only:
Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

Warnings and Safety Information

⚠️ CAUTION! Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

📌 NOTE! This symbol indicates valuable advice. This notes are important, but no related to hazardous situations.

Product overview (see figure)

- | | |
|------------------------------|------------------------------------|
| 1 G100 LED | 5 Slit metal specula, open |
| 2 HEINE Handle | 6 Short, closed specula |
| 3 Swivel lens | 7 HEINE UniSpec Disposable Specula |
| 4 Long, closed metal specula | 8 Adaptor for UniSpec Specula |

Setting up

To set up the HEINE devices, screw the instrument head into the HEINE battery handle or plug it on the HEINE rechargeable handle.

Verify that the lamp voltage complies with the supply voltage of the power source. The coloured marking on the bottom of the lamp shows you the lamp voltage:

- White ring = HEINE 2,5 V bulb
- Only for use with the HEINE BETA Battery handle.
- Black ring = HEINE LED 3,5 V bulb
- Only for use with the HEINE BETA rechargeable handle and HEINE wall unit.

Operation

The otoscope may only be inserted into the auditory canal when a tip or specula is fitted. Use an appropriate size of tip or specula that suits the animal or the anatomical structures.

Place the HEINE reusable specula or specula in such a way onto the instrument head that the guide pin fits into the bayonet slit. Tighten the tip by twisting it slightly to the right.

HEINE UniSpec disposable specula can be attached by using an adaptor.

Swivel Lens

The instrument head is fitted with a swivel lens of approx. 1.6 x magnification. The swivel lens may be used with all specula, except the 170 mm long, closed metal specula, as this speculum may affect the focusing

Application duration

📌 The G100 LED is intended for transient examination < 6 min with a break between applications of 10 minutes.

The setting up and operation of the HEINE power supplies are described in a separate instruction of use.

Hygienic reprocessing

- The instruction is available:
 - at www.heine.com
 - in a paper version which you can request from the address listed

Maintenance

Changing the light source

⚠️ Verify that the lamp voltage complies with the supply voltage of the handle.

Allow the device to cool down before changing the bulb.

Remove the instrument head from the handle. Remove the lamp by pulling on the narrow ridge of the lamp. Insert the new lamp into the guide tube as far as it will go.

Service

The device has no components serviceable by the end-user.

General Notes

📌 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

The expected life cycle, when the device is normal used and the warning and safety information as well as the maintenance instructions are observed, is up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

Note to the user:

All serious incidents that occur in connection with the product must be reported to HEINE Optotechnik GmbH & Co. KG and the member state's competent authority.

General Warnings

⚠️ Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anesthetic environments).

This product is not allowed to be entered or used in strong magnetic fields like MRI scanners.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not shine light of the instrument directly into the eyes.

Choose a tip/specula size suitable for the animal.

Check the reusable specula and speculum before each use for rough surfaces and edges.

In order not to exceed the surface temperature of 41 °C of the application part, keep within the application duration.

Remove the instrument head from the handle only outside the environment of the animal (at least 1.5 meters from the animal or animal support pursuant to IEC 60601-1).

Disposal

- 🗑️** The product must be recycled as separated electrical and electronic devices.
- 📌** Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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- Distúrbios eletromagnéticos – Requisitos e testes

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- Esclarecimento sobre os símbolos utilizados

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

	Inside professional healthcare facilities except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Statement for the operational environments	The supply voltage quality should be that of a typical hospital environment.
	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

Performance features of the ME system that have been determined to be essential to the performance	
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle	None

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the considered device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

EMC compatibility is only ensured if original HEINE spare parts, accessories and power sources are used. The EMC compatibility when using power sources from other manufacturers must be evaluated by the user.

Test	Compliance
RF emissions CISPR11	Group 1 Class B
Conducted emissions (EN 55011/CISPR 11)*	
Radiated emissions (EN 55011/CISPR 11)	
Harmonic current emissions (IEC 61000-3-2)*	Passed
Voltage changes, voltage fluctuations and flicker (IEC 61000-3-3)*	
Immunity	See attached immunity test levels

Immunity test levels		
Test	Test level	
	IEC 60601-1-2 test levels	Compliance test levels
Electrostatic Discharge (IEC 61000-4-2)	Contact discharge: ±8 kV <p>Air discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV</p>	
Radiated RF EM fields (IEC 61000-4-3)	3 V/m <p>80–2700 MHz</p> <p>80 % AM at 1kHz.</p>	
Electrical fast transients bursts (IEC 61000-4-4)*	± 2 kV <p>100 kHz repetition frequency</p>	G100 LED with a 2,5 V bulb: not applicable, because of the internally powered mode.
Surges (IEC 61000-4-5)*	± 0.5 kV, ± 1 kV	Otoscope with a 3,5 V bulb: Compliance test levels are identically to the IEC 60601-1-2 test levels
Conducted distur-bances induced by RF fields (IEC 61000-4-6)*	3 V <p>0.15 MHz – 80 MHz</p> <p>6 V in ISM bands between 0.15 MHz and 80 MHz</p> <p>80 % AM at 1 kHz</p>	
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse modulation: 18 Hz; 27 V/m <p>450 MHz, FM: ±5 Hz deviation: 1 kHz sine; 28 V/m</p> <p>710, 745, 780 MHz; Pulse modulation: 217 Hz; 9 V/m</p> <p>810, 870, 930 MHz; Pulse modulation: 18 Hz; 28 V/m</p> <p>1720, 1845, 1970 MHz; Pulse modulation: 217 Hz; 28 V/m</p> <p>2450 MHz; Pulse modulation: 217 Hz; 28 V/m;</p> <p>5240, 5500, 5785 MHz; Pulse modulation: 217 Hz; 9 V/m</p>	
Power frequency magnetic fields (IEC 61000-4-8)	30 A/m; 50Hz or 60 Hz	
Voltage dips (IEC 61000-4-11)*	0 % U _r ; 0.5 cycle; at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° <p>0 % U_r; 1 cycle and 70 % U_r; 25 30 cycles single phase: at 0°</p>	G100 LED with a 2,5 V bulb: not applicable, because of the internally powered mode.
Short interruptions (IEC 61000-4-11)*	0 % U _r ; 250 300 cycle	G100 LED with a 3,5 V bulb: Compliance test levels are identi-cally to the IEC 60601-1-2 test levels

*n.a.: „Not applicable" in the internally powered mode

Technical specification	
Environmental conditions for operation	+10 °C to +35 °C <p>30 % to 75 % rel. humidity</p> <p>700 hPa to 1060 hPa</p>
Environmental conditions for storage	+5 °C to +45 °C <p>45 % to 80 % rel. humidity</p> <p>500 hPa to 1060 hPa</p>
Environmental conditions for transport	-20 °C to +50 °C <p>45 % to 80 % rel. humidity</p> <p>500 hPa to 1060 hPa</p>
Nominal voltage	3.0 V – 3.7 V
Nominal current	3,5 V: max. 350 mA <p>2,5 V: max. 450 mA</p>
Protection class	Internally powered
Applied part	Type BF
Device classification according to IEC 62471	Exempt
IP	40

Accessories	
Closed metal specula <p>57 mm long 4 mm dia.</p> <p>65 mm long 6 mm dia</p> <p>90 mm long 7 mm dia.</p> <p>170 mm long 10 mm dia.</p>	G-000.21.330 <p>G-000.21.331</p> <p>G-000.21.332</p> <p>G-000.21.350</p>
Slit metal specula, open <p>65 mm long 6 mm dia.</p> <p>90 mm long 7 mm dia.</p>	G-000.21.341 <p>G-000.21.342</p>
Short, closed specula <p>1 set = 5 piece</p> <p>40 mm long 2,2 dia.</p> <p>40 mm long 2,8 dia.</p> <p>40 mm long 3,5 dia.</p> <p>40 mm long 4,5 dia.</p> <p>40 mm long 5,5 dia.</p>	G-000.21.316 <p>G-000.21.310</p> <p>G-000.21.311</p> <p>G-000.21.312</p> <p>G-000.21.313</p> <p>G-000.21.314</p>
UniSpec disposable specula <p>1000 pcs. 2,5 mm dia.</p> <p>1000 pcs. 4 mm dia.</p>	B-000.11.242 <p>B-000.11.241</p>
Tip Adaptor	G-000.21.302

Spare parts	
Swivel lens	G-000.21.209

Spare bulb		
Type	2.5 V	3.5 V
G100 LED	X-028.88.132	X-008.88.133

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