Instruction for use

HEINE BETA X VET

Ophthalmoscope

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Please read and follow these instructions for use and keep them for future reference.

Intended use

BETA X VET Ophthalmoscope

The HEINE BETA X VET direct ophthalmoscope is intended for transient examination of the ocular media (cornea, aqueous humor, lens, vitreous humor) and the retina. The instrument features an optical examination system and an illumination unit powered by a rechargeable battery.

The product may only be used by qualified veterinarians in closed rooms. The examination is limited to small animals (e.g., dogs, cats, pets) and farm animals (e.g., cattle, pigs, sheep, goats, horses).

X CS1 VET Single and X CS1 VET Dual

The HEINE X CS1 VET is a charging station intended exclusively for recharging the internal power supply of the HEINE X-series instruments. The product may only be used indoors by qualified veterinarians or trained personnel authorized to operate such products.

E4- USB and E12-USB

The HEINE E4-USB / E12-USB is intended for charging exclusively with the HEINE USB power supplies. The product may only be used indoors by qualified veterinarians or trained personnel authorized to operate such products.

For U.S. only:

Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

Warnings and Safety Information

CAUTION! This symbol indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent.
(Background color yellow; foreground color black)

NOTE! This symbol indicates valuable advice. This notes are important, but no related to hazardous situations.

Product overview

- | | |
|--|---|
| 1 Spectacle rest | 8 Filter lever (polarisation filter, without filter, red-free filter) |
| 2 Viewing aperture | 9 On/Off switch and brightness control |
| 3 Lens wheel | 10 Charging compartment |
| 4 4 Diopter readout
(plus = green, minus = red) | 11 Charge status indicator of the instrument |
| 5 Switch: HEINE focusCONTROL | 12 USB-C socket for charging the instrument |
| 6 Viewing window | 13 USB-C socket on the bottom for setting up with the power supply |
| 7 Aperture wheel | |

Lens wheel (3) correction values:
in 1 diopter steps: -36dpt to +38dpt

Setting up

Before using the device for the first time, remove all protective films and fully charge the device.

Operation

Use the rotary push button (9) to switch the device on and off. You can also adjust the brightness by turning the control button. Turning the button quickly results in a large brightness adjustment. Turning it slowly, on the other hand, produces a fine adjustment.

The aperture wheel (7) provides you with 7 different apertures for your examinations, which can be used additional to the polarisation filter or the red-free filter.

The lens wheel (3) is used to adjust the refraction compensation. You can set this in the range of -36 diopters to +38 diopters. The respective correction values are visible via the display (4) (negative values are shown in red, positive values in green).

With the focusCONTROL toggle switch (5), you can easily switch between focus on the anterior segment of the eye (C - cornea) and fundus (F - fundus).


With the filter lever (8) you can switch on the polarisation filter or the red free filter.

Possible aperture settings (7):

From left to right: blue filter, fixation star with polar coordinates, large aperture, medium aperture, MicroSpot, slit aperture, hemispace

Possible filter settings (8):

From left to right: red free filter, without filter, polarisation filter

HEINE focusCONTROL (5):

Above: Setting for focussing on the fundus (F - fundus)


Below: Setting for focussing on the anterior segment of the eye (C - cornea)

Hard-Reset:

Should a reset be necessary, it can be reset by pressing the on | off button for a period of 10 seconds


Charging

Charge the device outside of the patient environment (1.5 m).

Charging via USB-C

The mains socket adapter of the power supply is interchangeable, this enables worldwide operation. The pictograms describe the procedure of changing the adapters:

Connect the USB cable to the USB-C port on the ophthalmoscope and to the XCS1 charger and HEINE power supply unit. Place the power supply in such a position where it is accessible and easily unplugged. The green LED on the E12-USB power plug indicates that the power supply is ready for operation.

After you checked the mains voltage, connect the switch mode power supply to the mains supply. When taking the instrument out of service, proceed in reverse order.

Charging via X CS1 VET Single/Dual

The HEINE X CS1 VET charger is intended to be placed on a table or a comparable stable surface. In addition, the charger provides an M4 thread on the bottom. This can be used to mount the charger if required.

To set up the charger, a 5 V power supply (not included in delivery) is required, which meets the necessary safety requirements from the IEC 60601-1 Medical electrical equipment. The HEINE E12-USBC power supply is recommended, which can be purchased separately.

For setting up, connect the power supply to the USB-C port on the bottom of the charger (13). Make sure that the USB-C cable is firmly connected to the charger.

To ensure all-pole disconnection from the mains at all times, the ME device must be installed so that the switch mode power supply is accessible and disconnectable.

To put the charger out of operation, disconnect the charger from the power supply.

For charging the X Series instrument, place it completely in the charging compartment (10) of the charger as pictured in the product overview. Make sure that the charge status indicator of the instrument (11) is visible. Once the instrument has been correctly inserted into the charger, the instrument is automatically charged.

Make sure that there are no objects in the charger.
Make sure that only hygienic reprocessed X series instruments will be inserted into the charger.

The charge status indicator on the X series instrument (11) displays that it is being charged. It also displays the charge status.

If the charge status indicator on the X Series instrument does not light up after placing the instrument in the charger, please check whether the power supply is plugged into the mains and the charger and that the instrument is correctly inserted into the charger. If necessary, reset the charging process by removing the instruments from the charger for a few seconds and then reinserting it. If you do not notice any improvement, please contact the manufacturer.

Battery indicator

Instrument	Charger	Battery	Charge indicator segments		
			Green - 1	Green -2	Orange - 3
switched off	connected	100 %	on	on	on
		Full (99 % to 66 %)	flashing	on	on
		sufficiently charged (66 % - 33 %)	off	flashing	on
		empty (≤ 33 %)	off	off	flashing
switched on	not connected	full (≥ 66 %)	on	on	on
		sufficiently charged (66 % - 33 %)	off	on	on
		empty (33 % - 10 %)	off	off	on
		Residual capacity (≤ 10 %)	off	off	flashing

Hygienic Reprocessing

BETA X VET Ophthalmoscope, XCS1 VET Single and X CS1 VET Dual

The instruction is available:

- at www.heine.com
- in a paper version which you can request from the address listed

E4-USB and E12-USB

Clean the cable and the power supply unit only if necessary and only with a dry cloth.

Maintenance

There is no regular maintenance required.

Service

The device has no components serviceable by the end-user.

General Notes

The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

The expected life cycle amounts with designated use and the observation of warning and safety information as well as the maintenance instructions up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

Charge your device during longer storage periods to protect the battery from deep discharge.

For the examination please use as little light as possible. Please ensure to switch off the light after each examination.

The instrument automatically switches off after 5 minutes if the brightness control is not operated.

The device only indicates the year of manufacture. The service life starts from the first of January of the year of manufacture.

Note to the user:

All serious incidents that occur in connection with the product must be reported to HEINE Optotechnik GmbH & Co. KG and the member state's competent authority.

General Warnings

Ophthalmoscope: Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or flashing light or unusually bright or dimmed lightning.

Charger: Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not place metallic objects (e.g. paper clips, coins) in the charging cradle.

Do not use the device in fire- or explosive risk area.

Do not modify the device.

Do not open the enclosure.

This product is not allowed to be entered or used in strong magnetic fields like MRI scanners.

Repairs shall only be carried out by qualified persons.

Use only CE marked USB power supplies (5 V) from reputable manufacturers that conform with the safety requirements from the IEC 60601-1 Medical electrical equipment.

Do not use power supplies or cables if you notice any damage.

Do not use the power supply outdoors.

Only use the power supply in a dry environment. Avoid splashing water. Never immerse in liquids.

Use only original HEINE parts, spare parts, accessories and power sources.

Lay cables in a protected position and ensure that no-one can fall over them.

Heating during operation is normal and harmless. Make sure that the device isn't covered during operation, that it isn't placed in close proximity to a heating and that isn't exposed to direct sunlight.

The device must not be used with appliance outlets of extension leads.

The power supply is only approved for the use with ungrounded loads (Class II).

Before operating the device, verify that the mains voltage corresponds to the value indicated on the power supply.

Use this power supply only for the supply of the devices it is intended for.

Warning: Danger to life due to electric shock:

Only connect the power supply unit to an appropriate socket.

Never carry or pull the power supply unit by the electrical cable.

Only pull on the plug to disconnect a plug connection, never on the cable.

Lay cables in a protected manner and ensure that nobody can fall over them.

Hazard from Light

Prolonged intense light exposure can damage the retina, so the use of the device for eye examination should not be unnecessarily extended, and the brightness should not be set higher than necessary for clear observation of the target structure. The exposure dose for photochemical retinal hazard is the product of irradiance and exposure time. If the irradiance is reduced by half, the exposure time may be twice as long to reach the maximum limit.

While no acute optical radiation hazards have been identified for direct or indirect ophthalmoscopes, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

ISO 15004-2:2007: Group 2

Warning – The light from this instrument may be harmful. The risk of eye damage increases with exposure time. Exposure with this instrument at maximum intensity for longer than 2 minutes 26 seconds exceeds the recommended threshold for hazard.

Disposal

The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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Electromagnetic disturbances – Requirements and tests	
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such environments.	
Statement for the operational environments	Inside professional healthcare facilities except for near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
	The supply voltage quality should be that of a typical hospital environment
	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Performance features of the ME system that have been determined to be essential to the performance	None
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle	
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	EMC compatibility is only ensured if original HEINE spare parts, accessories and power sources are used. The EMC compatibility when using power sources from other manufacturers must be evaluated by the user.
Test	Compliance
RF emissions CISPR 11	Group 1 Class B
Harmonic Emissions	Passed*
Voltage fluctuations/Flickers	
Immunity	Passed*

*n/a: „Not applicable“ in the internally powered mode

Technical specifications

BETA X VET Ophthalmoscope

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+ 5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	- 20 °C to +50 °C 10 % to 90 % rel. humidity 500 hPa to 1060 hPa
Nominal voltage	5 V
Nominal current	500 mA
Protection class	Internally powered Class II (charging)
Device classification according to IEC 62471	Group 2
Cable	USB-cable max. 2 m

X CS1 VET Single and X CS1 VET Dual

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+ 5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	- 20 °C to +50 °C 10 % to 90 % rel. humidity 500 hPa to 1060 hPa
USB-C (input)	5 V; max. 1.2 A
USB-C (output)	5 V DC; max. 0.5 A per each output
Power consumption	≤ 15 W
Protection class	II (through external power supply)

E4-USB

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+ 5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	- 20 °C to +50 °C 10 % to 90 % rel. humidity 500 hPa to 1060 hPa
Nominal input voltage	100 – 240 V AC / 50 – 60 Hz
Power consumption	< 15 W
Output voltage	5 V
Output current	max. 1.2 A
Protection class	II
Protection provided by enclosures (IP code)	IP20

E12-USB

Environmental conditions for operation	0 °C to +40 °C 0 % to 95 % rel. humidity up to 5000 m
Environmental conditions for storage	- 40 °C to +80 °C 0 % to 95 % rel. humidity up to 5000 m
Environmental conditions for transport	- 40 °C to +80 °C 0 % to 95 % rel. humidity up to 5000 m
Nominal input voltage	100 – 240 V AC / 50 – 60 Hz / 0.45 A
Power consumption	max. 15 W
Output voltage	5 V
Output current	3 A
Protection class	II
Protection provided by enclosures (IP code)	IP20

Accessories

The corresponding accessories are listed in the following table:

USB cord with E4-USB plug-in power supply	X-000.99.300
X CS1 Single	X-002.99.215
X CS1 Dual	X-002.99.216

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