

HEINE Cube Quick Reference Guide

1. Connect your HEINE Cube using the supplied power cable. The HEINE Cube is started by pressing the ON button.



2. Connect your device to the HEINE Cube WLAN, which starts automatically after approx. two minutes.



SSID: HEINE-Cube
Password: heinecube

3a. On your
Mac / Windows / Tablet



open one of the
commonly used browsers



open the HEINE Cube address
192.168.88.1

Log in

3b. You have an iPhone* version 5s
or higher



download the HEINE Cube App



open the app and enter the address
192.168.88.1

Log in

LOGIN

Please login to unlock the application

User

Password

Login!

