## HEINE Otoscopes

HEINE Optotechnik GmbH & Co. KG  
Dornierstr. 6 · 82205 Gilching · Germany  
E-Mail: info@heine.com · www.heine.com

MED 112680 2020-06-02

V200.001.51

## DEUTSCH

## HEINE Otoscopes

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## HEINE Otoscopes

These instructions apply to the following products of the HEINE Otoscope series:  
BETA200 LED, BETA200, BETA400 LED, BETA400, K180, mini3000 F.O., mini3000 F.O. LED, BETA100, K100, mini3000.

**⚠️** Please read and follow these instructions for use and keep them for future reference.

## Intended Use

HEINE otoscopes are intended for the examination of the ear, the external auditory canal and the tympanum. The instruments have optical magnifying equipment and battery or rechargeable battery powered illumination equipment. The device may only be used by qualified medical personnel and in professional healthcare facilities. Patients can be examined as often as necessary without any further restrictions. The devices are intended for short term treatment and can also be used for instrumentation with non-active medical devices.

## For U.S. only:


Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

## Warnings and Safety Information

**⚠️** CAUTION! Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

**NOTE!** This symbol indicates valuable advice in terms of set up, operation, maintenance or repair, as applicable. Notes are important, but not related to hazardous situations.

## Product overview

## Setting up

To set up the HEINE devices, screw the instrument head into the HEINE battery handle or plug it on the HEINE rechargeable handle.

Verify that the lamp voltage complies with the supply voltage of the handle. The coloured marking on the bottom of the lamp shows you the lamp voltage:

White ring = HEINE XHL® 2.5 V bulb

Red ring = HEINE XHL® 3.5 V bulb

Black ring = HEINE LED illumination

only for use with the HEINE BETA® rechargeable handle and HEINE® wall unit

This does not apply to the otoscope of the mini3000 series. Their lamp voltage is always 2.5 V independently of the handle.

## Operation

The HEINE otoscope may only be inserted into the auditory canal when a tip is fitted.

Place the HEINE tip on the otoscope so that the projection inside the tip fits into the bayonet slot (1). Tighten the tip by twisting it slightly to the right.

The magnifying lens is inside the viewing window (2). To facilitate instrumentation, the viewing window can be swivelled to both sides or folded up.

The otoscopes have a connection port (3) for the optionally available insufflation bulb (6,7) (with the exception of the mini3000).

The test of tympanic mobility can be carried out when the viewing window is closed using the insufflation bulb. Apply the pressure carefully with the insufflation bulb.

If, after switching on the unit, a sharp drop in brightness or a flashing of the light is noticed, new batteries should be inserted or rechargeable batteries recharged.

## Application duration

The otoscopes are intended for transient examination < 1 min with a break between applications of 10 minutes.

The maximum temperature at the application point at an ambient temperature of 35 °C and in continuous operation is 44,5 °C.

The setting up and operation of the HEINE handles are described in a separate instruction of use.

## Hygienic Reprocessing

The instruction is available:

- at [www.heine.com](http://www.heine.com)
- in a paper version which you can request from the address listed

## Maintenance

## Changing the light source

**⚠️** Verify that the lamp voltage complies with the supply voltage of the handle.

Allow the device to cool down before changing the bulb.

## BETA200 LED / BETA400 LED / mini3000 F.O. LED

With these otoscopes, the LED cannot be changed.

## BETA200 / BETA400 / K180 / mini3000 F.O.

- Detach the otoscope from the handle and pull the bulb out of the guide tube of the connector plug (4) holding it at the narrow collar.
- Insert the new lamp all the way into the lamp guide.

## BETA100 / K100

- Remove the tip from the otoscope and unscrew the lamp. Remove the sleeve from the oil lamp.
- Refit the new bulb into the socket.
- Slide the sleeve on as far as it will go. This reduces stray light and reflections from the inside of the specula and should always be fitted.

## min3000

With this otoscope, swivel the viewing window (2) to the side and pull out the bulb backwards. Insert the new bulb as far as possible into the socket.

## Service

The device has no components serviceable by the end-user.

## General Notes

**⚠️** The warranty for the entire product is invalidated if non-original HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

The expected life cycle, when the device is normal used and the warning and safety information as well as the maintenance instructions are observed, is up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

## General Warnings

**⚠️** Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

This product is not allowed to enter or be used in areas with strong magnetic fields e.g. MRI scanners.

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