**HEINE** QUALITY  
MADE IN GERMANY

## ∴ HEINE iC1

CE

**HEINE Optotechnik GmbH & Co. KG**  
Dornierstr. 6 · 82205 Gilching · Germany  
E-Mail: [info@heine.com](mailto:info@heine.com) · [www.heine.com](http://www.heine.com)

MED 112655      2020-07-14

V-200.00.037

# HEINE iC1

<b>DEUTSCH</b>	<b>3</b>
<b>ENGLISH</b>	<b>7</b>
<b>FRANÇAIS</b>	<b>11</b>
<b>ESPAÑOL</b>	<b>15</b>
<b>ITALIANO</b>	<b>19</b>
<b>SVENSKA</b>	<b>23</b>
<b>NEDERLANDS</b>	<b>27</b>
<b>DANSK</b>	<b>31</b>
<b>SUOMI</b>	<b>35</b>
<b>PORTUGUÊS</b>	<b>39</b>

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# HEINE iC1

Please read and follow these instructions for use of and keep them for future reference.

## Intended use

### *iC1 dermatoscope*

The iC1 dermatoscope is an internally powered medical examination light. It is an reflected-light microscope for non-invasive, visual inspection of intact skin by a healthcare professional. The diagnostic is not limited to any patient demography.

### *iC1 dermatoscope with mounting case smartphone option*

The optional mounting case for the iC1 dermatoscope enables an integration of a smartphone which optimizes the macrophotography. The digitized image display takes place within the technical capabilities of the smartphone and in its range of application and is solely for documentation and training purposes.

## For U.S. only:

 Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

## Warnings and Safety Information

 **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** Note indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

## Product overview

- 1 Contact plate with scale
- 2 Filter insert
- 3 Toggle switch
- 4 Micro USB-charging port

- 5 Charge indicator
- 6 On/Off switch
- 7 Mounting case smartphone

## Setting up

### *Assembly of filter insert and contact plate*

Attach the insert (2) to the contact plate (1) and rotate the insert to lock the bayonet connection. To disassemble the filter (2) from the contact plate, reverse the procedure.

### *Attachment of contact plate*

The contact plates (1) are attached by a bayonet connection. To remove the contact plate turn it counterclockwise and pull it away from the dermatoscope. To attach, reverse the procedure. Always check whether the bayonet is safely locked.

## Setting up smartphone

Insert the smartphone into the appropriate mounting casing so it is firmly seated. Make sure that the camera is in the correct position.

To avoid any possible malfunctions, please deactivate all services (flight mode).

Images for documentation purposes can be taken with the standard camera app on your smartphone. Also, please deactivate the internal flash of the phone.

## Operation

The iC1 allows for a rapid change from polarized to non-polarized examination mode over a toggle switch. The difference is that no immersion fluid is required in the polarized examination mode. The default setting is the polarized mode.

Turn the dermatoscope on by using the on/off switch (6). Gently place the instrument so that the lesion is in the center of the contact plate.

### *For non polarized mode*

Prepare the skin by moistening with HEINE dermatoscopy-oil (use a cotton swab) or disinfectant spray.

 Always use the device in combination with the filter inserts. Only use contact plates from HEINE.

Switch off the device after each examination.

### *Brightness control*

The electronic voltage regulation in the dermatoscope guarantees a constant brightness.

 An electrical conducting connection between smartphone, PC and a mains power source is not permissible during operation.

## Charging

### *Charge indicator*

First insert the power supply into the wall socket or PC and then the micro USB into the charging port (4).

When the iC1 is charged via USB, the LED will flash green.

When the battery is fully charged, the LED will be constant green.

When the battery reaches a critical charge state, the LED will light up orange. The device must now be recharged.

Charge the device outside of the patient environment (1.5 m).

## Hygienic reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

Classification according to KRINKO: non-critical

Spaulding Classification USA: non-critical

 In the event of suspected contamination, carry out hygienic preparation of the instrument.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the resources and procedures named in these instructions for use.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the instructions of the manufacturer of the reprocessing media.

Do not use spray or immersion disinfection, dripping wet or heavily foaming tissues.

Do not use ultrasonic reprocessing and don't sterilize.

Do not use reprocessing media including alcohol.

The contact plates have to be cleaned and/or disinfected after each use.

## Procedure

Prior to hygienic reprocessing, turn the power off, remove the USB cable and extract the smartphone.

Ensure that moisture does not penetrate the area of the USB connector.

## Dermatoscope

Clean and disinfect the dermatoscope manually (clean and disinfect through wiping).

### Recommended agents

Cleaning agent: Neodisher MediClean

Disinfectant agent: quaternary ammonium compounds (e.g. Cleanisept Wipes or Microbac Tissues)

## Contact plate

Clean and disinfect the contact plate manually after removing it from the instrument (clean and disinfect through wiping).

Before cleaning or disinfection you must remove the filter insert.

### Recommended agents

Cleaning agent: Neodisher MediClean

Disinfectant agent: quaternary ammonium compounds (e.g. Microbac Tissues)

## Changing the light source

The LEDs cannot be changed.

## Maintenance and service

The instrument does not require maintenance or service.

## General warnings

Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not look directly into the light source to avoid dazzle from the intense light. The dermatoscope is not suitable for eye examination.

## General notes

The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE.

For more information, please visit [www.heine.com](http://www.heine.com).

## Disposal

The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

## Electromagnetic compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

 The iC1 dermatoscope, in connection with a smartphone, should not be used within the vicinity of life-sustaining medical devices.

This is a device in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site.

The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission and reduced electrical immunity of the medical equipment.

The device (iC1 + smartphone) may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

The appendix contains following tables

- Guidance and manufacturer's declaration – Electromagnetic immunity
- Technical specification
- Explanation of the used symbols

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### Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.

Emission test	Compliance	Electromagnetic environment – Guidelines
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	Symmetrical three-phase devices and other devices.
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Passed	

### Guidance and manufacturer declaration - Electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guideline
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	± 2 kV for mains cables ± 1 kV for input and output lines	The supply voltage quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT, (>95% dip in UT) for 1/2 period 40% UT, (60% dip in UT) for 5 periods 70% UT, (30% dip in UT) for 25 periods <5% UT, (>95% dip in UT) for 5 seconds	< 5% UT, (>95% dip in UT) for 1/2 period 40% UT, (60% dip in UT) for 5 periods 70% UT, (30% dip in UT) for 25 periods <5% UT, (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	10 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Comment: UT is the a.c. supply voltage prior to application of the test level.			

### Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	10 V eff	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	10 V/m	<p><b>Recommended separation distance:</b>  <math>d = 3,5/3 * \text{SQRT} (P/W)</math>  <math>d = 3,5/3 * \text{SQRT} (P/W)</math> 80 MHz to 800 MHz  <math>d = 7/3 * \text{SQRT} (P/W)</math> 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).                      Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup>                      Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

Note 1: At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

### Recommended separation distances for portable and mobile RF communication equipment and the device

The device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3,5/3 * \text{SQRT} (P)$	80 MHz to 800 MHz $d = 3,5/3 * \text{SQRT} (P)$	800 MHz to 2,5 GHz $d = 7/3 * \text{SQRT} (P)$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Technical specification

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Nominal charging voltage	5 V
Nominal current	max. 500 mA in charge mode
Protection class	internal power supply
Device classification according to IEC 62471	Group 2
Applied part	Type BF (for contact plate)
Operational time	typ. 3.5 h
Rechargeable battery charging time	typ. 2.5 h

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### Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

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