



# The new HEINE DERManager.

## One piece of software. Two operating options.

In simple terms, the dermatoscope and iPhone\* are used to capture the images, and the computer is used for the observation in detail.

Both devices can naturally be used to do much more, but here is an overview:



[Download the DERManager app](#)

[Access the DERManager web](#)

Function/Feature	DERManager app	DERManager web
Access patients	✓	✓
Create new patients	✓	✓
Manage work lists	✓	✓
Capture dermatological images	✓	
Capture series images	✓	
Capture overview images	✓	
Placement on bodymap	✓	
Placement on overview image	✓	
Checkup assistant	✓	
Evaluation of lesion + comment	✓	✓
Enlarged image		✓
Image comparison (1:1, 1:2, 1:4)		✓
Add more users / remove users		✓
Change / delete patient data		✓
Share patient data		✓
Delete images		✓
Upload locally stored images		✓


**HEINE** Digital documentation from the inventor of the dermatoscope.  
**DERManager**

Most of the functions and options that the HEINE DERManager offers are self-explanatory. For answers to almost all questions, "how to" videos and lots more, check out the [HEINE Knowledge Base](#).

\* Apple iPhone from version 5s upwards or iPod touch 6th generation. Apple, the Apple logo, iPhone and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

