



**If there's something  
to discover, you'll be  
sure to see it.**

# You could buy it because of its looks. Or because of its optics.

Reliable diagnostics thanks to a crystal-clear image, no matter whether you prefer to work close-up or further away. The HEINE DELTA 30 would not only win any beauty pageant; it stands out due to a very good imaging and a true field of view of 30 millimetre. So even the details of the details stand out.

Extremely sharp image  
with the achromatic  
HEINE optical system

Large field of view with  
real 30mm (lens 32mm) Ø

Fast change from contact  
to non-contact examination

Excellent colour  
rendering due to LED<sup>HO</sup>

Ergonomic handling due  
to angled product design



Digital documentation  
with the exclusive  
HEINE DERManager

The HEINE diSTANCE working  
ring enables the user to examine  
and manage a lesion under  
magnification without contact



[01]



**DELTA 30** with contact plate with scale, USB cord with  
medical approved plug-in power supply, case

**K-230.28.305**

**Charger 30** table charging station for DELTA 30

**X-002.99.212**

**HEINE diSTANCE** working ring [01]

**K-000.34.103**

Accessories and further articles under [heine.com/en/DELTA30](http://heine.com/en/DELTA30)

**HEINE Optotechnik GmbH & Co. KG**  
Dornierstr. 6, 82205 Gilching, Germany  
[www.heine.com](http://www.heine.com)



We reserve the right to change  
specification without notice.

10/23. A-000.00.228 EN

