

**Mobility and quality have now
become *one*.
The HEINE DELTA*one* Dermatoscope.**



Small on the outside. Big on the inside.

The HEINE DELTAone Dermatoscope.

Extremely sharp image with the achromatic HEINE optical system

Fast change from contact to non-contact examination

Fits into any pocket



Polarisation eliminates reflections and glare

LEDHQ
LED IN HEINE QUALITY.

Excellent colour rendering due to LEDHQ

Digital documentation with the exclusive HEINE DERManager*



- ∴ **Crystal clear image**, no matter if you prefer a closer or greater working distance
- ∴ **True field of view of 21 mm (lens 22 mm) Ø**
- ∴ **Brightness control**: 3 levels of light intensity
- ∴ **Easily cleaned with spray disinfectant**
- ∴ **10x magnification**
- ∴ **No cross-contamination**: autoclavable contact plate with scale
- ∴ **Dioptre compensation**
- ∴ **Switch between polarised and non-polarised illumination mode at the touch of a button** = toggle function. For more information, please visit heine.com/en/Toggle

- ∴ **Rechargeable battery with Li-ion technology**: no memory effect; safely recharge the battery regardless of its state of charge
- ∴ **Small contact plate** for difficult-to-access pigmented lesions
- ∴ **Possibilities of digital documentation**:
 - Image and data management convenient and intuitive with the HEINE DERManager.* More information on our website: www.heine.com/en/digital
 - Conventional digital images with many popular smartphones

DELTAone Standard Edition (black-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case

K-210.28.305

DELTAone White Edition (white-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case

K-211.28.305

∴ **Accessories and further articles under heine.com/en/DELTAone**

* Check on our website in which countries the DERManager is currently available.