

HEINE DELTAone

The most important points at a glance

Real 21 mm Ø
field of view

Switch between
polarised and
non-polarised light at
the touch of a button =
toggle function. For more
information, please visit
heine.com/en/Toggle


3 light intensity
levels

On/Off



HEINE DELTAone

The most important points at a glance



A film says more than a thousand words



The "Application" video shows the DELTAone in action.



Remember each and every mole.

The HEINE DERManager* is the ideal complement to your DELTAone. This software solution makes digital management of dermatoscopic image data simple.

*Check on our website in which countries the DERManager is currently available.

