


Straight Laryngeal Mirror



- ∴ A range of accessories make this instrument particularly versatile.
- ∴ **Enhanced XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- ∴ **Solid metal construction.** Durable and tough.

Straight Laryngeal Mirror	2.5V	3.5V
complete with plane mirror 20 mm dia., without handle	B-001.12.101	B-002.12.101
without mirror, without handle	B-001.12.109	B-002.12.109
XHL Xenon Halogen spare bulb	X-001.88.057	X-002.88.063
Plane mirror 20 mm dia.		B-000.12.103
Concave mirror 20 mm dia. (magnifying)		B-000.12.104
Metal tongue-depressor [01]		B-000.12.105
Transilluminator [02]		B-000.12.106
Protective cap for bulb (for mouth illumination) [03]		B-000.12.107

Curved Laryngeal Mirror



Assists examination with improved view of the larynx.

Curved Laryngeal Mirror	2.5V	3.5V
complete with plane mirror 20 mm dia., without handle	B-001.12.201	B-002.12.201
without mirror, without handle	B-001.12.208	B-002.12.208
XHL Xenon Halogen spare bulb	X-001.88.037	X-002.88.062
Plane mirror 20 mm dia.		B-000.12.203
Concave mirror 20 mm dia. (magnifying)		B-000.12.204
Metal tongue-depressor [01]		B-000.12.205

HEINE mini 3000® Laryngeal Mirror



Can only be used with the mini 3000 handle system.

- ∴ A range of accessories make this instrument particularly versatile (see above).
- ∴ **Enhanced XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- ∴ **Solid metal construction.** Durable and tough.

mini 3000 laryngeal mirror	2.5V
with plane mirror 20 mm dia, without handle	D-001.77.101
XHL Xenon Halogen spare bulb	X-001.88.057

Tongue-Blade Holder

For wooden tongue-blades



- **XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- **Light-guide for bulb.** Avoids dazzle.

Tongue-Blade Holder	2.5V	3.5V
for wooden tongue-blades, without handle	B-001.12.302	B-002.12.302
XHL Xenon Halogen spare bulb	X-001.88.037	X-002.88.049
Replacement light-guide for bulb [01]		B-000.12.308

Tongue-Blade Holder

With Blade Ejector for use with HEINE Disposable Tongue-Blades



- **XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- **Light-guide for bulb.** Avoids dazzle.
- **Blade ejector.** No need to touch used blades.
- **Only for HEINE Disposable Tongue-Blades.**
- **Disposable plastic tongue blades.** Strong and hygienic.

Tongue-Blade Holder	2.5V	3.5V
with blade ejector, without handle	B-001.12.305	B-002.12.305
XHL Xenon Halogen spare bulb	X-001.88.037	X-002.88.049
Pack of 100 disposable plastic tongue-blades		B-000.12.304

HEINE mini 3000® Tongue-Blade Holder

With blade ejector for HEINE Disposable Tongue-Blades



Can only be used with the mini 3000 handle system.
Available in black or blue.

- **Solid metal construction.** Durable and tough.
- **Finger-tip pressure releases the used blade.** No need to touch the blade after use.
- **Enhanced XHL Xenon Halogen Technology.** Very bright white light. Concentrated illumination of oral cavity and pharynx.
- **2-piece handle and head.** Easy to maintain, flexible. Compatible with other mini 3000 instrument heads.
- **High-quality handle: Chrome-finish upper section/refined plastic.** Shockproof, sturdy, non-slip.
- **Attachment clip with integrated on/off switch.** Secure. Switches off automatically when replaced in the pocket. 20,000 switch cycles guaranteed.
- **Only for HEINE Disposable Tongue-Blades.** Strong and hygienic.
- **Replaceable batteries.** Size AA, or optional rechargeable battery with mini NT charger*.

mini 3000 Tongue-Blade Holder	2.5V
complete with mini 3000 battery handle and 5 disposable blades, with batteries	D-001.74.118
Tongue-blade holder head with bulb, without handle	D-001.74.100



XHL Xenon Halogen spare bulb	X-001.88.037
Pack of 100 disposable plastic tongue-blades	B-000.12.304



*To upgrade to a rechargeable instrument with mini NT charger, please see page 148.



HEINE mini 3000® pocket instruments – a complete range of quality diagnostic instruments for your pocket. Available in black or blue. Please specify the colour you prefer when ordering. If not otherwise specified, we will supply black instruments.

Multi-purpose Illuminator



- **XHL Xenon Halogen Technology.** Bright, concentrated light.
- **Wide range of applications.** Used as a diagnostic light or tongue-depressor for wooden or disposable blades.

Multi-purpose Illuminator	2.5V	3.5V
Multi-purpose Illuminator	B-001.12.320	B-002.12.320
XHL Xenon Halogen spare bulb	X-001.88.077	X-002.88.078

Fiber Optic Nasal Illuminator



Accessory for Multi-purpose Illuminator (see above).

- **Fiber Optic Nasal Illuminator.** For specula and spreadable nasal speculum.

Fiber Optic Nasal Illuminator	B-000.12.321
-------------------------------	---------------------

Fiber Optic Nasal Illuminator



Fiber Optic Nasal Illuminator	2.5V	3.5V
complete with nasal speculum	B-001.12.322	B-002.12.322
same, with spreadable nasal speculum	B-001.12.323	B-002.12.323
XHL Xenon Halogen spare bulb	X-001.88.077	X-002.88.078
Nasal speculum for above [01]		B-000.11.220
Spreadable nasal speculum for above [02]		B-000.11.231

HEINE ENT Diagnostic Sets



Set complete with: **BETA 100 Diagnostic Otolaryngoscope**

- 3 reusable specula (B-000.11.222)
- tongue-depressor for disposable tongue-blades
- straight laryngeal mirror without accessories
- spreadable nose speculum
- one spare bulb each
- hard case

ENT Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V
BETA battery handle	B-244.10.118		
BETA R rechargeable handle for mains socket		B-244.20.376	B-244.27.376
BETA TR rechargeable handle with plug-in transformer		B-244.20.384	B-244.27.384
BETA NT rechargeable handle and NT 300 charger		B-244.20.420	
BETA L rechargeable handle and NT 300 charger			B-244.29.420



Set complete with: **K100 Diagnostic Otolaryngoscope**

- 3 reusable specula (B-000.11.222)
- tongue-depressor for wooden tongue-blades
- one spare bulb for each instrument
- hard case


ENT Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V
BETA battery handle	B-241.10.118		
BETA R rechargeable handle for mains socket		B-241.20.376	B-241.27.376
BETA TR rechargeable handle with plug-in transformer		B-241.20.384	B-241.27.384
BETA NT rechargeable handle and NT 300 charger		B-241.20.420	
BETA L rechargeable handle and NT 300 charger			B-241.29.420

: Rechargeable handle with charge status display in the bottom insert.

A company makes history and gives a name to outstanding product quality: HEINE.

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

www.heine.com



■ HEINE SUBSIDIARY COMPANIES
■ HEINE REPRESENTATIVES

GERMANY

HEINE Optotechnik GmbH & Co. KG
82211 Herrsching
Kientalstr. 7
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com

NORTH AMERICA

HEINE USA LTD.
10 Innovation Way
Dover, NH 03820
Tel. (603) 7 42-71 03
Fax (603) 7 42-72 17
Toll Free (800) 367-4872
E-Mail: service@heine-na.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
P.O. Box 717 Brookvale 2100
New South Wales
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobelackerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: heine-schweiz@schaffhausen.ch

Your distributor

10/11. A-000.00.100 e

