


**NEW! HEINE GAMMA® STETHOSCOPIES**

- ∴ **Solid chest piece** and ideal weight avoid stray noise from manipulation and ensure excellent resonance. Hygienic matt chrome finish.
- ∴ **Chrome-plated binaurals.** Individually-adjustable.
- ∴ **High-quality Y- piece tubing** for perfect acoustics and sound reproduction. Length 56 cm.
- ∴ **Tubes:** Latex-free.
- ∴ **5 year warranty** for materials and workmanship.

**HEINE GAMMA 3.1® Pulse Stethoscope**
**For fine acoustic performance.**


- ∴ **Stethoscope with solid flat chest piece.** The chest piece can be easily inserted under the sphygmomanometer cuff.
- ∴ **High-quality membrane,** 44mm.
- ∴ **Single-tube construction** with excellent sound transmission.
- ∴ **Includes ear olives** (small/hard) for sound insulation.
- ∴ **Optional available:** Ear olives, large/soft (M-000.09.945).

GAMMA 3.1

M-000.09.941

**HEINE GAMMA 3.2® Acoustic Stethoscope**
**For excellent acoustic performance.**


Technical data as GAMMA 3.1, but with a double chest piece.

- ∴ **Combined double chest piece** with two different frequency ranges to choose from: A high-quality membrane on one side and a conical acoustic head on the other side.
- ∴ **High-quality membrane,** 44mm diameter for excellent acoustic performance and sound transmission.
- ∴ **Cone diameter 30mm.**

GAMMA 3.2

M-000.09.942

**HEINE GAMMA 3.3® Acoustic Stethoscope**
**For the examination of children.**


Technical data as GAMMA 3.2.

- ∴ **Reduced contact area** of 31 mm makes GAMMA 3.3 ideal for examining children.

GAMMA 3.3

M-000.09.943

**HEINE GAMMA C3® Cardio Stethoscope**
**Specially-developed for cardiology.**


- ∴ **Exceptional acoustic performance,** particularly in the cardiac frequency range.
- ∴ **Combined double chest-piece** with two different frequency ranges. Instant selection of a double membrane on one side or a conical acoustic head on the other.
- ∴ **High-performance double membrane,** 45 mm diameter gives ideal sound transmission and perfect skin contact for exceptional sound reproduction.
- ∴ **Two-in-one tube system** for enhanced performance. Very supple, yet kink-resistant with exceptional acoustic reproduction.
- ∴ **Includes ear olives** (large/ soft) and one pair of replacements (small/ hard) for sound insulation.


GAMMA C3

M-000.09.944

**A company makes history and gives a name to  
outstanding product quality: HEINE.**

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

[www.heine.com](http://www.heine.com)



■ HEINE SUBSIDIARY COMPANIES  
■ HEINE REPRESENTATIVES

**GERMANY**

HEINE Optotechnik GmbH & Co. KG  
82211 Herrsching  
Kientalstr. 7  
Tel. +49 (0) 81 52-38 0  
Fax +49 (0) 81 52-3 82 02  
E-Mail: [info@heine.com](mailto:info@heine.com)

**NORTH AMERICA**

HEINE USA LTD.  
10 Innovation Way  
Dover, NH 03820  
Tel. (603) 7 42-71 03  
Fax (603) 7 42-72 17  
Toll Free (800) 367-4872  
E-Mail: [service@heine-na.com](mailto:service@heine-na.com)

**AUSTRALIA**

HEINE AUSTRALIA PTY. LTD.  
Unit 9, 98 Old Pittwater Road  
P.O. Box 717 Brookvale 2100  
New South Wales  
Tel. +61 (0) 2-99 38 95 00  
Fax +61 (0) 2-99 39 23 05  
E-Mail: [info@heine.com.au](mailto:info@heine.com.au)

**SWITZERLAND**

HEINE (Schweiz) AG  
Tobelackerstr. 9  
CH-8212 Neuhausen am Rheinfall  
Tel. +41 (0) 52-6 72 22 66  
Fax +41 (0) 52-6 72 63 77  
E-Mail: [heine-schweiz@schaffhausen.ch](mailto:heine-schweiz@schaffhausen.ch)

Your distributor

10/11. A-000.00.100 e

