

HEINE HSL 150® Hand-held Slit Lamp

For the examination of the anterior segment of the eye



∴ Continuous slit adjustment

∴ Blue filter

∴ Ideal for domiciliary visits

Compact and lightweight, the portable HSL 150 is ideal for situations where a full size slit lamp is not available or convenient. For the examination of the anterior segment in Human and Veterinary medicine.

- ∴ **Slit size adjustable** from 0.2 mm x 10 mm to 4 mm x 14 mm.
- ∴ **6x magnification.**
- ∴ **Blue interference filter (FITC)**, can be selected for corneal examination.
- ∴ **Multi-coated optics for maximum light transmission.**
- ∴ **2.5V or 3.5V XHL Xenon Halogen Technology** for bright, white light comparable with the brightness of a classic slit lamp.
- ∴ **The HSL 150 instrument head weighs just 70g**, the lightest slit lamp of its kind.

HSL 150 for BETA or BETA SLIM handles	AV	TL
HSL 150, head only	C-002.14.602	C-002.14.602 TL
XHL Xenon Halogen Spare Bulbs	X-002.88.099	X-002.88.099
Replacement filter attachment for HSL 150		C-000.14.605

HEINE HSL 10x® Loupe Attachment

Loupe attachment with 10x magnification



For HSL 150. With precision magnification optics for crisp, clear images.

- ∴ **10x magnification** when combined with HSL 150 Loupe.
- ∴ **Increased working distance** from examiner to patient.
- ∴ **Large field of view.**
- ∴ **With an additional detachable eyepiece** to steady the instrument on the brow and avoid stray light. Can be replaced with a soft cap, suitable for spectacle wearers.

HSL 10x Loupe Attachment	C-000.14.606
--------------------------	--------------

HEINE HSL 150® Hand-held Slit Lamp Sets






Set complete with: **HSL 150**
BETA SLIM handle (compact)


Hand-held Slit Lamp Sets	3.5V
BETA SLIM rechargeable handle and NT 300 charger, spare bulb, in hard case	C-266.20.470
same, with HSL 10x Loupe Attachment, in hard case	C-268.20.470



Set complete with: **HSL 150**,
BETA handle (standard size),
hard case

Hand-held Slit Lamp Sets	3.5V (NiMH)	3.5V 
BETA R rechargeable handle for mains socket	C-265.20.376	C-265.27.376
same, with HSL 10x Loupe Attachment	C-267.20.376	C-267.27.376
BETA NT rechargeable handle and NT300 charger	C-265.20.420	
same, with HSL 10x Loupe Attachment	C-267.20.420	
 BETA L rechargeable handle and NT 300 charger		C-265.29.420
same, with HSL 10x Loupe Attachment		C-267.29.420

: Rechargeable handle with LED charge status display in the bottom insert.


: To order your kit with the TL "Twist Lock" connection style, simply add the suffix "TL" to the end of the item number. e.g. B-141.20.376 TL for the set with Twist Lock.

110V is the standard mains voltage. If a different voltage is required, i.e. 220V, please specify when ordering.

A company makes history and gives a name to outstanding product quality: HEINE.

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

www.heine.com



■ HEINE SUBSIDIARY COMPANIES
■ HEINE REPRESENTATIVES

NORTH AMERICA

HEINE USA LTD.
HEINE INSTRUMENTS CANADA LTD
10 Innovation Way
Dover, NH 03820
Tel. (603) 742-7103
Fax (603) 742-7217
Toll Free (800) 367-4872
E-Mail: service@heine-na.com

GERMANY

HEINE Optotechnik GmbH & Co. KG
82211 Herrsching
Kientalstr. 7
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
P.O. Box 717 Brookvale 2100
New South Wales
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobeläckerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: heine-schweiz@schaffhausen.ch

Your distributor

11/11. A-000.00.100 e

