


### Straight Laryngeal Mirror




- **A range of accessories make this instrument particularly versatile.**
- **Enhanced XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- **Solid metal construction.** Durable and tough.

Straight Laryngeal Mirror	AV	TL
complete with plane mirror 20mm dia, without handle	B-002.12.101	B-002.12.101 TL
without mirror, without handle	B-002.12.109	B-002.12.109 TL
XHL Xenon Halogen spare bulb	X-002.88.063	X-002.88.063
Plane mirror 20mm dia.		B-000.12.103
Concave mirror 20 mm dia. (magnifying)		B-000.12.104
Metal tongue depressor [01]		B-000.12.105
Transilluminator [02]		B-000.12.106
Protective cap for bulb (for mouth illumination) [03]		B-000.12.107

\*Special Order - 4 to 6 weeks delivery time.

### Curved Laryngeal Mirror



Assists examination with improved view of the larynx.

Curved Laryngeal Mirror	AV
complete with plane mirror 20mm dia. without handle	B-002.12.201
same, without mirror, without handle	B-002.12.208
XHL Xenon Halogen spare bulb	X-002.88.062
Plane mirror 20mm dia.	B-000.12.203
Concave mirror 20 mm dia.	B-000.12.204
Metal tongue depressor [01]	B-000.12.205

\*Special Order - 4 to 6 weeks delivery time.

### HEINE mini 3000® Laryngeal Mirror



Can only be used with the mini 3000 handle system.

- **A range of accessories make this instrument particularly versatile (see above).**
- **Enhanced XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- **Solid metal construction.** Durable and tough.

mini 3000 laryngeal mirror	2.5V
with plane mirror 20mm dia, without handle	D-001.77.101
XHL Xenon Halogen spare bulb	X-001.88.057

\*Special Order - 4 to 6 weeks delivery time.

## Tongue Blade Holder

### For wooden tongue blades




- ⌘ **XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- ⌘ **Light guide for bulb.** Avoids dazzle.

<b>Tongue Blade Holder</b>	<b>AV</b>
for wooden tongue blades without handle	<b>B-002.12.302</b>
XHL Xenon Halogen spare bulb	<b>X-002.88.049</b>
Replacement Light guide for bulb [01]	<b>B-000.12.308</b>

\*Special Order - 4 to 6 weeks delivery time.

## Tongue Blade Holder

### With Blade Ejector for use with HEINE Disposable Tongue Blades



- ⌘ **XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- ⌘ **Light guide for bulb.** Eliminates dazzle.
- ⌘ **Blade ejector.** No need to touch used blades.
- ⌘ **Only for HEINE Disposable Tongue Blades.**
- ⌘ **Disposable plastic tongue blades.** Strong and hygienic.

<b>Tongue Blade Holder</b>	<b>AV</b>
with blade ejector without handle	<b>B-002.12.305</b>
XHL Xenon Halogen spare bulb	<b>X-002.88.049</b>
Pack of 100 disposable plastic tongue blades	<b>B-000.12.304</b>

\*Special Order - 4 to 6 weeks delivery time.

## HEINE mini 3000® Tongue Blade Holder


### With blade ejector for HEINE Disposable Tongue Blades



Can only be used with the mini 3000 handle system.  
Available in black or blue.

- : **Solid metal construction.** Durable and tough.
- : **Finger-tip pressure releases the used blade.** No need to touch the blade after use.
- : **Enhanced XHL Xenon Halogen Technology.** Very bright white light. Concentrated illumination of oral cavity and pharynx.
- : **2 piece handle and head.** Easy to maintain, flexible. Compatible with other mini3000 instrument heads.
- : **High quality handle: Chrome finish upper section/refined plastic.** Shockproof, sturdy, non-slip.
- : **Attachment clip with integrated on/off switch.** Secure. Switches off automatically when placed in the pocket. 20,000 switch cycles guaranteed.
- : **Only for HEINE Disposable Tongue Blades.** Strong and hygienic.
- : **Replaceable batteries.** Size AA, or optional rechargeable battery with mini NT charger\*.

<b>mini 3000 Tongue Blade Holder</b>	<b>2.5V</b>
*Complete with mini 3000 battery handle and 5 disposable blades, with batteries	<b>D-001.74.118</b>
*Tongue blade holder head with bulb, <b>without handle</b>	<b>D-001.74.100</b>



XHL Xenon Halogen spare bulb	<b>X-001.88.037</b>
Pack of 100 disposable plastic tongue blades	<b>B-000.12.304</b>

\*Special Order - 4 to 6 weeks delivery time.



\*To upgrade to a rechargeable instrument, with mini NT charger, please see page 146.



HEINE mini 3000® pocket instruments – a complete range of quality diagnostic instruments for your pocket. Available in black or blue. Please specify the color you prefer when ordering. If not otherwise specified, we will supply black instruments.

## Multi-purpose Illuminator



- ∴ **XHL Xenon Halogen Technology.** Bright, concentrated light.
- ∴ **Wide range of applications.** Used as a diagnostic light or tongue depressor for wooden or disposable blades.

Multi-purpose Illuminator	AV	TL
Multi-purpose Illuminator	<b>B-002.12.320</b>	<b>B-002.12.320TL</b>
XHL Xenon Halogen spare bulb	<b>X-002.88.078</b>	<b>X-002.88.044</b>

## Fiber Optic Nasal Illuminator



Accessory for Multi-purpose Illuminator.

- ∴ **Fiber Optic Nasal Illuminator.** For specula and spreadable nasal speculum.

Fiber Optic Nasal Illuminator	<b>B-000.12.321</b>
-------------------------------	---------------------

## Fiber Optic Nasal Illuminator, complete



[01]

[02]

Fiber Optic Nasal Illuminator	AV	TL
complete with nasal speculum	<b>B-002.12.322</b>	<b>B-002.12.322TL</b>
same, with spreadable nasal speculum	<b>B-002.12.323</b>	<b>B-002.12.323TL</b>
XHL Xenon Halogen spare bulb	<b>X-002.88.078</b>	<b>X-002.88.044</b>
Nasal speculum for above [01]		<b>B-000.11.220</b>
Spreadable nasal speculum for above [02]		<b>B-000.11.231</b>

HEINE ENT Diagnostic Sets

Set complete with: **BETA100 Diagnostic Otoscope**

- 3 reusable specula (B-000.11.222)
- tongue depressor for disposable tongue blades
- straight laryngeal mirror without accessories
- spreadable nose speculum
- hard case



ENT Diagnostic Sets	3.5V (NiMH)	3.5V
BETA R rechargeable handle for mains socket	<b>B-244.20.376</b>	<b>B-244.27.376</b>
BETA NT rechargeable handle and NT 300 charger	<b>B-244.20.420</b>	
BETA L rechargeable handle and NT 300 charger		<b>B-244.29.420</b>

\*Special Order - 4 to 6 weeks delivery time.

Set complete with: **K100 Diagnostic Otoscope**

- 3 reusable specula (B-000.11.222)
- tongue depressor for wooden tongue-blades
- hard case



ENT Diagnostic Sets	3.5V (NiMH)	3.5V
BETA R rechargeable handle for mins socket	<b>B-241.20.376</b>	<b>B-241.27.376</b>
BETA NT rechargeable handle and NT 300 charger	<b>B-241.20.420</b>	
BETA L rechargeable handle and NT 300 charger		<b>B-241.29.420</b>

\*Special Order - 4 to 6 weeks delivery time.


: Rechargeable handle with LED charge status display in the bottom insert.

: To order your kit with the TL "Twist Lock" connection style, simply add the suffix "TL" to the end of the item number. e.g. B-141.20.376 TL for the set with Twist Lock.

**A company makes history and gives a name to outstanding product quality: HEINE.**

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

[www.heine.com](http://www.heine.com)



■ HEINE SUBSIDIARY COMPANIES  
■ HEINE REPRESENTATIVES

**NORTH AMERICA**

HEINE USA LTD.  
HEINE INSTRUMENTS CANADA LTD  
10 Innovation Way  
Dover, NH 03820  
Tel. (603) 742-7103  
Fax (603) 742-7217  
Toll Free (800) 367-4872  
E-Mail: [service@heine-na.com](mailto:service@heine-na.com)

**GERMANY**

HEINE Optotechnik GmbH & Co. KG  
82211 Herrsching  
Kientalstr. 7  
Tel. +49 (0) 81 52-38 0  
Fax +49 (0) 81 52-3 82 02  
E-Mail: [info@heine.com](mailto:info@heine.com)

**AUSTRALIA**

HEINE AUSTRALIA PTY. LTD.  
Unit 9, 98 Old Pittwater Road  
P.O. Box 717 Brookvale 2100  
New South Wales  
Tel. +61 (0) 2-99 38 95 00  
Fax +61 (0) 2-99 39 23 05  
E-Mail: [info@heine.com.au](mailto:info@heine.com.au)

**SWITZERLAND**

HEINE (Schweiz) AG  
Tobeläckerstr. 9  
CH-8212 Neuhausen am Rheinfall  
Tel. +41 (0) 52-6 72 22 66  
Fax +41 (0) 52-6 72 63 77  
E-Mail: [heine-schweiz@schaffhausen.ch](mailto:heine-schweiz@schaffhausen.ch)

Your distributor

11/11. A-000.00.100 e

