

Veterinary ENT diagnostic instruments

RELIABLE, LONG-LASTING INSTRUMENTS

The HEINE veterinary range offers instruments of the highest quality and functionality as well as HEINE XHL® Xenon halogen bulbs and innovative power sources.

HEINE VET instruments are subjected to ongoing quality assurance checks. They stand up to the demands of the busiest veterinary practice and assure optimum levels of diagnosis and treatment.





HEINE Veterinary ENT diagnostic instruments



Slit Illumination Head



BETA200 F.O.
VET Otoscope



BETA 100
VET Otoscope

	Slit Illumination Head	BETA200 F.O. VET Otoscope	BETA 100 VET Otoscope
Housing construction			
Metal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Swivel lens - Magnification			
Glass/Magnification 2x	<input type="checkbox"/>		
Glass/Magnification 2,5x		<input type="checkbox"/>	<input type="checkbox"/>
Illumination			
F.O.	<input type="checkbox"/>	<input type="checkbox"/>	
Direct			<input type="checkbox"/>
Instrument size			
Professional	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Available power sources			
Battery 2.5V	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Rechargeable battery 3.5V	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tips			
Long, closed specula	<input type="checkbox"/>		
Slit metal specula, open	<input type="checkbox"/>		
Short, closed specula	<input type="checkbox"/>		
UniSpec VET disposable tips	<input type="checkbox"/>		
Tip adaptor	<input type="checkbox"/>		
VET reusable specula		<input type="checkbox"/>	<input type="checkbox"/>
AllSpec disposable tips (with adaptor)			<input type="checkbox"/>
Insufflation port		<input type="checkbox"/>	<input type="checkbox"/>
Page	033	035	035

* with Adaptor

Slit Illumination Head and accessories



- ∴ **Solid metal construction.** Will last a lifetime.
- ∴ **Slit instrument head.** Simplifies the use of instruments.
- ∴ **Fiber optics illumination.** Bright, homogeneous light.
- ∴ **Multi-coated, scratch-resistant swivel lens with 1.6x magnification.**
High-resolution image.
- ∴ **Simple open construction.** Easy to clean and disinfect with conventional solutions.
- ∴ **Double bayonet lock for specula.** Specula mounted firmly by double lock.

Slit Illumination Head	2.5V	3.5V
with swivel lens	G-001.21.301	G-002.21.301
XHL Xenon Halogen spare bulb	X-001.88.077	X-002.88.078
Swivel lens	G-000.21.209	

Long, closed specula

for Slit Illumination Head



For use with the Slit Illumination Head.

- ∴ **Metal construction, chrome-plated.** Easy to clean, disinfect and sterilize.
- ∴ **Matt internal surface.** Avoids reflexes.
- ∴ **Double bayonet lock for specula.** Firm mounting of specula to head.

Closed specula for Slit Illumination Head	
57 mm long / 4 mm dia.	G-000.21.330
65 mm long / 6 mm dia.	G-000.21.331
90 mm long / 7 mm dia.	G-000.21.332
150 mm long / 10 mm dia.	G-000.21.350

57 mm 65 mm 90 mm 150 mm
 Ø4 mm Ø6 mm Ø7 mm Ø10 mm

Slit metal specula, open

for Slit Illumination Head



For use with the Slit Illumination Head.

- ∴ **High-quality chrome-plated metal construction.** Easy to clean, disinfect and sterilize.
- ∴ **Matt internal surfaces.** Reduced reflexes.
- ∴ **60° opening with rounded edges.** Ideal for instrumentation, atraumatic insertion.
- ∴ **Double bayonet fitting.** Positive fit.

Slit metal specula for Slit Illumination Head	
65 mm long / 6 mm dia.	G-000.21.341
90 mm long / 7 mm dia.	G-000.21.342

65 mm long 90 mm long
 Ø 6 mm Ø 7 mm

Short, closed specula

for Slit Illumination Head and BETA 100 VET Otoscope



Ø5,5 mm Ø 4,5mm Ø3,5mm



Ø2,8 mm Ø 2,2mm



- ⌘ **Made of high-density SANALON S with a smooth external surface.** Dust and germs are easily removed. Disinfection and autoclaving at up to 134°C is possible.
- ⌘ **Ergonomic shape.** Painless atraumatic insertion.
- ⌘ **Integral metal ring with 2 bayonet locks.** Firm mounting of specula to head.
- ⌘ **40mm long.**

Short, closed SANALON S specula

for Slit Illumination Head and BETA 100 VET Otoscope

1 set = 5 veterinary specula in SANALON S	G-000.21.316
40mm long / 2.2mm dia.	G-000.21.310
40mm long / 2.8mm dia.	G-000.21.311
40mm long / 3.5mm dia.	G-000.21.312
40mm long / 4.5mm dia.	G-000.21.313
40mm long / 5.5mm dia.	G-000.21.314

HEINE UniSpec® disposable ear specula



Ø2,5mm

Ø4 mm



[01]

The adaptor illustrated [01] is necessary when using UniSpec specula with the Slit Illumination Head.

- ⌘ **No cleaning necessary.** Saves time.
- ⌘ **Hygienic and safe.** Eliminates cross-infection.
- ⌘ **Packed in units of 1000.** Assures a competitive price.
- ⌘ **Original HEINE design and manufacture.** Guaranteed good fit to the instrument head.

UniSpec disposable ear specula (grey)

Pack of 1000 UniSpec disposable specula 2.5 mm dia.	B-000.11.242
Pack of 1000 UniSpec disposable specula 4 mm dia.	B-000.11.241

Adaptor for use with Slit Illumination Head [01]	G-000.21.302
--	---------------------

⊗: Single use only.

Spreadable speculum for Slit Illumination Head



This versatile instrument is widely used as a universal speculum for small animals.

- ⌘ **For use with the Slit Illumination Head.**
- ⌘ **Distal jaw opens up to 28 mm with 4 click-stops.** Useful for examination and therapy.
- ⌘ **High-quality, chrome-plated metal construction.** Easy to clean, disinfect and sterilize.
- ⌘ **Matt internal surfaces.** Avoids reflexes.
- ⌘ **Oval branches with rounded edges.** Ideal for instrumentation, atraumatic.
- ⌘ **2 bayonet locks.** Firm mounting of spreadable Speculum to head.

Spreadable speculum for Slit Illumination Head (95mm long)	G-000.21.360
--	---------------------

HEINE BETA200® Fiber Optic (F.O.) VET Otoscope

Veterinary Otoscope with Fiber Optic illumination



Solid, all-metal construction. A quality instrument for the Vet's office with Fiber Optic (F.O.) illumination. Insufflation port. The swivel window features an airtight seal.

- ∴ **Distal Fiber Optic illumination.** Perfectly even illumination of auditory canal and tympanum.
- ∴ **Large, scratch-resistant glass viewing window with 2.5x magnification.** Unobstructed view, distortion-free image.
- ∴ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.
- ∴ **Insufflation port.** Reliable test of tympanic mobility.
- ∴ **Durable metal body, matt chrome-plated.** Unbreakable, long-life.
- ∴ **XHL Xenon Halogen Technology.** Extremely bright, white light.
- ∴ **For use with reusable SANALONS VET Specula (G-000.21.213).**
- ∴ **Instrument head matt-black inside.** Eliminates reflexes.



Fiber Optic (F.O.) illumination

BETA 200 VET Otoscope	2,5V	3,5V
without specula, without handle	G-001.21.250	G-002.21.250
XHL Xenon Halogen spare bulb	X-001.88.077	X-002.88.078

HEINE BETA 100® VET Otoscope

Veterinary Otoscope with direkt illumination



Otoscope with all-metal construction and direct Xenon Halogen illumination. Suitable for use with short, closed specula. With the use of the VET adaptor, it is possible to fit reusable VET specula for ear examinations.

- ∴ **All-metal housing.** Exceptionally-tough and long-lasting.
- ∴ **XHL Xenon Halogen Technology.** Extremely bright, white light.
- ∴ **Large, scratch-resistant glass viewing window with 2.5x magnification.** Unobstructed view, distortion-free image.
- ∴ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.
- ∴ **Insufflation port.** Reliable test of tympanic mobility.
- ∴ **Can be used with short, closed specula G-000.21.316 (no adaptor required).** Versatile, useful for very small animals.
- ∴ **Instrument head matt-black inside.** Eliminates reflexes.



Direct Xenon Halogen illumination

BETA 100 VET Otoscope	2,5V	3,5V
without specula, without handle	G-001.21.260	G-002.21.260
XHL Xenon Halogen spare bulb	X-001.88.037	X-002.88.049
Adaptor for reusable SANALON S VET specula [01]		G-000.21.214
Adaptor for AllSpec tips [02]		B-000.11.306

[01]



[02]

Reusable VET specula



- ⚡ **Made of high-density SANALON S with a smooth external surface.** Dust and germs are easily removed. Disinfection and autoclaving at up to 134° C is possible.
- ⚡ **Ergonomic shape.** Atraumatic.
- ⚡ **Convenient fixation system.** Positive fit to otoscope head.

Reusable SANALON S specula for BETA200 VET Otoscope

1 set = 3 reusable SANALON S specula	G-000.21.213
61 mm long / 4 mm dia.	G-000.21.210
61 mm long / 6 mm dia.	G-000.21.211
61 mm long / 9 mm dia.	G-000.21.212

HEINE AllSpec® Disposable tips



- ⚡ **For single use.** Hygienic, no cross-infection.
- ⚡ **Original HEINE design.** Firm fit to the otoscope head.
- ⚡ **Assured quality.** No sharp edges.

AllSpec Disposable tips (grey) for BETA 100VET Otoscope

Pack of 1000 pcs. 2.5 mm dia. (infants) [03]	B-000.11.128
Pack of 1000 pcs. 4 mm dia. (adults) [03]	B-000.11.127
Case of 10 dispenser packs of 250 tips, 2.5 mm dia (infants) [01]	B-000.11.151
Case of 10 dispenser packs of 250 tips, 4 mm dia (adults) [01]	B-000.11.150
AllSpec Tip-dispenser [02]	B-000.11.148

AllSpec Tip-dispenser [02]

The right size, visible and instantly available. A clean new tip for each and every patient. Contents: 50 pcs. Ø 2.5 mm and 50 pcs. Ø 4 mm.

ⓧ: Single use only.

Insufflation Bulb



For BETA200VET and BETA100VET Otoscope.

- ⚡ **For pneumatic testing.**

Insufflation bulb	B-000.11.240
-------------------	---------------------

HEINE Veterinary Diagnostic Sets

Set complete with: **Slit illumination head**
 1 set long, closed specula (57 x 4 mm, 65 x 6 mm, 90 x 7 mm)
 1 spare bulb
 hard case



Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V
BETA battery handle	G-112.10.118		
BETA R rechargeable handle for mains socket		G-112.20.376	G-112.27.376
BETA TR rechargeable handle with plug-in transformer		G-112.20.384	G-112.27.384
BETA NT rechargeable handle and NT 300 charger		G-112.20.420	
BETA L rechargeable handle and NT 300 charger			G-112.29.420

Set complete with: **Slit illumination head, BETA 200 Ophthalmoscope**
 1 set long, closed specula (57 x 4 mm, 65 x 6 mm, 90 x 7 mm)
 1 spare bulb for each instrument
 hard case





Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V
BETA battery handle	G-148.10.118		
BETA R rechargeable handle for mains socket		G-148.20.376	G-148.27.376
BETA TR rechargeable handle with plug-in transformer		G-148.20.384	G-148.27.384
BETA NT rechargeable handle and NT 300 charger		G-148.20.420	
BETA L rechargeable handle and NT 300 charger			G-148.29.420

: Rechargeable handle with charge status display in the bottom insert.

HEINE Veterinary Diagnostic Sets





Set complete with: **Slit illumination head, BETA 200 Ophthalmoscope**
 5 short, closed SANALON S specula
 10 UniSpec disposable ear specula with adaptor
 1 set = 4 long, closed specula (57 x 4 mm, 65 x 6 mm, 90 x 7 mm, 150 x 10 mm)
 2 slit specula
 spreadable speculum
 1 spare bulb for each instrument
 hard case

Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V 
BETA battery handle	G-150.10.118		
BETA R rechargeable handle for mains socket		G-150.20.376	G-150.27.376
BETA TR rechargeable handle with plug-in transformer		G-150.20.384	G-150.27.384
BETA NT rechargeable handle and NT 300 charger		G-150.20.420	
 BETA L rechargeable handle and NT 300 charger			G-150.29.420

Set complete with: **BETA 200 Ophthalmoscope, BETA 200 VET Otoscope**
 3 reusable SANALONS VET specula
 1 spare bulb for each instrument
 zipper case




Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V 
BETA battery handle	G-151.10.118		
BETA R rechargeable handle for mains socket		G-151.20.376	G-151.27.376
BETA TR rechargeable handle with plug-in transformer		G-151.20.384	G-151.27.384
BETA NT rechargeable handle and NT 300 charger		G-151.20.420	
 BETA L rechargeable handle and NT 300 charger			G-151.29.420

: Rechargeable handle with charge status display in the bottom insert.

A company makes history and gives a name to outstanding product quality: HEINE.

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

www.heine.com



■ HEINE SUBSIDIARY COMPANIES
■ HEINE REPRESENTATIVES

GERMANY

HEINE Optotechnik GmbH & Co. KG
82211 Herrsching
Kientalstr. 7
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com

NORTH AMERICA

HEINE USA LTD.
10 Innovation Way
Dover, NH 03820
Tel. (603) 7 42-71 03
Fax (603) 7 42-72 17
Toll Free (800) 367-4872
E-Mail: service@heine-na.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
P.O. Box 717 Brookvale 2100
New South Wales
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobeläckerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: heine-schweiz@schaffhausen.ch

Your distributor

10/11. A-000.00.100 e

