

HEINE HSL 150® Hand-held Slit Lamp

For the examination of the anterior segment



∴ Continuous slit adjustment

∴ Blue filter

∴ Ideal for domiciliary visits

Compact and lightweight, the portable HSL 150 is ideal for situations where a full-size slit lamp is not available or convenient. For the examination of the anterior segment in Human- and Veterinary medicine.

- ∴ **Slit size adjustable** from 0.2mm x 10mm to 4 mm x 14 mm.
- ∴ **6x magnification.**
- ∴ **Blue interference filter (FITC)**, can be selected for corneal examination.
- ∴ **Multi-coated optics for maximum light transmission.**
- ∴ **2.5V or 3.5V XHL Xenon Halogen Technology** for bright, white light comparable with the brightness of a classic slit lamp.
- ∴ **The HSL 150 instrument head weighs just 70g**, the lightest slit lamp of its kind.

HSL 150 for BETA or BETA SLIM handles	2.5V	3.5V
HSL 150, head only	C-001.14.602	C-002.14.602
XHL Xenon Halogen spare bulb	X-001.88.098	X-002.88.099
Replacement filter attachment for HSL 150	C-000.14.605	

HEINE HSL 10x® Loupe Attachment

Loupe attachment with 10x magnification



For HSL 150. With precision magnification optics for crisp, clear images.

- ∴ **10x magnification** when combined with HSL 150 Loupe.
- ∴ **Increased working distance** from examiner to patient.
- ∴ **Large field of view.**
- ∴ **With an additional detachable eyepiece** to steady the instrument on the brow and avoid stray light. Can be replaced with a soft cap, suitable for spectacle-wearers.

HSL 10x Loupe Attachment	C-000.14.606
--------------------------	--------------



HEINE HSL 150® Hand-held Slit Lamp Sets



[01]





[02]

Set complete with: **HSL 150 Hand-held Slit Lamp**
BETA SLIM handle (compact)

Hand-held Slit Lamp Sets	2.5V	3.5V
BETA SLIM battery handle, in soft pouch, without spare bulb [01]	C-252.10.105	
same, with HSL 10x Loupe Attachment, spare bulb, in hard case [02]	C-268.10.105	
BETA SLIM rechargeable handle and NT 300 charger, spare bulb, in hard case		C-266.20.470
same, with HSL 10x Loupe Attachment, in hard case		C-268.20.470

Set complete with: **HSL 150 Hand-held Slit Lamp**
spare bulb
BETA handle (standard size)
hard case




Hand-held Slit Lamp Sets	2.5V	3.5V (NiMH)	3.5V 
BETA battery handle	C-265.10.118		
same, with HSL 10x Loupe Attachment	C-267.10.118		
BETA R rechargeable handle for mains socket		C-265.20.376	C-265.27.376
same, with HSL 10x Loupe Attachment		C-267.20.376	C-267.27.376
BETA TR rechargeable handle with plug-in transformer		C-265.20.384	C-265.27.384
same, with HSL 10x Loupe Attachment		C-267.20.384	C-267.27.384
BETA NT rechargeable handle and NT 300 charger		C-265.20.420	
same, with HSL 10x Loupe Attachment		C-267.20.420	
 BETA L rechargeable handle and NT 300 charger			C-265.29.420
same, with HSL 10x Loupe Attachment			C-267.29.420

: Rechargeable handle with charge status display in the bottom insert.

A company makes history and gives a name to outstanding product quality: HEINE.

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

www.heine.com



■ HEINE SUBSIDIARY COMPANIES
■ HEINE REPRESENTATIVES

GERMANY

HEINE Optotechnik GmbH & Co. KG
82211 Herrsching
Kientalstr. 7
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com

NORTH AMERICA

HEINE USA LTD.
10 Innovation Way
Dover, NH 03820
Tel. (603) 7 42-71 03
Fax (603) 7 42-72 17
Toll Free (800) 367-4872
E-Mail: service@heine-na.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
P.O. Box 717 Brookvale 2100
New South Wales
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobelackerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: heine-schweiz@schaffhausen.ch

Your distributor

10/11. A-000.00.100 e

