

HEINE DELTA 30 dermatoscope

THE MOST IMPORTANT POINTS AT A GLANCE



HEINE DELTA 30 dermatoscope



charge state
66 - 100 %



charge state
33 - 66 %



charge state
10 - 33 %



flashes
charge state
< 10 %



For a secure fit, the contact plate snaps into place with a clearly audible click.



Diode lights up
blue = polarised
illumination

Diode does not
light up = non-polarised
illumination

A FILM SAYS MORE THAN A THOUSAND WORDS



The "Application" video shows the DELTA 30 in action.



Remember each and every mole.

The HEINE DERManager* is the ideal complement to your DELTA 30. This software solution makes digital management of dermatoscopic image data simple.

*Check on our website in which countries the DERManager is currently available.

