

**Mobility and quality have now
become *one*. The new
HEINE DELTA*one* Dermatoscope.**



30
YEARS
1st DERMATOSCOPE



Small on the outside. Big on the inside.

The new HEINE DELTAone Dermatoscope.



Extremely sharp image with the new achromatic HEINE optical system

Polarisation eliminates reflections and glare

Fast change from contact to non-contact examination

LEDHQ
LED NOW IN HEINE QUALITY.
Excellent colour rendering due to LEDHQ

Fits into any pocket

Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System*

- ∴ **Crystal clear image**, no matter if you prefer a closer or greater working distance
- ∴ **True field of view of 21 mm (lens 22 mm) Ø**
- ∴ **Brightness control**: 3 levels of light intensity
- ∴ **Easily cleaned with spray disinfectant**
- ∴ **10x magnification**
- ∴ **No cross-contamination**: autoclavable contact plate with scale
- ∴ **Dioptre compensation**
- ∴ **Switch between polarised and non-polarised illumination mode at the touch of a button** = toggle function. For more information, please visit heine.com/en/Toggle

- ∴ **Rechargeable battery with Li-ion technology**: no memory effect; safely recharge the battery regardless of its state of charge
- ∴ **Small contact plate** for difficult-to-access pigmented lesions
- ∴ **3 possibilities of digital documentation**:
 - With the exclusive HEINE DERM App: Bodymap, patient folder, e-mail function and much more
 - Wireless with the HEINE Cube System for up to 10 doctors and HEINE dermatoscopes
 - Conventional digital images with many popular smartphones

DELTAone Standard Edition (black-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-210.28.305
DELTAone White Edition (white-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-211.28.305

∴ **Accessories and further articles under heine.com/en/DELTAone**

*Not available in all countries.