

HEINE mini-c hand-held examination light

Pocket size. Rocket performance.
And an eye on safety.



Robust, durable examination light in a compact pocket size.
Bright, concentrated light – now with unique LED^{HQ}.

- ⚡ Safe examination of the eyes due to tested photobiological safety.
- ⚡ Comfortable and ergonomic.
- ⚡ Natural colours and virtually unlimited operating life of the LED^{HQ}.
- ⚡ Switches off automatically when placed in the pocket.
- ⚡ 20,000 on/off cycles guaranteed.
- ⚡ Replaceable batteries (size: AAA).

Ordering information

HEINE mini-c hand-held examination light with batteries

Art.-No.

D-008.73.109

