

HEINE QUALITY
MADE IN GERMANY

∴ **HEINE® EN 200-1 Wall Transformer
with DELTA® 20 T Dermatoscope**

Turns itself on and off. Always ready, never misplaced.



HEINE® EN 200-1

A space-saving and always ready-to-use solution at your fingertips.

EN 200-1 WALL TRANSFORMER

- ✦ **Built-in theft protection** due to the integrated locking screw. So the instrument is always in place.
- ✦ **Excellent energy efficiency:** Instruments turn on and off automatically. When they are removed and replaced. And after 15 minutes of non-usage when taken off, the instruments are also switched off. When the instrument is taken off, the most recently selected brightness level is enabled. This saves time and power.
- ✦ **Large operating radius** thanks to spiral cord with a length of up to 3.50 m.
- ✦ **Battery-free power supply.**



DELTA 20 T DERMATOSCOPE

- ✦ **Toggle function** for switching between polarised and non-polarised illumination mode.
- ✦ **High-resolution images** thanks to the high-quality achromatic optical system.
- ✦ **Detailed view of even the smallest structures;** 10 to 16x magnification is achieved with the focusing system.
- ✦ **Glare-free and reflection-free view thanks to polarisation,** even without immersion fluid.
- ✦ **With just one hand:** individual focusing, range of correction approximately -6 to +6 dpt.
- ✦ **Digital documentation:** using the HEINE SLR photo adaptor, the DELTA 20 T can be connected to the SLR cameras of major brands without any issues.

As a matter of course, a number of our instrument heads can be combined with the EN200-1's handle.

The EN 200-1 wall transformer meets the current standard IEC 60601-1:

medical electronic equipment – part 1: General safety requirements (edition 3.1), which have been intensified in some respects.



LED NOW IN HEINE QUALITY.

There's LED. And there's LED^{HQ}. HEINE has set a new standard which maintains that only the best is good enough, from the selection of materials to processing, from light intensity to dimmability, and from thermal management to a colour rendering index as high as possible. **That's what we call LED in HEINE Quality – or LED^{HQ}.**

EN200-1 Wall Transformer (Single Unit with one handle) as single module with long USB cord and E4-USB plug-in power supply **incl. DELTA 20 T Dermatoscope Head** with contact plate Ø 23 mm **with scale, dermatoscopy compendium**

K-010.28.252

For additional product configurations and further HEINE Quality Products please visit heine.com or give us a call!



HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7, 82211 Herrsching, Tel. +49 (0) 8152-38 0
Fax +49 (0) 8152-382 02, E-Mail: info@heine.com
www.heine.com

