## HEINE OMEGA 600 indirect ophthalmoscope

## A QUICK LOOK AT THE MOST IMPORTANT FEATURES.

We thank you for choosing HEINE OMEGA 600 and wish you a lot of success and fun working with your new instrument.





The video "How to set up your OMEGA 600" explains how to adjust the instrument to your personal needs.

## **HEINE OMEGA 600** indirect ophthalmoscope



## Stereoscopic adjustment lever

Best possible 3D imaging for any pupil size and peripheral views



Flip-up function of the optics unit For more convenient patient interactions



Eyepieces: To best accommodate your view, exchange the initially mounted eyepieces +2D with the neutral lenses 0D (supplied with the instrument)



Width adjustment

If you want to share your experiences or discoveries, we would be highly pleased to see it: #heine #iuseheine #OMEGA600











Teaching mirror TM2



Dust cover for front glass

