



Here comes the sun. Upgrade from XHL to the LED^{HQ} module.

CHANGING BULBS IS NOW A THING OF THE PAST!

With this new option, the sun is rising for many owners of HEINE instruments with XHL illumination.

- ▶ Virtually unlimited durability with LED^{HQ} in HEINE quality
- ▶ Maintenance-free and reliably ready for use
- ▶ Longer running times for both single-use batteries and rechargeable batteries thanks to lower power consumption
- ▶ Sustainable, as this is not a single-use product and there is less battery waste
- ▶ Tested spare part guarantees patient safety*

Instrument	Type	Power source	Previous XHL bulb	New LED ^{HQ} module
Otoscope	BETA 400 F.O. BETA 200 F.O. K180	Battery handle 2,5V	X-001.88.077	X-028.88.132
		Rechargeable handle 3,5V	X-002.88.078	X-008.88.133
Ophthalmoscope	BETA 200S BETA 200	Battery handle 2,5V	X-001.88.069	X-028.88.134
		Rechargeable handle 3,5V	X-002.88.070	X-008.88.135
Retinoscope	BETA 200	Battery handle 2,5V	X-001.88.087	X-028.88.130
		Rechargeable handle 3,5V	X-002.88.089	X-008.88.131
Veterinary	G100 slit illumination head BETA 200 VET otoscope	Battery handle 2,5V	X-001.88.077	X-028.88.132
		Rechargeable handle 3,5V	X-002.88.078	X-008.88.133



- ▶ In this short video we show you how the upgrade works and how to insert the new HEINE LED^{HQ} module instead of an XHL bulb.

*Only the original is as good as the original. Use HEINE ORIGINAL PARTS and not third-party products from other manufacturers. Only then will our guarantee apply and performance features such as colour rendering, function and safety of our examination instruments are always guaranteed.

