



**| For every intubation.
| Even the strange ones.**



The HEINE visionPRO video laryngoscope.

WITH VERY, VERY, VERY GOOD IMAGE QUALITY.

For fast and reliable intubation, nothing is more important than the best view of the glottis – with a bright and sharp image. This is true for every intubation – and even more so for difficult, unforeseen situations. That's why we've combined the best components in the HEINE visionPRO video laryngoscope, from the camera to the image processing to our HEINE allBRIGHT display.

EXTREMELY DURABLE. AFTER ALL, IT IS MADE IN GERMANY.

The HEINE visionPRO is extremely robust and reliable. The grip head and display frame are made of anodized aluminum, while the handle and camera arm are made of stainless steel. The cover glass of the allBRIGHT display is highly protective. After all, the high reliability of HEINE instruments is legendary. The same applies to our video laryngoscope: we make Made in Germany.

Exclusive:

The HEINE allBRIGHT display.

Our allBRIGHT display provides true color and lifelike imaging even in bright ambient light. No wonder, because HEINE is the first and only supplier of video laryngoscopes worldwide to offer a transflective display.*

STRAIGHT TO THE TARGET.

With its extra-large field of view in portrait format, the HEINE visionPRO shows a significantly larger section of the glottis and the relevant surrounding anatomy. This means that quick and easy orientation – even in time-critical situations – is always well supported. Intuitive to handle, HEINE visionPRO allows the user to focus fully on the intubation.


SAFE, FAST AND EASY TO CLEAN.

To ensure maximum patient safety, the display and handle can be easily separated for numerous validated reprocessing procedures and are Joint Commission compliant. The handle can be wipe disinfected or even sterilized by low-temperature gas plasma. The assembled instrument can be hygienically reprocessed with low-level disinfection. This is thanks to our smooth and easy-to-clean surfaces called HEINE smoothSURFACE.

* Data on file

allBRIGHT





As eco-friendly as possible.

SUSTAINABILITY BUILT IN.

The HEINE visionPRO features a rechargeable lithium battery, which we chose for its near infinite life. This reduces costs and, above all, the impact on the environment. Not only are single-use batteries expensive but are also highly problematic to manufacture and dispose of – and so not sustainable.

FIRST BLADE MADE FROM UPCYCLED MATERIAL.

And one more real benefit at the end: our visionPRO blades are made of upcycled material from discarded refrigerators.

By recycling these refrigerators in a controlled manner, we spare the environment around 200 metric tons of CO₂ annually already in the very first year. This is equivalent to 3.8 million miles of rail travel per person, or about 140 times around the world. What's more, we manufacture the blades for the visionPRO video laryngoscope on a machine that runs on greensourced electricity. Learn more about sustainability at www.heine.com/en/CSR.

Various national and international industrial property rights have been registered and granted for the smart design of our blades.

Quality made in Germany.

We live and breathe quality. In everything we do. Every product we make. Every lens we polish. Every screw we tighten. Every optic fiber we bundle. Always. Every day.

The durability of our primary diagnostic instruments is legendary. Every single instrument is the result of decades of product development, continuous improvements, production and a 100 % quality control in our own facilities. We make certain that our high expectations of quality are maintained.

It is how we ensure that our high-end products meet our claim of providing the most accurate and earliest possible diagnosis. No more, no less.



[01]

visionPRO set Mac [01]

visionPRO video laryngoscope consisting of visionPRO HR1 Mac handle and visionPRO DS1 display unit, visionPRO CS1 charging station with plug-in power supply

F-270.95.862

You can find ordering information and other products at
www.heine.com/laryngoscopes



Mac 2



Mac 3



Mac 4



HyMac 3

HEINE Optotechnik GmbH & Co. KG

Dornierstr. 6, 82205 Gilching, Germany
www.heine.com



We reserve the right to change
specification without notice.

02|25. A-000.00.328 EN-US

