

HEINE QUALITY
MADE IN GERMANY

∴ HEINE BETA 400 F.O. Otoscope

The highest magnification for the best diagnosis.



THE HEINE BETA 400 FIBER OPTIC OTOSCOPE. 4.2x MAGNIFICATION. EXCEPTIONAL DEPTH OF FIELD. ONE HAND OPERATION.

The HEINE BETA 400 F.O. Otoscope is the most advanced high quality otoscope making more details visible, while maintaining familiar handling and functionality. With a magnification factor of 4.2x and multi-coated precision optics for sharp spatial images, it truly provides for the best diagnosis. Combined with a large field of view the ear canal and tympanum are clearly visible at a glance. The unique opening mechanism (patent pending) can be operated with one hand and flips up the optical system. Application of instruments through the device becomes easily possible.



4.2x

4.2x MAGNIFICATION

More detail of anatomical structures and easy recognition of foreign bodies.



ONE HAND OPERATION

Easy thumb control: Optical system to be flipped open and out of the way, enabling instrumentation and ease of cleaning.



EXCEPTIONAL DEPTH OF FIELD

Precise image: Provides sharp spatial images in a convenient working range.

BETA 400 F.O. Otoscope with 1 set (4 pcs.) of reusable tips and 10 AllSpec disposable tips 4 mm Ø, spare bulb, **BETA battery handle 2.5V**, hard case [01]

2.5V: B-143.10.118

BETA 400 F.O. Otoscope with 1 set (4 pcs.) of reusable tips and 10 AllSpec disposable tips 4 mm Ø, spare bulb, **BETA4 USB rechargeable handle with USB cord and plug-in power supply**, hard case

3.5V: B-143.27.388

LED: B-143.28.388

BETA 400 F.O. Otoscope with 1 set (4 pcs.) of reusable tips and 10 AllSpec disposable tips 4 mm Ø, spare bulb, **BETA4 NT rechargeable handle with NT4 table charger**, hard case [02]

3.5V: B-143.23.420

LED: B-143.24.420

Further sets on request.

