

# : HEINE SIGMA® 250 Binocular Indirect Ophthalmoscope

Accurate diagnosis made easy.



# HEINE SIGMA<sup>®</sup> 250 BINOCULAR INDIRECT OPHTHALMOSCOPE Spectacle-mounted lightweight with a multitude of benefits.



High tensile carbon fiber construction is lightweight and has excellent stability.



# HEINE mPack mini – The compact and light power source for SIGMA 250.

Precision HEINE Optics.

The new, high performance power source HEINE mPack mini offers freedom from cables that tie the users to desk, wall or static power sources. With a weight of just 95 g and dimensions of 44 x 23 x 101 mm, this mobile power source is now approximately approx. 70% smaller and lighter than its predecessor.

NEW: Optimised mPack mini belt clip for simple attachment of the power source to the belt or waistband. The mPack mini can also be charged in the belt clip. Thanks to the Li-ion battery technology, the power source has an operating time of typ. 11 hours in continuous operation at maximum power.

#### SIGMA 250 Kit with S-FRAME

SIGMA250, S-FRAME, mPack mini with E4-USB plug-in transformer, Retaining cord, Cleaning fluid for S-FRAME, mPack mini Belt Clip

#### C-008.33.341

## SIGMA250 M2\* Kit with S-FRAME

SIGMA250 M2, S-FRAME, mPack mini with E4-USB plug-in transformer, Retaining cord, Cleaning fluid for S-FRAME, mPack mini Belt Clip

## C-008.33.346

#### SIGMA 250 Set with S-FRAME

SIGMA250, S-FRAME, mPack mini with E4-USB plug-in transformer, Scleral depressor small, Scleral depressor large, A.R. 20D Ophthalmology Lens, Teaching Mirror, Cobalt Blue and Yellow Filters, Diffusor attachment, Retaining cord, Cleaning fluid for S-FRAME, mPack mini Belt Clip, Carrying case

#### C-281.41.671

#### SIGMA 250 M2\* Set with S-FRAME

SIGMA250 M2, S-FRAME, mPack mini with E4-USB plug-in transformer, Scleral depressor small, Scleral depressor large, A.R. 20D Ophthalmology Lens, Teaching Mirror, Cobalt Blue and Yellow Filters, Diffusor attachment, Retaining cord, Cleaning fluid for S-FRAME, mPack mini Belt Clip, Carrying case

# C-279.41.671

\* The SIGMA 250 M2 has a micro spot aperture instead of small aperture Aperture SIGMA250: Aperture SIGMA 250 M2: 

Further information on HEINE products can be found at www.heine.com



 $\geq$  90, good index for red colors.

HEINE Optotechnik GmbH & Co. KG, Kientalstr. 7, 82211 Herrsching, Germany Tel. +49 (0) 8152/38-0, Fax +49 (0) 8152/38-202, E-Mail: info@heine.com, www.heine.com