



Reprocessing Agents for HEINE visionPRO

HEINE Optotechnik GmbH & Co. KG has tested and approved the following agents for material compatibility for disinfection of the HEINE visionPRO video laryngoscope.

Manufacturer	Product name
B.Braun	Meliseptol Wipes Sensitive
	Meliseptol Wipes Ultra
Bode Chemie / Hartmann	Bacillol Tissues
	Mikrobac Tissues
	Mikrobac viruzidal
Dr. Schuhmacher	Cleanisept Wipes
	Cleanisept Wipes forte
Ecolab Healthcare	Incides N
	Sani-Cloth Active
PDI	Sani-cloth AF3
Schülke	Mikrozid AF wipes

The validated reprocessing procedure is described in the reprocessing instruction MED235566.