

Page	Index	Page	Index
	A		H
086/142	Adaptors for NT4 Table Charger	055	Hand-held Slit Lamps
	B	055	- HSL 150
120–127	Binocular Loupes	060	Hand-held Slit Lamp Accessories
126	- C2.3/C2.3K	060	- HSL 10x Loupe Attachment
128/123	- HR		Handles
124/125	- HRP	140	- BETA Handle Module
128	Binocular Loupe Accessories	140	- BETA Battery Handle
128	- Spectacle frame S-FRAME	141	- BETA4 NT Rechargeable Handle
128	- Headband	141	- BETA4 USB Rechargeable Handle
129/130	- LoupeLight2	142	- BETA4 SLIM NT Rechargeable Handle
154–158	Bulbs, XHL Xenon Halogen Technology	140	- BETA SLIM Battery Handle
	C	084	- EasyClean LED Laryngoscope Handle SHORT
162/163	Cases	084	- EasyClean LED Laryngoscope Handle SLIM
	Chargers	084	- EasyClean LED Laryngoscope Handle Standard
148	- mini NT	085	- F.O.4 SHORT XHL/LED NT Laryngoscope Handle
142	- NT4 and Adaptor	085	- F.O.4 SLIM XHL/LED NT Laryngoscope Handle
028–029	Clip Lamps	086	- F.O. Angled Laryngoscope Handle
028	- mini-c	147	- Large Battery Handle
028	- mini 3000 Clip Lamp	147	- mini 3000 Battery Handle
029	- mini 3000 Combi Lamp	148	- mini 3000 Rechargeable Handle
030	Clip Lamps Accessories	085	- Standard F.O.4 XHL/LED NT Laryngoscope Handle
	D	083	- XP Disposable Laryngoscope Handle Shell
090–096	Dermatoscopes		Handles – Upgradekits
▶ 091	- DELTA 30	143	- BETA4
▶ 092	- DELTAone	086	- F.O.4
094	- DELTA 20T	131–138	Headlights
093	- iC1	132–133	- MicroLight2
096	- mini 3000	134–137	- ML4 LED HeadLight
096	- mini 3000 LED	134–136	- ML4 LED HeadLight UNPLUGGED
095	- NC2	138	- MD 1000 F.O.
097–099	Dermatoscope Accessories		L
097	- DELTA 30 Accessories		Lamps
097	- DELTAone Accessories	028	- mini-c Clip Lamp
098	- DELTA 20T Accessories	029	- mini 3000 Combi Lamp
097	- iC1 Accessories	078–088	Laryngoscopes
099	- mini 3000/mini 3000 LED Accessories		Laryngoscopes Classic+ F.O.
098	- NC2 Accessories	079	- Macintosh F.O.
144–146	Diagnostic Centre EN 200	080	- Miller F.O.
	E	080	- Paed F.O.
084	EasyClean LED Laryngoscope Handles	087	- Sets
116–117	Examination Lights	080	- WIS F.O.
116	- EL 10 LED		Laryngoscopes FlexTip+ F.O.
117	- EL3 LED	081	- FlexTip+ F.O.
	F	088	- Sets
160	Fiber Optic Cable for HEINE Projectors		Laryngoscopes XP, Disposable
058/061	Finoff Transilluminator	088	- Emergency Set
073	Fundus Charts	082	- Macintosh
059/061	Funduscope iC2	082	- Miller
			M
		148	mini NT table charger
		151	mPack
		153	mPack mini
		150	mPack UNPLUGGED

Page	Index
	O
057/061	Ophthalmic Examination Lamps
042–051	Ophthalmoscopes, direct
045	- BETA 200
045	- BETA 200 LED
043	- BETA 200S
043	- BETA 200S LED
048	- K180
050	- mini 3000
050	- mini 3000 LED
052	Ophthalmoscopes, direct; Accessories
052	- BETA 200 M2
062–071	Ophthalmoscopes, indirect
071	- BINOCULAR, hand-held
070	- MONOCULAR, hand-held
063/064	- OMEGA 500, binocular, head-worn
066	- OMEGA 500 UNPLUGGED
067	- OMEGA 500 with DV1
068	- SIGMA 250 binocular, spectacle-worn
072–075	Ophthalmoscopes, indirect; Accessories
072	- Trainer, Model Eye
073	- Ophthalmoscopy Lens A.R.
073	- Combi-Case / Hard case
	Otosopes
020	- BETA 100
013	- BETA 200 F.O.
013	- BETA 200 LED F.O.
011	- BETA 400 F.O.
011	- BETA 400 LED F.O.
036	- BETA 100 VET (Veterinary)
035	- BETA 200 LED VET F.O. (Veterinary)
035	- BETA 200 VET F.O. (Veterinary)
022	- K100
015	- K180 F.O.
016	- mini 3000 F.O.
016	- mini 3000 LED F.O.
023	- mini 3000
034	Otoscope, Operating
018/019	Otosopes with F.O. illumination Accessories
019	- Insufflation Bulb for BETA 200, BETA 400, K180
019	- Insufflation Bulb for mini 3000 F.O.
018	- Tips/Specula
024/025	Otosopes with direct illumination, Accessories
025	- Insufflation Bulb for BETA 100, K100
024	- Tips/Specula
038–040	Otosopes VET, Accessories
040	- Insufflation Bulb for BETA 100 VET, 200 VET
038/039	- Tips/Specula
	P
140–153	Power Sources
140–143	- BETA Handles
144–146	- EN 200 Wall Transformer
152	- EN 50 Wall/Table Transformer
150	- EN 50 UNPLUGGED Wall Transformer
149	- HC 50L Headband Rheostat
151	- mPack
153	- mPack mini
150	- mPack UNPLUGGED
147	- mini 3000 battery handle
148	- mini 3000 rechargeable handle
148	- mini NT table charger
142	- NT4 table charger
149	- Plug-in transformer for ML4 LED HeadLight/OMEGA 500

Page	Index
160	Projectors F.O.
160	- HK 7000
160	Projector Accessories
160	- Fiber Optic Cable
160	- Wheeled Stand
108–112	Proctology
108–110	- Anoscopes
108	- Instrument Head for UniSpec Tubes
109	- Proctoscope / Anoscope Illumination Head
111	- Proctoscope F.O.
112	- Sigmoidoscope
114	- Sigmoidoscope/Proctoscope kit RE 7000
114	- Sigmoidoscope/Proctoscope kit BETA4 USB
109	- Tubes with straight direction of view
109	- Tubes with 90° direction of view
108	- UniSpec Disposable Tubes
113	Proctology Accessories
113	- Insufflation Bulb
113	- Sponge Holder
	R
	Rechargeable Batteries
141	- Li-ion Rechargeable Battery
142	- M3Z Rechargeable Battery
148	- NiMH 2Z Rechargeable Battery
056/060	Retinometer LAMBDA 100
053	Retinoscope BETA 200
060	Retinoscope BETA 200, Accessories
072	- Trainer, Model Eye
	S
128	S-FRAME Spectacle Frame
073	Scleral-Depressor
032	Slit Illumination Head G100/G100 LED (Veterinary)
037–040	- Slit Illumination Head G100/G100 LED Accessories
	Specula
018	- mini 3000 F.O., mini 3000 with direct illumination
018	- Otoscopes with F.O. illumination
024	- Otoscopes with direct illumination
037–039	- Veterinary Diagnostic Instruments
102–104	Sphygmomanometers GAMMA latex-free
103	- GAMMA G5
102	- GAMMA G7
103	- GAMMA GP
104	- GAMMA XXL LF
105	Sphygmomanometer Accessories
105	- Cuffs GAMMA latex-free
106	Stethoscopes GAMMA latex-free
	T
025	Tip-Adaptor (for AllSpec)
026/027	Tongue-Blade Holder
	V
032	Veterinary, Slit Illumination Head

HEINE – Simply the best guarantee in the market!

HEINE offers you a guarantee which reflects the quality and reliability of our products. All products conform to the European Directive 93/42/EEC or the Medical Devices Regulation (MDR) concerning medical devices. They carry the **CE** mark.

Standard 5 Year Parts and Labor Guarantee

All HEINE products carry a 5 year parts and labor guarantee against defects in design, materials and workmanship.* Excluded from this are wear and tear, single-use items and consumables as well as rechargeable batteries and exposed and accessible control and supply cables with corresponding plug connection.



For details please visit www.heine.com.

Product-related guarantee conditions:

HEINE iC1 Dermatoscope and HEINE iC2 Funduscope. HEINE grants a 5 year guarantee only if the product has been purchased through an authorized dealer and is registered within 6 months from the date of purchase. To register please go here www.heine.com/en_GB/register/. The guarantee provides coverage against defects in design, materials and workmanship. Rechargeable battery and mobile phone case are not included in this guarantee.

HEINE Fibre Optic Laryngoscopes. HEINE Classic+ F.O. laryngoscope blades have an extended guarantee of 10 years on the faultlessness from material, processing and construction defects, for blades manufactured on or after 01/01/2020 or with serial number SN1105228000 or higher.

The 10-year guarantee ensures that output will not fall below 1,000 lux, even after repeated autoclave cycles (based on 4,000 cycles at 134 °C/4 min.)

The guarantee covers HEINE Classic+ laryngoscope blades which are cleaned, disinfected and/or sterilised in accordance with the instructions for reprocessing.

HEINE Classic+ F.O. laryngoscope blades manufactured prior to 01/01/2020 and Modular+ blades are subject to the same guarantee conditions provided above, with a 5-year guarantee period.

HEINE GAMMA G-Series. HEINE GAMMA G7 and HEINE GAMMA XXL LF Sphygmomanometers carry an extended guarantee of 10 years for the faultlessness from material, processing and construction defects. Wear parts and metrological control are not included.

HEINE EN 200 BP Automatic Digital Sphygmomanometer. HEINE grants a 5 year guarantee on the measurement system. Wear parts and metrological control are not included.

***HEINE Cube.** For the Cube HEINE grants a guarantee of 2 years.