

## The new HEINE Premium Service: Customised laryngoscope blades with an individual inscription.

HEINE blades are in a class of their own with their 10-year guarantee and smoothSURFACE design. To ensure that they always end up where they belong following hygienic reprocessing, all Classic+ Macintosh and Paed blades can now be inscribed with the name of the department, faculty or clinic.

Our laser marking is specifically adapted to the blades' surface:

- The hygienic reprocessing options remain unaffected.
- The 10-year guarantee still applies – naturally.



This is how our HEINE Premium Service works:

We offer **inscription for our best-selling blades**, the Classic+ Macintosh blades, plus our special Classic+ Paed children's blades, so we cover the entire range of blade sizes:

- Classic+ Macintosh blades (all sizes: 0, 1, 2, 3m, 3, 4 and 5)
- Classic+ Paed blades for children (all sizes: 0 and 1)

## Inscription options:

- Free selection of a single line of text with **30 characters** (incl. spaces) from the ASCII characters (international ISO 8859-1 standard; see Table 1).

**Table 1: ASCII characters**

[Blank space]	!	“	#	\$	%	&	'	(	)	*	+	,	-	.	/
0	1	2	3	4	5	6	7	8	9	:	;	<	=	>	?
@	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
P	Q	R	S	T	U	V	W	X	Y	Z	[	\	]	^	_
`	a	b	c	d	e	f	g	h	i	j	k	l	m	n	o
p	q	r	s	t	u	v	w	x	y	z	{		}	~	

- The **font** (Arial), **font size** (approx. 2 – 4 mm character height, depending on the blade size) and **position** permit good legibility and **are standardised** by HEINE (see image).
- For regulatory reasons, we cannot add **text with critical, misleading or protected content**. We cannot inscribe text with content e.g. relating to deviating cleaning instructions, a different intended use, company names, foreign brand and product names, etc.
- We offer **inscription only for new orders**. Unfortunately, we cannot offer this service for products that have already been delivered or used.

Your local HEINE contact person will be delighted to give you more details.

[www.heine.com/en/find-a-distributor](http://www.heine.com/en/find-a-distributor)