

If there's something to discover, you'll be sure to see it.

The new HEINE DELTA 30 Dermatoscope.



30
YEARS
1st DERMATOSCOPE



The new HEINE DELTA 30 Dermatoscope. You could buy it because of its looks. Or because of its optics.

Not only would our new high-end dermatoscope win any beauty pageant; it stands out due to the best imaging we've ever developed. Combined with a field of view measuring a veritable 30 millimetres and unprecedented colour rendering, in LED^{HQ} of course.

Polarisation provides glare-free and non-reflective working conditions. Above all else, it is comfortable to operate with its angled, ergonomic design and intuitive control. Feel free to take a look – you're in for a real discovery.

Crystal-clear image with the new achromatic HEINE optical system

Largest field of view of real 30 mm (lens 32 mm)

LED^{HQ}
LED NOW IN HEINE QUALITY.

Excellent colour rendering due to LED^{HQ}

Ergonomic handling thanks to the angled product design

Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System**




Table charger for DELTA 30

- **Glare-free and non-reflective working conditions** thanks to polarisation
 - **Custom light intensity:** 3 degrees of luminous intensity
 - **10x magnification**
 - **Switch between polarised and non-polarised light at the touch of a button** = toggle function.
- For more information, please visit heine.com/en/Toggle

- **No cross-contamination:** autoclavable contact plate with scale
- **Maintenance-free** due to our dustproof design
- **Dioptré compensation**
- **Rechargeable battery with Li-ion technology:** no memory effect; you can safely recharge the battery regardless of its charge state
- **Wireless connection** via iPhone* and WLAN with the HEINE Cube System**

DELTA 30 with contact plate with scale, USB cord with medical approved plug-in power supply, case

K-230.28.305

Table charger for DELTA 30

X-002.99.212

· **Accessories and further articles under heine.com/en/DELTA30**

*The Apple iPhone is not included in the scope of delivery. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. **Not available in all countries.

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7, 82211 Herrsching, Tel. +49 (0) 8152-38 0
Fax +49 (0) 8152-382 02, E-Mail: info@heine.com
www.heine.com



10/19. A-000.00.228 e



We reserve the right to change specification without notice.
© 2019 HEINE Optotechnik GmbH & Co. KG. All rights reserved.