

1102

Ophthalmoskop Trainer

Ophthalmoscope Trainer

Kat.-Nr./Cat.-No.: C-000.33.010

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
03	Z-300.00.190	FEDERRING SPRING WASHER
06	Z-300.00.241	LINSENENKSCHRAUBE M1,6X3 COUNTERSUNK SCREW M1,6X3
12	Z-300.00.191	PROBIERGLÄSERVORSATZ TRIAL LENS HOLDER
26	Z-300.00.189	RÄNDELSCHRAUBE KNURLED SCREW

