

1312

Standard F.O. 4 NT Laryngoskopgriff

Standard F.O. 4 NT Laryngoskop handle

Kat.-Nr./Cat.-No.: F-007.22.894

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	F-000.22.824	GRIFFHUELSE F.FO-LARYNG.-GR. HANDLE OUTER/FO-LARYNGO-HAND
02	X-002.88.044	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V
03	Z-300.00.177	KAPPE OBEN BATTERY INSERT HEAD ASSEMBLY
04	Z-300.00.176	EINSCHUBHÜLSE BATTERY INSERT SLEEVE
05	X-007.99.383	LADEBATTERIE LI-ION L RECH. BATTERY LI-ION L
06	F-002.22.884	BODENEINHEIT STANDARD F.O.4 NT BOTTOM INSERT STAND. F.O.4 NT

