

1317

F.O. 4 SLIM NT

F.O. 4 SLIM NT

Kat.-Nr./Cat.-No.: F-007.22.804

F-008.22.804

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	F-000.22.804	GRIFFHUELSE/KL.LARYN.GR. OUTER HANDLE/COMP.LARYN.HA.
02	X-002.88.044	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V
03	F-002.22.801	EINSCH. KL.F.O. LARYN. LADEGR. INSERT SMALL F.O.RECH.HANDLE
04	F-008.22.801	LED BATT.EINSCH.F.KL.F.O.LA.GR LED BATT.INS.F.SMALL F.O.LA.HA
05	X-007.99.104	LADEBATTERIE M2Z 4 NT LI-ION RECH. BATT. M2Z 4 NT LI-ION
06	X-007.99.380	LADEBATTERIE M3Z 4 NT LI-ION RECH. BATT. M3Z 4 NT LI-ION
07	F-000.22.854	BODENEINHEIT F.O. 4 SLIM NT BOTTOM INSERT F.O. 4 SLIM NT

