

Kat.-Nr./Cat.-No.: D-001.70.205-041

| Nr.<br>Pos | Bestell-Nr.<br>Order-No. | Bezeichnung<br>Description                                  |
|------------|--------------------------|---|
| 01         | Z-300.00.072             | LINSENKOPFSCHRAUBE M 2X4<br>SCREW                           |
| 02         | Z-300.00.299             | FEDERSCHEIBE<br>SPRING WASHER                               |
| 03         | Z-300.00.298             | LUPENDECKEL<br>VIEWING WINDOW                               |
| 04         | X-001.88.056             | XHL XENON HALOGEN LAMPE 2,5V<br>XHL XENON HALOGEN BULB 2,5V |

